

# CATALOGUE



## WELCOME TO THE BIOSAN PRODUCT CATALOGUE

#### Biosan: Reliable Solutions for an Unpredictable World

For more than 30 years, Biosan has stood for one core mission: to remain reliable and credible in the storm of hi-tech, rapidly changing environments, and the constant chaos of new technologies. While the world races ahead with innovation, we focus on what matters most-ensuring your laboratory results are consistent, trustworthy, and ready for whatever comes next.

#### **Solutions for Every Generation of Scientists**

Our catalogue offers a complete range of laboratory equipment to cover every stage of sample preparation–from essential mixing, shaking, centrifuging, and incubation to advanced cell cultivation, nucleic acid isolation, and analytical instruments for final results. Each solution is designed to meet the needs of both experienced professionals and the expectations of a new, digitally native generation.

#### What's New?

- Microspin 12 Plus High speed mini centrifuge higher speed, broader application, same reliability;
- FTA-U the most powerful laboratory aspirator to handle it all;
- ES-20/80C shaker-Incubator with integrated cooling and new slider platform for convenient use and broader applications;
- TS-100/C/Smart thermo-shakers with new blocks;
- PCR cabinets now have a touchscreen. Easy and seamless operations with the same trustworthy results;
- UIS-360 precise and consistent petri dish plating.

#### **Adaptable Workflows-Manual or Automated**

Biosan supports both hands-on and automated workflows. Choose from **Assist pipettes** and **MagSorb-16** magnetic racks for nucleic acid isolation, or upgrade to fully automated solutions, such as **Biomagpure 12 Plus**, for high-throughput purification.

#### **Data You Can Trust**

With **DEN-1** densitometers, **HiPo MPP-96** microplate photometers, and the **Bioquant 96** real-time PCR system, Biosan helps you achieve results that are accurate, reproducible, and ready for publication.

#### **Driven by Real-World Needs**

Our mission is more than technology – it's about supporting you through the challenges of modern science. We constantly improve our products based on evolving industry standards and, most importantly, your feedback. Each update addresses real problems:

- · Accurate handling of small volumes
- · Reliable temperature control
- · Precise mixing of micro-quantities
- Air decontamination for safer workspaces
- · Secure storage of sensitive materials

#### **Support That Moves at Your Speed**

Behind every Biosan product is a dedicated team–people who understand the realities of laboratory work and who are ready to support you before, during, and after purchase. We offer fast, responsive service, practical training, and a library of videos and guides at <a href="https://www.biosan.lv">www.biosan.lv</a>.

#### **Our vision**

Biosan's promise is to stay reliable and credible, no matter how fast the world changes. We are committed to empowering professionals of every generation, from seasoned experts to digital natives, with innovative, dependable, and time-proven solutions.

Thank you for choosing Biosan.

Team that trusts in you and that you can trust!

Sincerely.

Vadims, Juris, Mark and Louise – Board of Biosan











## A Dividing Cell - Growth and New Purpose

Today, Biosan is entering a more challenging but also advanced and purposeful stage of development. Together with significant changes like the addition of a new ISO standard and a much deeper focus on medical diagnostics, our new identity represents a forward-looking shift. The logo now represents a next stage – a living cell in the process of division. Metaphorical emphasis for expansion, adaptability, and a thriving ecosystem.



#### This evolution reflects our commitment to:

- Dynamic growth, visible in our expanding technologies and global partnerships;
- Biological vitality, emphasizing our deep connection to life sciences;
- Purposeful development, aligned with the evolving needs of the medical community.

#### Our Updated Color System: Calm, Clean, Connected

With this transformation comes an updated color identity, designed to visually represent Biosan's values and direction today.

The new identity features:

- "Laboratory blue" a continuation of our legacy, symbolizing scientific precision, stability, and trust;
- "Medical green" representing sustainability, sterility, and our close relationship with the medical field;
- A renewed font for "biosan" honoring and respecting our history and roots.

#### More Than a Visual — A Living Identity

Our new identity is more than a change in design. It is a reflection of who we are becoming: a company that grows like a living system - adaptive, responsive, and committed to supporting health and science across borders.

We believe this transformation will be embraced with enthusiasm by everyone who trusted and relied on us for over three decades.



#### Disclaime

During the logo renewal transition period, the products you receive may differ from those shown in the images and may still feature the previous Biosan logo.





## Leadership and Quality Evolution

Under the new management, as never before, we remain a strong, reliable, and trusted partner, proven by countless returning customers across the globe, certified quality system according to ISO 9001, 13486, and 14001, and, most importantly, our team - people who trusts in you and whom you can trust!

#### **Complete Molecular Diagnostic Workflow**

Amidst COVID-19 lockdowns and logistical chaos, we managed to launch the complete molecular portfolio and remained a trusted and delivering. partner in over 100 countries.





#### **National Recognition**

We hit the Top 25 Exporting Brands in Latvia.

#### **Innovation** in Cell Cultivation

We introduced the personal bioreactor RTS-1, featuring the unique Reverse Spin® technology. Something the world has never seen before.





#### **Design Evolution**

Under the new logo, we introduced a new, unified concept of "living cell" - a foundation of life

## 2004

Joined by Grant Instruments, we are now ready more than ever to go beyond





#### **International Debut**

First appearance at Biotechnica'95 and immediate recognition from Japan, Korea and Germany.

#### **Foundation**

Vasily and Svetlana Bankovsky had a dream that they made come true. On the remains of the former Union, they saw a potential, they united people who shared a passion for helping others. Together, they made a company with development potential for decades.



2015

2012

2010

1995

1992

2025

2020

**Strategic Partnership** 

#### 4

#### **NEW PRODUCTS, ANNOUNCEMENT AND UPGRADES**

## NEW PRODUCTS

High-speed Mini-centrifuge: Microspin-12 Plus

Minicentrifuge-Vortex for PCR plates: MSC-2P

Aspirator: FTA-U

**Eight Interchangeable Blocks for Thermo-Shakers** 

**UV-cleaner boxes:** UVC/T-AR, UVC/T-M-AR, UVT-B-AR, UVT-S-AR

## **NEW PRODUCTS**

**Universal Inoculation Spinner:** UIS-360 **Laboratory table:** LT-120, LT-150, LT-180

Platform Exchange Set for Shakers-incubators:

ES-20/60, ES-20/80, ES-20/80C: P-EX

## **GENERAL LAB EQUIPMENT**

ROCKERS, SHAKERS, ROTATORS, VORTEXES, HOMOGENIZER	13
THERMO-SHAKERS	35
MINICENTRIFUGES-VORTEXES         4           FV-2400, FVL-2400N, MSC-3000, MSC-6000, CVP-2, MSC-2P         4	43
CENTRIFUGES	<del>1</del> 9
DRY BLOCK THERMOSTATS	51
WATER BATHS	59
MAGNETIC STIRRERS, OVERHEAD STIRRERS	37
BIOSAFETY AIR, SURFACE – UV-cleaner boxes, UV-cleaner recirculators	95
WATER PURIFICATION SYSTEMS	)5

DENSITOMETERS, PHOTOMETER	
ASPIRATORS, PIPETTES	
WASHERS	
BIOPROCESSING	
PERSONAL BIOREACTORS	
SHAKERS-INCUBATORS	
$CO_2$ INCUBATOR	
LAB DIAGNOSTICS	
DNA/RNA PURIFICATION	
REAL TIME PCR DETECTION	
IMMUNODIAGNOSTICS	
APPLICATIONS	
Article links	
GENERAL INFORMATION	
GENERAL INFORMATION ABOUT BIOSAN POLICY	





## Microspin-12 Plus NEW

see page 49

High-speed Mini-centrifuge

The Microspin 12 Plus is a high-quality mini- centrifuge designed for efficient separation of various components in a range of applications, including RNA/DNA extraction, cell suspension separation, and other micro quantitative analyses. It's sleek, bioform design and compact footprint make it a space-saving solution to the laboratory workspace. The centrifuge can accommodate microtubes or strip tubes, reaching a maximum speed of 16,250 RPM or 15,588 RCF.

Safety features include metal protective inserts and senclosures inside the body and lid of the centrifuge, an automatic imbalance detection system with an automatic stop function, and a lid locking mechanism, providing secure operation throughout the speed range. A sound signal alerts when centrifugation is complete.



## **NEW FEATURE**

Quick spin button







**MSC-2P** is a compact sized digital centrifuge intended to collect droplets, mix reagents and collect once more for improved PCR yield in subsequent analysis. The combination of spin-mix functions ensures fast operation, thorough mixing and repeatable results. Centrifuge rotor can accommodate 2 PCR plates at the same time, thus saving time considerably.


#### MSC-2P is possible to operate in 4 independent modes:


- Centrifuge Max. 3,500 rpm
- Vortex up to 5 min
- Centrifuge/Vortex combined two motion types
- Spin-mix-spin algorithm up to 10 cycles


The spin-mix-spin algorithm (SMS-algorithm) is designed to collect (or reset) micro volumes of reagents to the bottom of the PCR plate tubes (the first centrifugation or spin), then vortexing (mix) and re-collecting reagents (repeated spin) from the walls and cover. This repetitive algorithm of operations, aimed at reducing sample preparation errors, we call the SMS algorithm.













## TS-100, TS-100C and TS-100C Smart

Thermo-Shakers

**TS-100 and TS-100C** thermo-shakers are designed for intensive mixing of various reactions vessels in a temperature controlled environment. The **TS-100C** model of thermo-shaker differs from **TS-100** in the possibility of cooling samples down to  $+4^{\circ}$ C.

We introduce 8 new thermoblocks.

## Interchangeable Blocks for TS-100

VP-8/5	8 × 5 ml conical tubes
VP-4	4 × 50 ml conical tubes
VP-8/15	8 × 15 ml conical tubes
VP-CV-20	20 × 10 mm cuvettes (0.2–4.5 ml)
VP-32	32 × 0.5 ml microtubes
VP-CL-24	24 × 3.6–4.5 ml cryotubes
VP-CS-24	24 × 1–1.8 ml cryotubes
VP-20	ø12 mm round bottom tubes
VP-18	$20 \times 0.5 \text{ ml} + 12 \times 1.5 \text{ ml microtubes}$
VP-18/02	$20 \times 0.2$ ml microtubes +12 × 1.5 ml microtubes
VP-24N	24 × 1.5 ml microtubes
VP-24	24 × 2 ml microtubes
VP-96A	96-well unskirted or semi-skirted microplate (0.2 ml) for PCR or 12 × 8–0.2ml strips or 96 tubes of 0.2 ml

## Similar Blocks available also for TS-100C

## **NEW BLOCKS**



VP-8/5 8 × 5 ml conical tubes



VP-4 4 × 50 ml conical tubes



VP-8/15 8 × 15 ml conical tubes



VP-CV-20 20 × 10 mm cuvettes (0.2-4.5 ml)



VP-32 32 × 0.5 ml microtubes



VP-CL-24 24 × 3.6–4.5 ml cryotubes



VP-CS-24 24 × 1–1.8 ml cryotubes



VP-20 20 × ø12 mm round bottom tubes

more info on page 38





FTA-U NEW see page 116

## Universal Aspirator with Trap Flask

FTA-U is designed for the precise aspiration and removal of alcohol, buffer solutions, and other liquids from reaction vessels, including applications such as DNA/RNA extraction, ELISA plate washing, and cell culture media exchange. The device is suitable for use with microtubes, plates, and flasks, and can be configured with either a 2 L or 4 L polypropylene trapping flask.

Suction power is regulated continuously between -200 mbar and -950 mbar, with a maximum aspiration speed of 12 l/min (air) and a suction flow rate up to 47 ml/s (liquid). Aspiration intensity is controlled by a linear regulator on the front panel, and operation status is clearly indicated with a light ring around the control knob. Overflow protection is provided by an invasive level sensor, which stops the pump and triggers sound and light alarms to prevent accidental spills.

**FTA-U** is equipped with a hydrophobic microbiological air filter (pore size  $0.027 \mu m$ ) that removes up to 99.99% of bacteria and viral particles from the airstream, ensuring contamination-free operation and laboratory safety. The standard set includes the MA-U universal adapter for single-use 200  $\mu L$  or 1000  $\mu L$  tips, with optional accessories such as an 8-channel aspiration tip (MA-8), HAS-1 hand operator for ergonomic handling, and a foot switch for hands-free operation.











**UVT-S-AR** 

## UVC/T-AR, UVC/T-M-AR, UVT-B-AR and UVT-S-AR see page 98

PCR cabinets

DNA/RNA UV-cleaner cabinets (UVC/T-AR, UVC/T-M-AR, UVT-B-AR, UVT-S-AR) are designed for dependable, contamination-free operations with DNA samples, supporting stable results in any laboratory environment. A new addition is the touchscreen control panel, which allows straightforward cabinet operations, programming of UV exposure cycles, access to the User Manuals, Service access and real-time monitoring of lamp life-making daily operation easy, precise, consistent, and reproducible

Biosan UV-cleaner cabinets are recommended for operations with DNA/RNA, providing a stable and trustworthy environment so you don't have to worry about direct and cross-contamination.

## LT-120, LT-150 and LT-180

Laboratory table NEW see page 104

LT-120, LT-150 or LT-180 laboratory table is the perfect solution for any research or laboratory facility. It features a chemically resistant high pressure laminate worktop, making it suitable for use in a variety of laboratory settings. The worktop is easy to clean and maintain, ensuring a hygienic work environment. The table is also adjustable in height, which allows you to customize the table to your preferred working height for added comfort and convenience. In addition to its standard features, this laboratory table also offers an optional drawer with an anti-slippery surface and a retractable block of electrical sockets. LT series table is an all-in-one solution for any laboratory or research facility, providing a sturdy, reliable and practical work surface that can stand up to the demands of daily use.



UVC/T-AR



UVC/T-M-AR



UVT-B-AR







P-EX, Platform Exchange Set

**ES-20/80C** NEW see page 130 Shaker-incubator with cooling

Orbital Shaker–Incubator with cooling **ES-20/80C** for biotechnological and pharmaceutical laboratories is a next-generation, professional category equipment.

**ES-20/80C** duplicates the functionality of ES-20/80 and also uses a Peltier element to cool the camera. A built-in heat-resistant brushless fan provides precise temperature distribution inside the chamber (from 12.5 °C below ambient up to +80 °C). Additionally, excellent sample temperature uniformity of  $\pm 0.2$  °C at 37 °C is achieved.







P-EX NEW

Platform Exchange Set (Slider + Tray)

P-EX, Platform Exchange Set allows for faster replacement of flat platforms. Thanks to the platform sliding system, it is now possible to install 4  $\times$  2l flasks in ES-20/60, ES-20/80, ES-20/80C shakers-incubators. The platform exchange set can accommodate UP-168, P-6/1000, P-9/500, P-15/250, P-30/100, HSP-6/1000, HSP-9/500, HSP-15/250, HSP-30/100 platforms.



# **UIS-360** NEW see page 33 Universal Inoculation Spinner

The **UIS-360** Universal Inoculation Spinner is a versatile laboratory device designed to streamline microbial plating and ensure consistent results. It provides smooth, uniform rotation for spreading samples across agar surfaces, improving both efficiency and reproducibility in microbiological workflows. With an adjustable speed range of 10–300 RPM – the broadest in its class – the **UIS-360** accommodates various protocols from gentle inoculum spreading to vigorous. Its compact footprint and lightweight build make it easy to handle or even take into the field (it runs on just 3.7 W, so it can be powered by a standard 12 V battery pack). The inclusion of a foot-switch for hands-free control and a protective shield for splash protection further enhance user comfort and safety.



## **Key Features and Advantages:**

- Universal Compatibility fits both square (125×125 mm) and round Petri dishes (up to Ø105 mm) out-of-the-box, and even up to Ø150 mm with an optional adapter.
- Widest Speed Range adjustable speed control from 10 to 300 RPM – offering the widest range among similar plate spinners.
- Continuous and Hands free Operation: supports continuous rotation for extended procedures, with an included foot-switch for hands-free control.
- Enhanced Safety equipped with a protective polycarbonate face shield that guards against accidental splashes and aerosol droplets.
- Compact & Portable Design measures only 170×190 mm on the lab bench, weighs ~0.85 kg, and can even be powered by a 12 V DC power bank.
- Delivers smooth, uniform rotation for even spreading of samples, improving reproducibility between plates.


## **Applications:**

- Microbial culture plating Uniform spreading of inoculum on agar surfaces for colony isolation.
- Diagnostic testing Preparation of plates for pathogen detection and antibiotic susceptibility testing.
- Quality control in food and beverages Microbial screening for contamination detection.
- Pharmaceutical research Sterility testing and inoculation preparation for drug development.
- Environmental monitoring Efficient spreading of soil, water, or air samples for microbial analysis.
- Educational training Teaching microbiological plating techniques with a reliable and user-friendly device.



## MIXING DEVICES:

## **ROCKERS, SHAKERS,** ROTATORS, VORTEXES, **HOMOGENIXER**



Programmable rotator

# DESCRIPTION

## MR-1, Mini Rocker-Shaker

Mini Rocker-Shaker **MR-1** provides regulated gentle rocking motion of the platform and is ideal for mini gel destaining after electrophoresis, conducting Northern, Southern and Western blot analysis.

Shaker is a compact, noiseless device designed for personal use. Drive and brushless motor allow continuous mixing up to 7 days and ensures reliable, trouble-free operation for more than 2 years.

Non-slip, temperature resistant, silicone mat located on the rocker's platform provides a stable position for vessels during shaking. Optional dimpled PDM mat fixes tubes of different sizes.

The unit is designed for operation in cold rooms, incubators (excluding  $CO_2$  incubators) and closed laboratory rooms at ambient temperature from  $+4^{\circ}C$  to  $+40^{\circ}C$  in a non-condensing atmosphere and maximum relative humidity 80% for temperatures up to 31°C decreasing linearly to 50% relative humidity at  $40^{\circ}C$ .

#### ACCESSORIES FOR THE STANDARD PLATFORM:

Optional dimpled mat PDM prevents different size tubes from rolling around the platform.



## MR-12, Rocker–Shaker

**MR-12** Rocker–Shaker provides both soft and intensive mixing of solutions or nutrient media in vessels or plastic bags placed on the platform. Adjustable speed and platform tilt angle allows setting parameters for optimal solution transfer and mixing.

The device is ideal for gel destaining after electrophoresis and homogenisation of bioextraction media. It is optimal for biomolecule hybridisation on strips and staining/destaining procedures. When installed inside a bioincubator it is ideal for growing cells and cell cultures in disposable plastic reactor-bags (working volumes up to 10 litres, media volumes up to 5 litres).

The unit is designed for operation in cold rooms, incubators (excluding  $CO_2$  incubators) and closed laboratory rooms at ambient temperature from +4°C to +40°C in a non-condensing atmosphere and maximum relative humidity 80% for temperatures up to 31°C decreasing linearly to 50% relative humidity at 40°C. Low voltage external power supply (12V) provides electrical safety in humid environment



## MR-1 and MR-12, Rocker-Shakers

	MR-1	MR-12	
Mixing frequency range	1–30 oscill./min	1–99 oscill./min (increment 1 oscill./min)	
Fixed tilt angle	7° (fixed)	0°-10° (increment 1°) (for 1-50 oscill./min) 10° (for 51-99 oscil./min)	
Max. continuous operation time	168	8 h	
Digital time setting	1 min-24 h/non-stop	1 min–99 h 59 min (increment 1 min)/non-stop	
Timer sound signal	-	yes	
Non-slip silicone mat is supplied as standard	215 × 215 mm	480 × 380 mm	
Maximum load	1 kg	5 kg	
Display	LED	LCD, 2 × 16 signs	
Platform working area	215 × 215 mm	480 × 380 mm	
Overall dimensions (W×D×H)	220 × 205 × 120 mm	430 × 480 × 210 mm	
Weight	2.1 kg	11.9 kg	
Input current/power consumption	12 V, 320 mA/3.8 W	12 V, 1.1A/13 W	
External power supply	Input AC 100–240 V, 50/60 Hz; Output DC 12 V		


PDM, dimpled mat



MR-1 with PDM dimpled mat



MR-12



## ORDERING INFORMATION:

Cat. number

BS-010152-AAG

MR-1 with standard platform Bio PP-4S

MR-12 with standard platform PP-480

BS-010130-AAI

Optional accessories: for MR-1:

PDM PDM, dimpled mat

## 3D, Sunflower Mini Shaker

"Sunflower" **3D** Mini–Shaker provides adjustable three-dimensional smooth rotation of the platform and is designed for mixing blood samples, minigel staining and destaining, sample washing, blot hybridisation reactions. Mini–Shaker is a compact device with low energy consumption. The use of direct drive and brushless motor allows continuous mixing up to 7 days and ensures reliable, trouble-free operation for many years. Non–slip, temperature resistant, silicone mat located on the shaker's platform provides a stable position for vessels during shaking. The platform is suitable for placing a versatile dimpled PDM mat for different size tubes.

Mini–Shaker can be used in cold rooms or incubators, operating at ambient temperature range +4°C to +40°C.



## Multi Bio3D,

## Programmable mini-shaker («Sunflower» type)

Programmable mini-shaker **Multi Bio 3D** is designed for various applications: hybridization reactions, cell growing, gel washing, soft extraction and homogenisation of biological components in solutions.

Multi Bio 3D provides realization of several types of motion in one module. This option of Biosan instruments essentially extends possibilities and enhances the efficiency of preparation of test samples as well as allows selecting the mixing type according to individual requirements.

Microprocessor control allows performing

Orbital 3D rotation of the platform and but also
Reciprocal 3D motion (of ping-pong type) and
Soft vibrating rocking. These three motion types can be performed separately, pairwise and in cycles, periodically repeating the sequence of three motion types. The shaker is designed for laboratories with increased demands for the quality of mixing, extraction and cell growing processes.

Non-slip, temperature resistant, silicone mat located on the shaker platform provides a stable position for vessels during shaking. Optional dimpled PDM mat fixes tubes of different sizes.

Programmable shaker can be used in cold rooms or incubators, operating at the ambient temperature range +4°C to +40°C.



**Product** 

video

# **3D** Mini-Shaker and **Multi Bio 3D**, Programmable 3D shaker («Sunflower» type)

	3D	Multi Bio 3D	
Speed control range (orbital and reciprocal motion)	5–60 rpm	1–100 rpm	
2 Turning angle (reciprocal motion)	-	0–360° (increment 30°)	
3 Rocking angle (vibro motion)	-	0–5° (increment 1°)	
Fixed tilt angle	7	70	
Orbit	-	22 mm	
Platform working area	215 × 215 mm		
Non-slip silicone mat is supplied as standard			
Maximum continuous operation time	168 h		
Time setting range for <b>1 2</b>	-	0-250 s	
Time setting range for <b>3</b>	-	0-5 s	
Number of cycles	-	0–125 times	
Timer sound signal	-	yes	
Maximum load	1 kg		
Overall dimensions (W $\times$ D $\times$ H)	235 × 235 × 140 mm		
Weight	1.2 kg	1.8 kg	
Input current/power consumption	12 V, 260 mA/3.1 W	12 V, 380 mA/4.6 W	
External power supply	Input AC 100-240 V, 50/60 Hz; Output DC 12 V		

## Accessories for the standard platform:

Optional dimpled mat PDM prevents different size tubes from rolling around the platform









#### ORDERING INFORMATION:

Cat. number 🔭

 ${f 3D}$  with stand. platform  ${f Bio\ PP\text{-}4S}$ 

BS-010151-AAG

Multi Bio 3D with stand. platform Bio PP-4S

BS-010125-AAG

Optional accessories:

PDM dimpled mat PDM

## PSU-10i, Orbital Shaker

Shaker **PSU-10i** provides regulated orbital motion of the platform and is designed for use both in small specialized biotechnological laboratories and in large multidisciplinary laboratories: a choice of five (5) interchangeable platforms provides the possibility of performing various procedures and techniques.

Shaker **PSU-10i** incorporates a direct drive system, a brushless motor with a guaranteed service life of up to 35,000 hours and an automatic loading balancing system. These innovations allow for continuous mixing up to 7 days, ensure reliable, trouble-free operation for more than 2 years and significantly expand the device performance range in both high and low limits.

The unit is designed for operation in cold rooms, incubators (excluding  $CO_2$  incubators) and closed laboratory rooms at ambient temperature from +4°C to +40°C in a noncondensing atmosphere and maximum relative humidity 80% for temperatures up to 31°C decreasing linearly to 50% relative humidity at 40°C.



## PSU-20i, Orbital Shaker

Shaker **PSU-20i** provides three motion types: **① Orbital, ② Reciprocal** and **③ Vibrating,** which can be performed separately, pairwise and sequentially in repeated cycles.

Shaker is designed for applications both in small specialized laboratories and in large multidisciplinary laboratories. **PSU-20i** is an ideal instrument for laboratories researching biopharmaceutics and biomedicine.

Shaker **PSU-20i** is noiseless and reliable in operation, incorporates a direct drive system and brushless motor with a guaranteed service life up to 35,000 working hours. Direct drive and brushless motor allows for continuous mixing for up to 7 days and ensures reliable operation for more than two years.

A choice of nine different interchangeable platforms provides the possibility of performing various procedures and techniques. Special attention should be paid to a multilevel platform, which allows accommodation of a large number of various microplates, Petri dishes, cultural bags and other low containers.

The unit is designed for operation in cold rooms, incubators (excluding  $CO_2$  incubators) and closed laboratory rooms at ambient temperature from +4°C to +40°C in a noncondensing atmosphere and maximum relative humidity 80% for temperatures up to 31°C decreasing linearly to 50% relative humidity at 40°C.



## PSU-10i and PSU-20i, Orbital Shakers

	PSU-10i	PSU-20i	
Multi-motion	-	yes	
Speed control range*	50-450* rpm (increment 10 rpm)	20–250* rpm (increment 5 rpm)	
Digital speed control	yes		
Max. continuous operation time	168 h		
Orbit	10 mm 20 mm		
Digital time setting	1 min-96 h/non-stop		
Timer sound signal	yes		
Maximum load	3 kg	8 kg	
Overall dimensions (W×D×H)	255 × 255 × 100 mm 410 × 410 × 130 mm		
Weight	3.4 kg 11.7 kg		
Input current/power consumption	12 V, 800 mA/9.6 W 12 V, 3.2 A/40 W		
External power supply	Input AC 100–240 V, 50/60 Hz; Output DC 12 V		

<sup>\* -</sup> max. speed depends on the load and vessels' shape

Platform Bio PP-4 for PSU-10i



Platform PP-20/4 for PSU-20i



Platform Bio PP-4 for PSU-10i



Platform P-6/250 for PSU-10i

ORDERING INFORMATION:



<b>PSU-10i,</b> Shaker without platform	BS-010144-AAN
<b>PSU-20i,</b> Shaker without platform	BS-010145-ACI

Cat. number 📜

PSU-20i motion types	Description	Speed range	Turning angle	Motion timer*	Digital time setting
① 〇 Orbital	Orbital motion with an option of shifting direction	20–250 rpm	-	0-250 s	
2 Reciprocal	Orbital motion with shifting direction of rotation	20–250 rpm	0–360° (30° increment)	0-250 s	1 min-96 h (increment 1 min) or non-stop
3 S Vibrating	High speed, low amplitude motion	-	0–5° (1° increment)	0-5 s	,

<sup>\* –</sup> for switching to the next motion in the cycle

## Platforms for PSU-10i and ES-20

Platform		Description	Dimensions (Working area)	Cat. number
UP-12 Used on PSU-10i, ES-20	The same of the sa	Universal platform with adjustable bars for different types of flasks, bottles and beakers with silicone mat	285 × 220 × 40 mm (270 × 185 × 40 mm)	BS-010108-AK
<b>Bio PP-4</b> Used on PSU-10i		Flat platform with silicone mat for Petri dishes, culture flasks, agglutination cards	255 × 255 mm (230 × 230 mm)	BS-010116-AK
PP-4 Used on ES-20, PSU-10i		Metallic flat platform with silicone mat for Petri dishes, culture flasks, agglutination cards	220 × 220 mm (215 × 215 mm)	BS-010108-BK
SPM		Double-sided adhesive mat as an alternative for regular flask clamps (for PP-4)	220 × 220 mm	BS-010111-BK
<b>P-12/100</b> Used on PSU-10i, ES-20		Platform with clamps for flasks, 100–150 ml (12 places)	250 × 190 mm (250 × 190 mm)	BS-010108-EK
<b>P-6/250</b> Used on PSU-10i, ES-20		Platform with clamps for flasks, 250–300 ml (6 places)	250 × 190 mm (250 × 190 mm)	BS-010108-DK
<b>P-16/88</b> Used on PSU-10i, ES-20		Platform with spring holders for up to 88 tubes up to 30 mm diameter (e. g. 10 ml, 15 ml, 50 ml tubes)	275 × 205 × 75 mm (275 × 205 × 75 mm)	BS-010116-BK

## SPML, SPM, Double-sided adhesive strips and mat

Convenient alternative to traditional steel holders, an easy way to fix tubes, plates, flasks and other laboratory wares on platforms for cultivation, incubation and mixing. Two size options are offered **SPML** can be used with UP-168 platform on Biosan orbital shaker PSU-20i and in ES 20/80, ES-20/80C, ES 20/60 orbital shakers.

**SPM** is compatible with PP-4 platform, which fits both on PSU-10i orbital shaker and in ES-20 Shaker-Incubator.

Made of polyurethane with adhesive, simple to clean and durable, able to withstand up to 1,000 times placement/removal or 12 months of use. Additional information about temperature, working volume and speed limitations available in the user manual.

SPM on PP-4 platform





SPML Size (L×W×H)	390 × 80 × 3 mm (double sided 1.5 mm PU with PET adhesive)
SPM Size (L×W×H)	$210 \times 210 \times 3$ mm (double sided 1.5 mm PU with PET adhesive)
Colour	transparent
Duration of use	up to 1,000 times placement/ removal or 12 months
Temperature range	+4C° to +80C°
Shaking speed	0–400 rpm

SPML on UP-168 platform





## ORDERING INFORMATION:

## Platforms for **PSU-20i** and **ES-20/60**

Platform	Description	Dimensions (Working area)	Cat. number
UP-330 Used on PSU-20i	Universal platform with adjustable bars for different types of flasks, beakers	345 × 430 × 105 mm (300 × 400 × 80 mm)	BS-010145-AK
<b>P-30/100</b> Used on PSU-20i, ES-20/60	Platform with 30 clamps for 100–150 ml flasks	360 × 400 mm (360 × 400 mm)	BS-010135-BK
<b>P-16/250</b> Used on PSU-20i, ES-20/60	Platform with 16 clamps for 250–300 ml flasks	360 × 400 mm (360 × 400 mm)	BS-010135-CK
<b>P-9/500</b> Used on PSU-20i, ES-20/60	Platform with 9 clamps for 500 ml flasks	360 × 400 mm (360 × 400 mm)	BS-010135-AK
<b>P-6/1000</b> Used on PSU-20i, ES-20/60	Platform with 6 clamps for 1,000 ml flasks	360 × 400 mm (360 × 400 mm)	BS-010135-DK
PP-400 Used on PSU-20i, ES-20/60, ES-20/80	Flat platform with non-slip silicone mat	360 × 400 mm (360 × 400 mm)	BS-010135-FK
<b>UP-168</b> Used on PSU-20i, ES-20/60, ES-20/80	Universal platform for different flasks (Clamps ordered separately)	360 × 400 mm (360 × 400 mm)	BS-010135-JK
FC-50 FC-100 FC-250 FC-500 FC-1000 FC-2000	Clamp for 50, 100, 250, 500, 1000, 2,000 ml flask (for UP-168)	Ø 50 mm Ø 65 mm Ø 85 mm Ø 105 mm Ø 130 mm Ø 165 mm	BS-010126-MK BS-010126-HK BS-010126-JK BS-010126-LK BS-010126-IK BS-010126-NK
SPML	Set of 3 double-sided adhesive strips as an alternative for regular flask clamps (for UP-168)	390 × 80 × 3 mm	BS-010135-MK
TR-21/50	Test tube rack for 50 ml with 21 drillings (for UP-168)	340 × 124 mm (2 per platform)	BS-010135-KK
TR-44/15	Test tube rack for 15 ml with 44 drillings (for UP-168)	340 × 124 mm (2 per platform)	BS-010135-LK
<b>PP-20/4</b> Used on PSU-20i	Four-level flat platform with non-slip rubber mat	380 × 480 × 440 mm (355 × 455 mm ×4)	BS-010126-EK
<b>PP-20/3</b> Used on PSU-20i	Three-level flat platform with non-slip rubber mat	380 × 480 × 300 mm (355 × 455 mm ×3)	BS-010126-DK
<b>PP-20/2</b> Used on PSU-20i	Two-level flat platform with non-slip rubber mat	380 × 480 × 160 mm (355 × 455 mm ×2)	BS-010126-CK
PP-20 Used on PSU-20i	One-level flat platform with non-slip rubber mat	380 × 480 mm (355 × 455 mm)	BS-010126-BK
P-EX For ES-20/60, ES-20/80, ES-20/80C	Platform Exchange Set (Slider + Tray).  Now possible to install 4 × 2l flasks.  Accommodate: UP-168, P-6/1000, P-9/500, P-15/250, P-30/100, HSP-6/1000, HSP-9/500, HSP-15/250, HSP-30/100 platforms.		BS-010173-CK

## MPS-1, High-Speed Multi Plate Shaker



High–Speed Multi Plate Shaker **MPS-1** can be used in virtually any application by providing adjustable mixing of reagents in microtest plates, PCR plates, deepwell plates and test tubes (shaking tubes 0.2 to 2 ml and vortexing any volume up to 50 ml).

The shaker is compact and user-friendly. The shaker is ideal for personal use.

MPS-1 features a head for vortexing a single tube.

The unit is designed for operation in cold rooms, incubators (excluding  $CO_2$  incubators) and closed laboratory rooms at ambient temperature from +4°C to +40°C in a non-condensing atmosphere and maximum relative humidity 80% for temperatures up to 31°C, decreasing linearly to 50% relative humidity at 40°C. Low voltage external power supply (12 V) provides electrical safety in a humid environment.

MPS-1 features Pulse Mode mixing function that works on the principle of giving a periodic impulse: the tube is accelerated to the set speed, holds it for 3 seconds and then drops the speed to zero. This motion is repeated until the timer runs out. This method provides a constant state particle resuspension inside a tube, as the acceleration is always changing. The advantage of this method is the high throughput of mixed samples compared to vortexing a single tube.



#### **Features**

- Speed control range 300-3,200 rpm
- · Stable mixing with 3 mm orbit
- · Five mixing presets
- · Pulse Mode mixing function
- · Quiet operation low noise at maximum speed
- Universal platform holder for Deepwell plates and Microtest plates
- Additional four platforms for semiskirted and unskirted PCR plates 200 µl as well as for tubes from 0.2 to 2 ml



**P-02/96,** latform for semi-/unskirted PCR plate 200 μl

## MPS-1, High-Speed Multi Plate Shaker

Vortexing a 50 ml tube



Vortexing a 15 ml tube



Deepwell plate 96/1000 µl



Microtest plate 200 µl



Deepwell plate 96/500 µl



Mixing Speed control range	300-3,200 rpm
Platform options:	
<ul> <li>For semi-\unskirted PCR plate or 96 microtest tubes 0.2 ml</li> </ul>	P-02/96
- For 24 microtest tubes 1.5-2 ml	P-2/24
- For 32 microtest tubes 0.5 ml	P-05/32
- For 24 microtest tubes 0.5 ml and 48 microtest tubes 0.2 ml	P-02/05
<ul> <li>Universal platform for deepwell plates, 96-well microtest plates (U, V or flat bottomed), 384-well microtest plates</li> </ul>	

	presets:

VORTEX	3,200 rpm
HARD	2,600 rpm
MEDIUM	1,800 rpm
SOFT	1,000 rpm
CUSTOM	adjustable rpm

Features a Pulse Mode mixing function	
Features a Vortex function	
Maximum load	0.3 kg
Mixing Orbit	3 mm
Acceleration time to maximum speed	5 s
Digital time setting	0-60 min (15 s increment)/non-stop
Timer sound signal	yes
Maximum continuous operation time	8 h
Noise level, not more	65 dB
Weight	5.1 kg
Overall dimensions (W×D×H)	225 × 215 × 150 mm
Input current/power consumption	12 V, 800 mA / 10 W
External power supply	Input AC 100-240 V 50/60 Hz; Output DC 12 V

ORDERING INFORMATION:	Ca	it. num

ORDERING INFORMATION:	Cat. number 👈
MPS-1, Multi Plate Shaker with built-in universal platform	BS-010216-A03
MPS-1, Multi Plate Shaker with built-in universal platform	
and set of 4 platforms (P-02/96, P-2/24, P-05/32, P-02/05)	BS-010216-A11

Optional platforms:		Cat. number
<b>1</b> P-02/96	For semi-/unskirted PCR plate or 96 microtest tubes 0.2 ml	BS-010216-CK
<b>2</b> P-2/24	For 24 microtest tubes 1.5–2 ml	BS-010216-AK
<b>3</b> P-05/32	For 32 microtest tubes 0.5 ml	BS-010216-BK
<b>4</b> P-02/05	For 24 microtest tubes 0.5 ml and 48 microtest tubes 0.2 ml	BS-010216-DK





2 Platform P-2/24



A	Platform	P-05/32



4 Platform P-02/05



SPECIFICATIONS

## PSU-2T, Mini-Shaker

Mini–Shaker **PSU-2T** is designed for immunoassays and provides adjustable mixing of reagents in microplates. The device ensures smooth movement of the platform even at low speeds.

Shaker is a compact and user-friendly device. It takes up little space on a desk and is ideal for personal use. Direct drive and brushless motor allow continuous mixing up to 7 days and ensures reliable, trouble-free operation for more than 2 years. Display of the device switches between time and speed readings.

The unit is designed for operation in cold rooms, incubators (excluding  $CO_2$  incubators) and closed laboratory rooms at ambient temperature from  $+4^{\circ}C$  to  $+40^{\circ}C$  in a non-condensing atmosphere and maximum relative humidity 80% for temperatures up to 31°C decreasing linearly to 50% relative humidity at  $40^{\circ}C$ .

Speed control range	150–1,200 rpm
Digital time setting	1 min-24 h/non-stop
Digital setting and control of time	and speed
Max. continuous operation time	168 h
Direct drive mechanism	
Orbit	2 mm
Overall dimensions (W×D×H)	220 × 205 × 90 mm
Weight	2 kg
Input current/ power consumption	12 V, 280 mA/3.4 W



External power supply

Output DC 12 V

Cat. number

Input AC 100-240 V, 50/60 Hz;

**PSU-2T** with standard platform IPP-2 BS-010155-AAG

## Optional platforms

**IPP-4** BS-010102-AK





A Platform IPP-2



B Platform IPP-4



## Platforms for microtest plates:

for 4 microtest plates

**184**  $\times$  132 mm for 2 microtest plates **1919-4** (optional platform) **266**  $\times$  170 mm



## Multi Bio RS-24 and Multi RS-60, rotators






It is possible to choose the position of tubes for rocking motion – horizontal or vertical. The platform does not make an additional revolution before stopping in the horizontal plane.

Product video Programmable Rotators performs several motion types in one module. Microprocessor control allows performing not only **1** Vertical overhead rotation of the platform, but also **2** Reciprocal rotation (rocking motion) as well as **3** Vibration. These three motion types can be performed separately, pairwise and in cycles, periodically repeating the sequence of three motion types. Multi–Rotation option of Biosan instruments substantially expands possibilities and enhances the efficiency of sample preparation for the examined materials and allows adjusting the mixing procedure according to the individual tasks.

Programmable Rotators can be used for variety of applications in modern life science laboratories: for hybridisation reactions, cell growing, soft extraction and homogenisation of biological components in solutions, as well as for reactions of binding and washing of magnetic particles.

**Multi Bio RS-24** and **Multi RS-60** are designed for operation in cold rooms, incubators (excluding  $CO_2$  incubators) and closed laboratory rooms at ambient temperature from +4°C to +40°C in a non-condensing atmosphere and maximum relative humidity 80% for temperatures up to 31°C decreasing linearly to 50% relative humidity at 40 °C. Low voltage external power supply (12 V / 24V) provides electrical safety in a humid environment.



Programmable Rotator provides 3 rotation types and Pause:

- 1 Vertical overhead rotation
- Reciprocal rotation (rocking motion)
- 3 Vibro
- Pause

## Multi Bio RS-24 and Multi RS-60, rotator

	Multi Bio RS-24	Multi RS-60		
Vertical overhead rotation:				
Speed control range 1–100 rpm (increment 1 rpm)				
Vertical rotation movement	36	0°		
Time setting range	0-25	50 s		
2 Reciprocal rotation (rocking motion):				
Speed control range	1–100 rpm (inc	rement 1 rpm)		
Tilt angle range	1–90° (incr	rement 1°)		
Time setting range	0-250 s			
3 Vibro:				
Tilt angle range	0–5° (increment 1°)			
Pause/Vibro time setting range	0-5 s			
GENERAL SPECIFICATIONS:				
Digital time setting	1 min-24 h/non-sto	p (increment 1 min)		
Timer sound signal	ye	es		
Maximum load	0.5 kg	0.8 kg		
Overall dimensions (W $\times$ D $\times$ H)	365 × 195 × 155 mm	430 × 230 × 230 mm		
Weight	1.7 kg	3.8 kg		
Input current/power consumption	12 V, 660 mA/8 W	24 V, 750 mA/18 W		
External power supply	Input AC 100–240 V, 50/60 Hz; Output DC 12 V	Input AC 100–240 V, 50/60 Hz; Output DC 24 V		

## Optional platforms for Multi RS-60



Standard:	Capacity	Tube Volume	Tube Diameter	Cat. number
<b>1</b> PRS-48	48	1.5–15 ml	10–16 mm	BS-010118-CK
Optional:				
2 PRS-8/22	8 and 22	up to 50 and 1.5–15 ml	20–30 and 10–16 mm	BS-010118-AK
<b>3</b> PRS-14	14	up to 50 ml	20-30 mm	BS-010118-BK









**3** PRS-14



## **Optional platforms for Multi Bio RS-24**

	Opti	ional platfor	ms for wu	ti Bio K5-24
Standard:	Capacity	Tube Volume	Tube Diameter	Cat. number
<b>1</b> PRS-26	26	1.5–15 ml	10–16 mm	BS-010117-GK
Optional				
<b>2</b> PRS-5/12	5 and 12	up to 50 and 1.5–15 ml	20–30 and 10–16 mr	n BS-010117-HK
<b>3</b> PRS-10	10	up to 50 ml	20-30 mm	BS-010117-IK
PRSC-22	22	15 ml	16 mm	BS-010117-LK
PRSC-10	10	50 ml	25-30 mm	BS-010117-JK
<b>6</b> M-8/50	8	50 ml	25-30 mm	BS-010117-PK
7 RP-8/15	8	15 ml	16 mm	BS-010117-DK
8 RP-8/50	8	50 ml	25–30 mm	BS-010117-KK
9 PRS-1DP		croplates and racks for tall Thermo 3741MTX, 3742MTX, 3		BS-010149-DK
PRS-26	2 PRS-5/12	<b>3</b> PRS-1	0	4 PRSC-22
PRSC-10	6 M-8/50	9 PRS-1	DP C	Elamps on PRSC-10
			D	
RP-8/15				
8 RP-8/50				

**PRS** series platforms are equipped with universal rubber clamps for different size tube fixation; **PRSC** series platforms have metal clamps able to hold heavier solutions (e.g. soil, sand).



Cat. number

ORDERING INFORMATION:

## Bio RS-24, Mini-Rotator

Mini-rotator Bio RS-24 provides vertical rotation of the platform. The rotator is an ideal instrument for preventing blood coagulation in tubes and fulfilling of procedures of biological components extraction.

The device is simple to operate; it is designed as a lowcost solution.

The unit is designed for operation in cold rooms, incubators (excluding CO2 incubators) and closed laboratory rooms at ambient temperature from +4°C to +40°C in a non-condensing atmosphere and maximum relative humidity 80% for temperatures up to 31°C decreasing linearly to 50% relative humidity at 40°C. Low voltage external power supply (12 V) provides electrical safety in a humidenvironment.

Speed control range	5–30 rpm
Vertical rotation movement	overhead, 360°
Digital time setting	1 min–24 h/non-stop (increment 1 min)
Timer sound signal	yes
Maximum continuous operation time 8 h	
Overall dimensions (W×D×H)	325 × 190 × 155 mm
Weight	1.4 kg
Recommended load	75% of the rated volume
Input current/power consumption	12 V, 110 mA/1.3 W
External power supply	Input AC 100–240 V 50/60 Hz; Output DC 12 V

PRS series platforms are equipped with universal rubber clamps for different size tube fixation;

PRSC series platforms have metal clamps able to hold heavier solutions (e.g. soil, sand).







Vertical rotation 360°

Bio RS-24 in operation



#### ORDERING INFORMATION:

Cat. number

Bio RS-24 BS-010133-AAG with standard platform PRS-22

Optional platforms:

PRS-4/12 BS-010117-AK PRSC-18 BS-010117-EK

Platform	Capacity	Tube Volume	Tube Diameter, Ø
1 PRS-22 (standard)	22	1.5–15 ml	10–16 mm
2 PRS-4/12 (optional)	4 and 12	up to 50 and 1.5–15 ml	20–30 mm and 10–16 mm
3 PRSC-18 (optional)	18	15 ml	16 mm

















## V-1 plus and V-32, Vortexes

**V-1 plus** vortex and **V-32** multi vortex are intended for intensive mixing of samples in tubes with an eccentric mechanism.

#### Vortex can be used for different operations:

- · Mixing tissue samples;
- · Suspending cell samples;
- · Mixing chemical samples;
- Mixing bacterial and yeast cells when washing from the culture medium;
- Extracting metabolites and enzymes from cells and cell cultures, etc.

Vortex can be used to perform various DNA/RNA operations, such as purification of low-molecular DNA/RNA fragments in PCR-diagnostics.

Vortex is applicable in all the fields of laboratory research in biotechnology, microbiology and medicine.

#### Vortexes have two operation modes:

- · Continuous operation;
- Impulse operation. (V1 plus pressure activated)


Model **V-1 plus** is a personal vortex with a fluoroplastic head for single tube vortexing.

Model **V-32** is a universal vortex multipurpose device with different accessories. It is supplied with a 32-socket universal platform PV-32 for Eppendorf type tubes up to 2 ml (2–1.5/0.5/0.2 ml–16/8/8 sockets) and a PL-1 head for vortexing a single tube up to 50 ml. An optional 6-socket platform PV-6/10 for 10 ml tubes (maximum tube diameter 15 mm) or a platform PV-48 for six strips of eight 0.2 ml microtubes can be supplied on request.



Product video

Platform PL-1 for V-32



## V-1 plus and V-32, Vortexes

	V-1 plus	V-32
Mixing principle	Vibro Eccentric	
Speed control range	500–3,000 rpm	
Acceleration time	2 s	3 s
Maximum continuous operation time	24 h	
Mixing module for tubes	from 0.2 to 50 ml	from 0.2 to 10 ml
Maximum mixing volume	30 ml	45 ml
Maximum load	30 g	70 g
Orbit	4 mm	2 mm
Dimensions (W×D×H)	90 × 150 × 80 mm	120 × 180 × 100 mm
Weight	0.8 kg	1.5 kg
Input current/power consumption	12 V, 320 mA/3.8 W	
External power supply	Input AC 100-240 V, 50/60 Hz; Output DC 12 V	




V-1 Plus



Platform PV-6/10 for V-32





Platform PV-48 for V-32

V-32 with platform PV-48



## ☐ ORDERING INFORMATION:

Cat. number

V-1 plus

BS-010203-AAG

V-32 with standard platforms PL-1 and PV-32

Optional platforms for V-32:

PV-6/10 platform for 6–10 ml tubes (max. Ø 15 mm)

BS-010207-BK

PV-48, platform for 6–8 × 0.2ml strips or 48 tubes of 0.2 ml

BS-010207-GK

## MSV-3500, Multi Speed Vortex

Product Class

MSV-3500
with platform SV-8/15

Make Spend Nurtex
MSS-3500





MSV-3500 with all platforms

MSV-3500 without platform

BS-010210-TAH

BS-010210-AAH

Multi Speed Vortex **MSV-3500** is designed for soft or intensive mixing of reagents in different size and type plastic tubes (0.2 to 50 ml).

It is designed for operation in life science laboratories working in biochemistry, cell and molecular biology.

Unit has four types of interchangeable platforms: for Eppendorf type microtest tubes, 10/15/50 ml tubes (diameter 12/16/30 mm). Platforms can be ordered separately or as one set with MSV-3500.

Speed and time are under microprocessor control. LCD display indicates two lines of values: the set and actual values of speed and time.

Unit provides high maximum speed of platform rotation efficiently mixing microvolumes (less than 5  $\mu$ l) of samples.

Speed control range	300-3,500* rpm
Digital time setting	0–60 min/non-stop (increment 1 min)
Timer sound signal	yes
Display	LCD, 2 × 16 signs
Orbit	4 mm
Maximum load	0.2 kg
Maximum continuous ope	eration time 8 h
Dimensions (W×D×H)	180 × 170 × 145 mm
Weight	2.6 kg
Input current/power cons	umption 12 V, 1 A / 12 W
External power supply	Input AC 100–240 V, 50/60 Hz, Output DC 12 V
* – Maximum speed depends on load	

Optional platforms:		Cat. number
① SV-16/8	Platform for 16 × 1.5 ml + 8 × 0.5 ml + 8 × 0.2 ml microtubes, Ø 11/8/6 mm	BS-010210-CK
2 SV-10/10	Platform for 10 × 10 ml tubes 12 mm diameter	BS-010210-BK
<b>3</b> SV-8/15	Platform for 8 × 15 ml tubes 16 mm diameter	BS-010210-DK
4 SV-4/30	Platform for 4 × 50 ml tubes 30 mm diameter	BS-010210-AK









## RCP-24, Homogenizer

Reciprocal Homogenizer RCP-24, a bench-top mechanical device designed for mixing, grinding, homogenizing and emulsifying biological objects in microtubes by vigorously mixing by reciprocal motion with various beads for sample preparation for subsequent academic, pharmaceutical, biotechnological or biomedical studies.

Homogenizer facilitates the formation of a supernatant containing nucleic acids and proteins suitable for subsequent purification, extraction or analysis. The device is optimized for extracting proteins, DNA, RNA or tRNA from various tissue sources, but it can also be used for other applications. **RCP-24** performs efficient homogenization of mammalian tissue, plant tissue or other biomaterials.

Test tubes capacity	up to 24
Test tubes	2 ml
Speed control range	500–2,000 rpm (increment 100 rpm)
Digital time setting	1–15 min (increment 1 min)
Oscillation amplitude	44 mm, vertical
Dimension (W×D×H)	285 × 400 × 440 mm
Input current / power consumption	230 V, 50 Hz/ 220 W (1.3 A)
Weight	19.1 kg











ORDERING INFORMATION:

RCP-24, Homogenizer

Cat. number

BS-010701-A02

## **NEW UIS-360, Universal Inoculation Spinner**

Speed control range

Weight







#### **Key Features and Advantages:**

- Universal Compatibility fits both square (125×125 mm) and round Petri dishes (up to Ø105 mm) out-of-the-box, and even up to Ø150 mm with an optional adapter.
- Widest Speed Range adjustable speed control from 10 to 300 RPM – offering the widest range among similar plate spinners.
- Continuous and Hands free Operation supports continuous rotation for extended procedures, with an included foot-switch for hands-free control.
- Enhanced Safety equipped with a protective polycarbonate face shield that guards against accidental splashes and aerosol droplets.
- Compact & Portable Design measures only 170×190 mm on the lab bench, weighs ~0.85 kg, and can even be powered by a 12 V DC power bank.
- Delivers smooth, uniform rotation for even spreading of samples, improving reproducibility between plates.

The **UIS-360** Universal Inoculation Spinner is a versatile laboratory device designed to streamline microbial plating and ensure consistent results. It provides smooth, uniform rotation for spreading samples across agar surfaces, improving both efficiency and reproducibility in microbiological workflows. With an adjustable speed range of 10–300 RPM – the broadest in its class – the **UIS-360** accommodates various protocols from gentle inoculum spreading to vigorous. Its compact footprint and lightweight build make it easy to handle or even take into the field (it runs on just 3.7 W, so it can be powered by a standard 12 V battery pack). The inclusion of a foot-switch for hands-free control and a protective shield for splash protection further enhance user comfort and safety.

	1 , , ,
Operation modes	Continuous operation or foot-switch activation (on/off control)
Maximum continuo	uous 168 h
Plate compatibilit	up to 125 × 125 mm; up to Ø105 mm (platform available for up to Ø150 mm, on request)
Safety features	Removable splash shield, non-slip rubber feet, automatic shutoff on foot-switch release
Dimensions (W×D×H)	$170 \times 190 \times 110$ mm (without shield) $200 \times 190 \times 240$ mm (with shield)

Input current 12 V, 310 mA / 3.7 W /power consumption

External power Input AC 100–240 V, 50/60 Hz, supply Output DC 12 V

USB-C to power socket 12V adapter **FS-1,** Foot switch





10-300 rpm (adjustable)

ORDERING INFORMATION:

BS-010177-A01

Cat. number

0.85 kg

Optional accessories:

USB-C to power socket 12V adapter

FS-1, Foot switch

**UIS-360** 

BS-000001-S27 BS-010177-AK



## MIXING DEVICES:

# THERMO-SHAKERS



TS-DW

Thermo-Shaker for Deep Well Plates

# PST-60HL, PST-60HL-4 and PST-100HL, Thermo-Shakers

**PST-60HL, PST-60HL-4** and **PST-100HL** Thermo-shakers are designed for shaking standard 96-well microtiter plates in the thermal regulation mode. Models **PST-60HL** and **PST-100HL** hold 2 plates, model **PST-60HL-4** has four plates.

A multisystem principle used in design of the Thermo-Shaker, allows operating it as three independent devices:

- · Incubator;
- · Microplate shaker;
- · Thermo-Shaker.

A distinctive feature of Biosan Plate Thermo–Shakers is effective **Two-Side Microplates Heating**, which allows achieving full correspondence of the set and actual temperature in the microplate wells.

Standard versions of Thermo-shakers provide heating up to 60°C, sufficient for carrying out ELISA tests.

Thermo-shaker **PST-100HL** with the ability to stabilize the temperature up to 100°C is specially designed for hybridisation reactions.

### Plate Thermo-Shakers provide:

- · Soft or intensive sample shaking;
- · Rotation speed regulation, stabilization and indication
- Even rotation amplitude throughout the Thermo Automatic
- Setting and indication of the required temperature on the platform
- Automatic fault diagnostics (temperature sensor, platform heating, lid heating etc.)
- With the help of the temperature calibration function, the user can calibrate the unit to compensate for differences in the thermal behaviou of plates from different manufacturers.

# Application fields:

# PST shakers can be used in various applications such as:

- Immunochemistry Enzyme-Linked Immuno Sorbent Assay (ELISA). Unique bottom and top heating, while shaking, ensures the most efficient linkage of the target, thus providing the most reliable results;
- Molecular biology Micro and Macro array applications incubation with shaking provides more efficient hybridization of target nucleic acid with on the surface of Micro and Macro chip printed probes (Specific holder is required)








# PST-60HL, PST-60HL-4 and PST-100HL, Thermo-Shakers

	PST-60HL	PST-60HL-4	PST-100HL	
Temperature setting range	+25°C	+25°C +100°C		
Temperature control range	+5°C above an	nbient +60°C	+5°C above ambient +100°C	
Temperature setting resolution		0.1°C		
Temperature stability		±0.1°C		
Temperature uniformity @ +37°C	±0.2	25°C	±0.2°C	
Temperature calibration coefficient range		0.936-1.063 (±0.063)		
Heating	Two-side microplate he	Two-side microplate heating (platform and lid) + double heating contour of the platform		
Orbit		2 mm		
Speed regulation range	250	–1,200 rpm (increment 10 r	pm)	
Digital time setting	1 min	–96 h/non-stop (increment	1 min)	
Timer sound signal		yes		
Display		LCD, 2 × 16 signs		
Max. height of microtest plate		18 mm		
Number of microtest plates	2	4	2	
Weight	6.1 kg 8.8 kg		5.9 kg	
Platform dimensions (W×D)	250 × 150 mm	290 × 210 mm	250 × 150 mm	
Overall dimensions (W×D×H)	270 × 260 × 125 mm	380 × 390 × 140 mm	270 × 260 × 125 mm	
Input current/power consumption	12 V DC, 3.3 A/40 W 12 V DC, 4.15 A/50 W 12 V, 5 A/60 W			
External power supply	Input AC 100-240 V 50/60 Hz, Output DC 12 V			

# Heat up time PST-60HL and PST-60HL-4:



# Heat up time PST-100HL:

from 25°C 60 min to 100°C

PST-60HL-4 spring holders

ORDERING INFORMATION:





Cat. number

**PST-60HL** BS-010119-AAI

**PST-60HL-4** BS-010128-AAI

**PST-100HL** BS-010142-AAI

# TS-100 and TS-100C, TS-100C Smart Thermo-Shakers

TS-100 and TS-100C, TS-100C Smart thermo-shakers are designed for intensive mixing of various reaction vessels in a temperature control environment. The TS-100C model of thermo-shaker differs from TS-100 in the possibility of cooling samples down to +4°C.

# Features of thermo-shakers meet the highest expectations of users according to many parameters:

- · Fast reaching of specified mixing speed and maintenance of equal amplitude of rotation throughout the thermo-shaker block;
- Stability of maintaining the set temperature in a wide range throughout the block surface of thermo-shakers;
- With the help of the temperature calibration function, the user can calibrate the unit approximately ±6% of the selected temperature to compensate differences in the thermal behaviour of tubes from different manufacturers;
- LCD display indicates pre-set and current values of temperature, speed and time of operation;
- · Quiet motor operation, compact size, prolonged service life. Functions of heating and mixing can be performed either simultaneously or independently, which allows using the unit as three independent devices:
- · Thermostat;
- Shaker:
- · Thermo-shaker.

We offer 13 thermostating platforms for TS-100 and 13 for TS-100C, TS-100C Smart. Within one model of thermoshaker, the blocks are mutually interchangeable and can be easily installed.

# TS-100C Smart model allows you to control the device in the following modes:

- 1. Manual using the front panel interface.
- 2. Through a computer program using Bluetooth® technology. The software allows you to manage following

# parameters:

- · Rotation speed
- Temperature
- Time
- · Sound signal
- Creating Profiling programs using controlled parameters
- · Visualization of temperature vs time and speed vs time graphs
- · Data export to Excel and CSV formats
- · Error messages/Fault diagnostics

Possibility of control up to seven units from PC. Independent parameter setting allows performing different tasks simultaneously on several units.







Ø2mm



connection

# **TS-100** and **TS-100C, TS-100C Smart** Thermo-Shakers

	TS-100	TS-100C, TS-100C Smart		
Temperature setting range	+25°C +100°C +4°C +100°C			
Temperature control range	5°C above ambient +100°C	15°C below ambient +100°C		
Temperature setting resolution	0.1	1°C		
Temperature stability @ +37°C	±0.	1°C		
Calibration option	ye	es		
Temperature calibration coefficient range	0.936-1.06	53 (±0.063)		
Speed control range	250–14	00 rpm		
Speed setting resolution	10 :	rpm		
Orbit	2 n	nm		
Display	LCD, 2 ×	16 signs		
Digital time setting	1 min-96 h (1 r	min increment)		
Timer sound signal	ye	es		
Overall dimensions with thermoblock (W×D×H)	220 × 240	× 130 mm		
Weight	3.7 kg			
Input current/power consumption	12 V, 3.5 A/42 W 12 V, 4.9 A/60 W			
External power supply	Input AC 100–240 V, 50/60 Hz; Output DC 12 V			
PC software	– only for TS-100C Smart			







ORDERING INFORMATION:	Cat. number 🗔
TS-100 without block	BS-010120-AAI
TS-100C without block	BS-010143-AAI

**TS-100C Smart** with software, without block BS-010171-A01

Photos and descriptions of all blocks can be found on next page

# Interchangeable Blocks for **TS-100**

ų.	ORI	DERING INFO	PRMATION			Cat. number
		cional cks:	Description	Maximal RPM	Maximal temperature	
	0	VP-8/5	8 × 5 ml conical tubes	1200	100	BS-010175-SK
	0	VP-4	4 × 50 ml conical tubes	1000	80	BS-010175-GK
	3	VP-8/15	8 × 15 ml conical tubes	1100	80	BS-010175-HK
NEW	4	VP-CV-20	20 × 10 mm cuvettes (0.2–4.5 ml)	1100	80	BS-010175-IK
Z	6	VP-32	32 × 0.5 ml microtubes	1400	100	BS-010175-JK
	6	VP-CL-24	24 × 3.6–4.5 ml cryotubes	1300	100	BS-010175-KK
	Ø	VP-CS-24	24 × 1–1.8 ml cryotubes	1400	100	BS-010175-LK
	8	VP-20	20 × ø12 mm round bottom tubes	1400	100	BS-010175-TK
	9	SC-18	20 × 0.5 ml +12 × 1.5 ml microtubes	1400	100	BS-010120-AK
	1	SC-18/02	20 × 0.2 ml microtubes +12 × 1.5 ml microtubes	1400	100	BS-010120-CK
	0	SC-24N	24 × 1.5 ml microtubes	1400	100	BS-010120-GK
	Ð	SC-24	24 × 2 ml microtubes	1400	100	BS-010120-EK
	ß	SC-96A	96-well unskirted or semi-skirted microplate (0.2 ml) for PCR or 12 × 8–0.2ml strips or 96 tubes of 0.2 ml	1400	100	BS-010120-FK







# Interchangeable Blocks for **TS-100C** and **TS-100C Smart**

Opt Bloc	ional cks:	Description	Maximal RPM	Maximal temperature	
0	VP-8/5C	8 × 5 ml conical tubes	1200	100	BS-010176-
0	VP-4C	4 × 50 ml conical tubes	1000	80	BS-010176-0
3	VP-8/15C	8 × 15 ml conical tubes	1100	80	BS-010176-I
4	VP-CV-20C	20 × 10 mm cuvettes (0.2–4.5 ml)	1100	80	BS-010176
0	VP-32C	32 × 0.5 ml microtubes	1400	100	BS-010176
0	VP-CL-24C	24 × 3.6–4.5 ml cryotubes	1300	100	BS-010176-
Ð	VP-CS-24C	24 × 1–1.8 ml cryotubes	1400	100	BS-010176-
8	VP-20C	$20 \times \text{ø}12 \text{ mm}$ round bottom tubes	1400	100	BS-010176-
9	SC-18C	$20 \times 0.5 \text{ ml} + 12 \times 1.5 \text{ ml microtubes}$	1400	100	BS-010143-
0	SC-18/02C	$20 \times 0.2$ ml microtubes +12 × 1.5 ml microtubes	1400	100	BS-010143-
0	SC-24NC	24 × 1.5 ml microtubes	1400	100	BS-010143-0
Ø	SC-24C	24 × 2 ml microtubes	1400	100	BS-010143-
B	SC-96AC	96-well unskirted or semi-skirted microplate (0.2 m for PCR or $12 \times 8$ –0.2ml strips or 96 tubes of 0.2 m		100	BS-010143-
0	VP-8/5C	2 VP-4C 3 VP-	-8/15C	<b>4</b> VP-CV	-20C
(	VP-8/5C		-8/15C	4 VP-CV	
5			-CS-24C		

SPECIFICATIONS

# TS-DW, Thermo-Shaker for deep well plates

**TS-DW** Thermo-Shaker is designed for shaking and incubating deep well plates.

A multisystem principle, used in the Thermo-Shaker design, allows operating it as three independent devices: Incubator, Plate shaker and Thermo-Shaker.

**TS-DW** provides excellent temperature uniformity across the plate due to effective two-sided heating of the block and the lid, contour heating of the block and close proximity of heating elements to plate walls.

There is a number of interchangeable blocks to suit different plates such as Eppendorf® 96/1,000 µl, Sarstedt® Megablock 96/2,200 µl, Porvair® 96/2,000 µl, Axygen® 96/2,200 µl. Also, we can manufacture a customized block on request.



The block for deepwell plate is mountable, thus a custom plate module can be manufactured on request

Temperature setting range

Temperature control range	5 °C above ambient +100 °C
Temperature setting resoluti	on 0.1°C
Temperature uniformity @ +	37 °C ±0.1 °C*
Temperature accuracy @ +37	7 °C ±0.5 °C*
Temperature calibration coefficient range	0.936-1.063 (±0.063)
Time of platform heating from	m +25 °C to +37 °C 6 min*
Speed control range	250-1,400 rpm
Orbit	2 mm
Display	LCD, 2×16 signs
Digital time setting	1 min-96 h (1 min increment)
Timer sound signal	yes
Overall dimensions (W×D×H)	240 × 260 × 160 mm
Weight	5.1 kg
Input current/power consum	nption 12 V, 4.8 A/58 W
External power supply	Input AC 100-240 V 50/60 Hz;

\* – For B-2E block



Cat. number

Output DC 12 V

+25 °C ... +100 °C

**TS-DW** without block BS-010159-A02





Product video

### **Deep Well Plate Thermo-Shaker provides:**

- · Soft or intensive sample shaking;
- Rotation speed regulation, stabilisation and indication;
- Even rotation amplitude throughout the Thermo-Shaker platform;
- Exceptional temperature uniformity across the plate:
- · Required operation time setting and indication;
- Automatic stopping of the platform movement after expiration of the set time;
- Setting and indication of the required temperature on the platform;
- A variety of changeable blocks that can accommodate most popular deepwell plates;
- Automatic fault diagnostics (temperature sensor, platform heating, lid heating etc.).

# Separate blocks to accommodate different deepwell plates will be released. For example:

- Deep Well Plates NUNC® 96/2,000 μl
- Deep Well Plates Eppendorf® 96/0.5 ml

### **Application fields:**

- Cytochemistry for in situ reactions;
- **Immunochemistry** for immunofermentative reactions;
  - **Biochemistry** for enzyme and protein analysis;
- Molecular biology for nucleic acid extraction.

# **TS-DW,** Thermo-Shaker for deep well plates



ORDERING INFORMATION			Cat. number
Interchangeable Blocks:	Volume	D-w plate cat. no.	

interchai	igeable blocks.	volume	D-w plate cat. no.	
B-PCR	Block for one unskirted, semi-skirted, skirted PCR microplate	96/200 μl	-	BS-010159-YK
AP-96	Block for one deep-well plate Eppendorf®	96/500 µl	0030501101	BS-010159-QK
B-05PO	Block for one deep-well plate PlateOne®	96/0.5 ml	S1896-5000	BS-010159-OK
B-06A	Block for one deep-well plate Corning Axygen®	96/600 μΙ	P-DW-500C	BS-010159-KK
B-08AB	Block for one deep-well Abgene™ Storage Plate	96/800 µl	AB0765, AB0859	BS-010159-MK
B-12AB	Block for one deep-well plate Abgene™	96/1.2 ml	AB0564, AB0787	BS-010159-SK
B-1R	Block for one deep-well plate Riplate®	96/1,000 µl	43001-0101	BS-010159-UK
B-2A	Block for one deep-well plate Axygen®	96/2,200 µl	P-2ML-SQ-C	BS-010159-FK
B-2BBI	Block for one deep-well strips with rack SSIbio®	96/1.2 ml	703B00, 713B00	BS-010159-JK
B-2KF	Block for one PowerMag® Glass Bead 96 Plate	96/-	27600-4-KF-BP	BS-010159-LK
B-2M	Block for one KingFisher™ deep-well 96 Plate	96/-	95040450	BS-010159-RK
B-2N	Block for one deep-well plate Nunc®	96/2,000 µl	278743, 278752	BS-010159-DK
B-2PO	Block for one deep-well plate PlateOne®	96/2 ml	S1896-2000	BS-010159-NK
B-2R	Block for one deep-well plate Riplate®	96/2,000 μΙ	43001-0103	BS-010159-TK
B-2S	Block for one deep-well plate Sarstedt® Megablock	96/2,200 µl	82.1972.002	BS-010159-CK
B-2SD	Block for one deep-well plate Slicprep™ 96 Device	96/2.2 ml	AB0932	BS-010159-PK
B-2SL	Block for one deep-well plate Starlab®	96/1,200 µl	E2896-0120	BS-010159-IK
B-2E	Block for one deep-well plate Eppendorf®	96/1,000 µl	0030505204, 0030506200, 0030502205	BS-010159-AK
B-2P	Block for one deep-well plate Porvair <sup>®</sup> , Nunc <sup>®</sup> , Biotage <sup>®</sup>	96/2,000 μl	219009, 95040452, 121-5203	BS-010159-EK



# MINICENTRIFUGES VORTEXES, MINI CENTRIFUGE, CENTRIFUGES



**CVP-2**Centrifuge/Vortex for PCR plates

# FV-2400, Microspin and FVL-2400N, Combi-Spin

Minicentrifuges-Vortexes Microspin **FV-2400** and Combi-Spin **FVL-2400N** is specially designed for genetic engineering research (for PCR-diagnostics experiments). Units can be used in biomedical and biotechnological laboratories.

Minicentrifuges-Vortexes provide simultaneous mixing and separation of 12 samples, using centrifuge and mixing modules located on the common spin-module. Sequential combination of these operations allows you to collect all material at the bottom of the tube.

**FV-2400** is an "open type" centrifuge (without lid), that increases the speed of centrifugation and resuspension operations.

**FVL-2400N** has a bioform design and equipped with a transparent protective lid accompanied by a protection mechanism that stops the rotor motion when the lid is opened.













	FV-2400	FVL-2400N	FV-2400	FVL-2400N
Rotation speed (fixed)	2,800	) rpm	3,500 rpm	
Max. RCF	500×g		700×g	
Continuous and impulse operation	n modes			
Safety		Stop at open lid		Stop at open lid
Overall dimensions (W×D×H)	120 × 170 × 120 mm	190 × 235 × 125 mm	120 × 170 × 120 mm	190 × 235 × 125 mm
Weight	1.4 kg	1.7 kg	1.4 kg	1.7 kg
Nominal operating voltage	230 V, 50 Hz	230 V, 50 Hz	120 V, 60 Hz	120 V, 60 Hz
Power consumption (120 / 230 V)	30 W (	0.13 A)	30 W (	0.27 A)

SPECIFICATIONS

# Rotors for FV-2400 and FVL-2400N

## ORDERING INFORMATION:

Cat. number

FV-2400 white with standard rotor R-1.5M and R-0.5/0.2M

BS-010201-AAA

FVL-2400N with standard rotors R-1.5 and R-0.5/0.2

BS-010202-AAA

Optional rotors: see table below

Rotors for FV-2400	):	Capacity	Туре	Cat. number
1 R-0.5/0.2M	12 × 0.5 ml and 12 × 0.2 ml microtubes	24	Standard	BS-010201-BK
2 R-1.5M	12 × 1.5/2 ml microtubes	12	Standard	BS-010201-AK
3 R-2/0.5	8 × 1.5/2 ml and 8 × 0.5 ml microtubes	16	Optional	BS-010205-CK
4 R-2/0.5/0.2	6 × 1.5/2 ml. 6 × 0.5 ml and 6 × 0.2 ml microtubes	18	Optional	BS-010205-DK
5 SR-16	Two 8-section strips for 0.2 ml microtubes	16	Optional	BS-010202-AK
6 SR-64*	Eight 8-section strips for 0.2 ml microtubes	64	Optional	BS-010201-EK

\* – For any type of strips including paired





**3** R-2/0.5











**4** R-2/0.5/0.2







Rotators for FVL-24	100N:	Capacity	Type	Cat. number
1 R-0.5/0.2	12 × 0.5 ml and 12 × 0.2 ml microtubes	24	Standard	BS-010205-BK
2 R-1.5	12 × 1.5/2 ml microtubes	12	Standard	BS-010205-AK
3 R-2/0.5	8 × 1.5/2 ml and 8 × 0.5 ml microtubes	16	Optional	BS-010205-CK
4 R-2/0.5/0.2	6 × 1.5/2 ml. 6 × 0.5 ml and 6 × 0.2 ml microtubes	18	Optional	BS-010205-DK
<b>5</b> SR-16	Two 8-section strips for 0.2 ml microtubes	16	Optional	BS-010202-AK
<b>6</b> SR-32*	Four 8-section strips for 0.2 ml microtubes	32	Optional	BS-010205-FK

\* – Not compatible with Combi-Spins produced before 2015















# MSC-3000 and MSC-6000, Multi-Spins

Centrifuge/vortex Multi-Spins MSC-3000 and MSC-6000 are products of extensively evolving Spin-mix-Spin technology that is intended for collecting micro volumes of reagents on the microtube's bottom (first centrifugation spin), following mixing (mix) and collecting the reagents again from the walls and cap of the microtube (second spin). Aim of this repetitive algorithm of operation is to reduce the mistakes during sample preparation for PCR analysis. We named it "sms-algorithm".

Multi–Spin is a fully automatic device for reproducing sms–algorithm for 12 tubes at one time, thus saving time considerably. A must-have instrument for PCR and DNA analyses laboratory.

### Multi Spin is four devices combined in one:

1. Centrifuge - Maximum RCF:

**MSC-3000**: up to  $800 \times g$  **MSC-6000**: up to  $2,350 \times g$ 

- 2. Vortex (3 mixing modes soft, medium, hard; regulated time; Vortexing regulation timer 1–20 s)
- 3. Centrifuge/Vortex;
- **4.** SMS-cycler for realisation of the "sms-algorithm".







Saving time with multi-spin Multi-Spin allows considerable time saving compared to Combi-Spin by automatically performing a cycling program of sample mixing and spinning according to the set spin-mix-spin cycle for 12 microtubes			
simultaneously.	FVL-2400N	MSC-3000	MSC-6000
Speed control max.	2,800 rpm	3,500 rpm	6,000 rpm
RCF max.	500×g	800×g	2,350×g
Number of tubes vortexing	1 individually	12 simult	aneously
Time for completing "spin-mix-spin" cycle:			
for 2 microtubes	60 s	25 s	15 s
for 12 microtubes	5–6 min	90 s	60 s
for 100 microtubes	60 min	15 min	10 min
Unit price ratio	1 ×	1.5 ×	1.6 ×

# MSC-3000 and MSC-6000, Multi-Spins

	MSC-3000	MSC-6000	
Speed regulation range (increment 100 rpm)	1,000–3,500 rpm	1,000-6,000 rpm	
RCF max.	800×g	2,350×g	
Spin timer	1 s-99 min	1 s-30 min	
Timer sound signal	yes		
Vortexing intensity	Soft, medium, hard		
Vortexing time	0–20 s (increment 1 s)		
SMS-cycle regulation	1–999 cycles		
Display	LCD, 2×16 signs		
Safety	Autostop at open lid	Lid lock	
Overall dimensions (W×D×H)	190 × 235 × 125 mm		
Weight	2.1 kg	2.5 kg	
Input current/power consumption	12 V, 11 W (0.9 A )	24 V, 24 W (1 A)	
External power supply	Input AC 100–240 V 50/60 Hz; Output DC 12 V	Input AC 100–240 V 50/60 Hz; Output DC 24 V	

ORDERING INFORMATION:

Cat. number

MSC-3000 with standard rotors R-1.5, R-0.5/0.2

BS-010205-AAN

MSC-6000 with standard rotors R-1.5, R-0.5/0.2

BS-010211-AAL

### Optional rotors: see table below

Rotor R-1.5



MSC-3000




Optional rotors:		Capacity	Туре	Cat. number
1 R-0.5/0.2	12 × 0.5 ml and 12 × 0.2 ml microtubes	24	Standard	BS-010205-BK
2 R-1.5	12×1.5/2 ml microtubes	12	Standard	BS-010205-AK
3 R-2/0.5	$8 \times 1.5/2$ ml and $8 \times 0.5$ ml microtubes	16	Optional	BS-010205-CK
4 R-2/0.5/0.2	$6 \times 1.5/2$ ml, $6 \times 0.5$ ml and $6 \times 0.2$ ml microtubes	18	Optional	BS-010205-DK
<b>5</b> SR-16	Two 8-section strips for 0.2 ml microtubes	16	Optional	BS-010202-AK
6 SR-32*	Four 8-section strips for 0.2 ml microtubes	32	Optional	BS-010205-FK

\* - Not compatible with Multi-Spins produced before 2015









4 R-2/0.5/0.2



**5** SR-16



**6** SR-32











# CVP-2, Centrifuge vortex for PCR plates

After many years of Combined Centrifuge/Vortex concept success, we are proud to introduce the long-awaited Centrifuge vortex for PCR plates, **CVP-2**, to the sample preparation market.

The Spin-Mix-Spin technology is intended to spin-down micro volumes of reagents on the well's bottom (first centrifugation spin), following mixing (mix) and spin-down the reagents again from the walls and cap of the well (second spin). Aim of this repetitive algorithm of operation is to reduce the mistakes during sample preparation for PCR analysis.

**CVP-2** is a fully automatic device for reproducing sms-algorithm for 2 PCR plates at the same time, thus saving time considerably. A must-have instrument for PCR and DNA analyses laboratory.

### CVP-2 is 4 devices combined in 1:

- 1. Centrifuge Maximum RCF: 225 × g (1,500 rpm)
- 2. Vortex (300–1,200 rpm; Vortexing regulation timer 0–60 sec)
- 3. Centrifuge vortex

Speed regulation range

4. SMS-cycler for realization of the "sms-algorithm"

### Tested plate types for use with CVP-2 centrifuge:

- Full-skirted 96-well standard micro-plates (without adapter)
- Half-skirted 96-well standard micro-plates (with adapter AP-96)
- Unskirted 96-wel standard I micro-plates (with adapter AP-96)
- Applied Biosystems<sup>™</sup> MicroAmp<sup>™</sup> Optical 96-well reaction plate (with adapter AP-96)
- Applied Biosystems<sup>™</sup> MicroAmp<sup>™</sup> Optical 384-well reaction plate (with adapter AP-384)

300-1,500 rpm

Output DC 12 V

• For specific plate usage, please contact us for evaluation.

opeca regulation range	500 1/500 p
Min. RCF at 1,500 rpm	185 × g
Vortex regulation range	300-1,200 rpm
Setting resolution	100 rpm
Plate type:	
Without adapter:	96-well skirted PCR plates,
	PCR strips in a frame;
• With adapter <b>AP-96</b> :	96-well semi-skirted
	and non-skirted PCR plates;
• With adapter AP-384:	384-well PCR plates;
Display	LCD, 2×16 signs
Centrifugation mode time range	0-30 min
Centrifugation mode time increment	1 s; after 1 min–1 min
Vortex mode time range	0-60 s
Timer sound signal	yes
Number of programmable cycles	1-999
Chamber diameter	210 mm
Overall dimensions (W×D×H)	285 × 350 × 190 mm
Weight	6.15 kg
Input current/power consumption	12V,1.5 A/18 W
External power supply	Input AC 100-240 V 50/60 Hz;





### ORDERING INFORMATION

Cat. number



# CVP-2

BS-010219-A02

With rotor for two PCR plates, protection lid and adapters AP-96\* (a set of 2 adapters for 96-well semi-skirted and unskirted PCR plates)

### Optional accessories:

### AP-384\*

BS-010219-EK

A set of 2 adapters for 384-well PCR plates

<sup>\* -</sup> Adapters are made of Ertacetal® C and are autoclavable

# Microspin 12 Plus, High-speed Mini-centrifuge



# SPECIFICATIONS

- SECTICATIONS	
Rotor imbalance automatic diagnostics	emergency stop, IMBALANCE indication
Speed control range	1000 – 16,250 rpm (increment 100 rpm)
Relative centrifugal force control	range 53 – 15,588 <i>g</i>
Digital time setting	15 s – 90 min
Time setting resolution	15 s (<1 min) 1 min (>1 min)
Quick spin	+
Acceleration time to 16,250 rpm, accurate within ±8s	20 s
Slowdown time to 0 rpm, not mo	re 10 s
Display	LCD
Rotors available	2
Max. capacity (MSR-12)	12 × 1.5/2 ml tubes
Max. capacity (MSR-16) 1	6 × 0,2ml tubes/8-well strips
Tilt angle (MSR-12)	30°
Tilt angle (MSR-16)	45°
Maximum tube height (MSR-12)	≤ 47 mm
Maximum tube height – above the	e rotor (MSR-12) ≤ 8 mm
Maximum rotor load	36 g
Dimensions	200 × 240 × 125 mm
Weight, accurate within ±10%	3.5 kg
Input current/ power consumption	24 V, 2.5 A / 60 W
External power supply	Input AC 100–240 V 50/60 Hz; Output DC 24 V

The Microspin 12 Plus is a high-quality minicentrifuge designed for efficient separation of various components in a range of applications, including RNA/DNA extraction, cell suspension separation, and other micro quantitative analyses. It's sleek, bioform design and compact footprint make it a space-saving solution to the laboratory workspace. The centrifuge can accommodate microtubes or strip tubes, reaching a maximum acceleration of 16,250 RPM or 15,588

The centrifuge comes with a selection of two durable aluminum fixed angle rotors:

- 1. Microtube rotor MSR-12 has the capacity to hold up to 12 × 2 ml vessels (microtubes, spin columns), with adapters provided for 20.5 ml and 10.2 ml microtubes. The lid enables centrifugation as mictotubes as an spin columns.
- 2. Optional strip rotor MSR-16 (the lid MSL-16 is included) can hold up to two 8 × 0.2 ml PCR strip rows or 16 individual 0.2 ml PCR tubes.

Equipped with microprocessor control, the centrifuge ensures precise control of actual and set parameters, offering a user-friendly interface with simple setup. Brushless motor provides quiet, vibration-free performance even at high speeds, ensuring a long product service life.

Safety features include metal protective inserts and enclosures inside the body and lid of the centrifuge, an automatic imbalance detection system with an automatic stop function, and a lid locking mechanism, providing secure operation throughout the speed range. A sound signal alerts when centrifugation is complete.

## 1 A-02 Adapters





### **ORDERING INFORMATION:**

Cat. number BS-010218-AA1



# Microspin 12 Plus

Built-in rotor MSR-12 (12 places for microtubes 1.5/2 ml) with protection lid MSL-SC and adapters A-02, A-05 (autoclavable)

### Additional/replacement parts:

MSL-SC, protection lid for rotors BS-010213-EK

- A-02, 12 pieces for microtubes 0.2 ml BS-010213-BK
- **2** A-05, 12 pieces for microtubes 0.5 ml BS-010213-AK MSR-16, rotor for two 8 x 0.2 ml BS-010218-AK PCR strip rows or 16 individual 0.2 ml PCR tubes.

# MSC-2P, Minicentrifuge-Vortex for PCR plates NEW

**MSC-2P** is a compact sized digital centrifuge intended to collect droplets, mix reagents and collect once more for improved PCR yield in subsequent analysis. The combination of spin-mix functions ensures fast operation, thorough mixing and repeatable results. Centrifuge rotor can accommodate 2 PCR plates at the same time, thus saving time considerably.

# MSC-2P is possible to operate in 4 independent modes:

- · Centrifuge Max. 3,500 rpm
- Vortex up to 5 min
- · Centrifuge/Vortex combined two motion types
- · Spin-mix-spin algorithm up to 10 cycles

The spin-mix-spin algorithm (SMS-algorithm) is designed to collect (or reset) micro volumes of reagents to the bottom of the PCR plate tubes (the first centrifugation or spin), then vortexing (mix) and re-collecting reagents (repeated spin) from the walls and cover. This repetitive algorithm of operations, aimed at reducing sample preparation errors, we call the SMS algorithm.

Capacity	2 × 96 well PCR microplates, semi-skirted or non-skirted, PCR tube strips
Centrifugation mode	500-3,500 rpm
Relative centrifugal force	Up to 610 g
Centrifugation timer	1 s – 10 min
Vortexing mode	Full rotor
Vortexing timer	0–5 min
SMS-cycle regulation	1–10 cycles
Protection	Lid lock
Display	LCD, 2×16 symbols
Dimensions, (W×D×H)	165×220×230 mm
Weight	2.7 kg
Input current/power consu	mption 12 V, 1.5 A/18W
External power supply	Input AC 100–240 V 50/60 Hz; Output DC 12V










ORDERING INFORMATION:


Cat. number

**MSC-2P** BS-010204-A02

Optional accessories:

# MSC-2P, Minicentrifuge-Vortex for PCR plates





# LMC-3000, Laboratory Centrifuge

LMC-3000 is a modern low-speed bench-top centrifuge designed for operation with microtest plates and centrifuge tubes up to 50 ml, Gel Cards. This device is widely used in biomedical profile laboratories.





### **Features:**

- · Soft start and run-down of the rotor;
- · User-friendly setting of centrifugation parameters and simultaneous display of both set and actual values;
- · Safe operation at any speed is provided by metal protection chamber and case cover, automatic stop at imbalance and a lock keeping the lid closed while the centrifuge is running;
- Low noise level:
- Rotor selection;
- Setting rotor speed in RPM or RCF (Relative Centrifugal

Force);	
<ul> <li>Multiple accelerations (Slow, Normal, Fadeceleration (0, Slow, Normal, Fast) moto switch off forced braking;</li> </ul>	
Wide choice of accessory rotors.	
Speed regulation range for centrifuge tubes	100–3,000 rpm (1,610 × g)
Speed regulation range for microtitre plates	100–2,000 rpm (560 × g)
Setting resolution	100 rpm
Rotor imbalance diagnostics (automatic stop, "IMBALANCE" warning)	

(automatic stop, INIBALANCE	warning)
Display	LCD, 2×16 signs
Digital time setting	1-90 min (increment 1 min)
Timer sound signal	yes
Chamber diameter	340 mm
Overall dimensions (W×D×H)	425 × 495 × 235 mm
Weight	11.8 kg
Nominal operating voltage	230 V, 50/60 Hz or 120 V, 50/60 Hz
Power consumption (230/120	V) 110 W (0.5 A)/ 120 W (1 A)









ORDERING INFORMATION:

Cat. number

LMC-3000 without rotors

BS-010208-AAA

mnorature central range

# LMC-4200R, Laboratory Refrigerated Centrifuge







### **Features:**

- Effective way of acceleration and deceleration: Run-up time 20 sec; Run-down time, not more 30 sec;
- Efficient rate of chamber refrigeration: under 10 min;
- Maintenance of stable temperature during operation;
- User-friendly setting of centrifugation parameters (speed, temperature, time) and simultaneous display of both set and actual values;
- Safe operation is provided by a metal protection chamber and a case cover, automatic stop at imbalance (emergency shutdown, "IMBALANCE" displayed) and a lock keeping the lid closed while the centrifuge is running;
- · Low noise level;
- · Possibility to switch off forced braking;
- · Wide choice of accessory rotors;
- · Rotor selection;
- Setting rotor speed in RPM or RCF (Relative Centrifugal Force);
- Multiple accelerations (Slow, Normal, Fast) and deceleration (0, Slow, Normal, Fast) modes and possibility to switch off forced braking;

Laboratory bench-top centrifuge with refrigeration LMC-4200R provides temperature control of biomaterial during centrifugation. Temperature control of the so-called "cold-shelf" is a gold standard for enzymologists and cell biologists because it ensures conditions necessary for reproducibility of the sample preparation stage. Temperature control absence at this stage can cause unpredictable results.

**LMC-4200R** is a modern centrifuge designed for operation with microtest plates, Gel Cards and tubes from 2 to 50 ml.

Temperature control range	-10°C +25°C
Stable temperature maintenance range	25°C below ambient to +25°C
Temperature setting resolu	tion 1°C
Speed regulation range for centrifuge tubes	100-4,200 rpm (3,160×g)
Speed regulation range for microtitre plates	100-2,000 rpm (560 × g)
Speed setting resolution	100 rpm
Rotor imbalance diagnostic (automatic stop, "IMBALANG	
Slowdown time, not more	30 s
Display	LCD, 2 lines
Digital time setting	1-90 min (increment 1 min)
Timer sound signal	yes
Chamber diameter	360 mm
Dimensions ( W×D×H)	635 × 580 × 335 mm
Weight	56 kg
Nominal operating voltage	230 V, 50 Hz
Power consumption (230 V)	990 W (4.3 A)

Rotor R-24/10



**ORDERING INFORMATION:** 

Cat. number

LMC-4200R without rotors

BS-010212-AAA

Rotors description and pictures can be found on next pages



# HOW TO CHOSE ROTOR

# Interchangeable Rotors for LMC-3000 and LMC-4200R











Only for LMC-4200R

		Rotor <b>R-12/10</b>	Rotor <b>R-24/10</b>	Rotor <b>R-6</b>	Rotor <b>R-6P</b>
Rotor type		Swing-out			
Dimensions		Ø16 × 105 mm Ø29 × 115 mm			
Capacity		12	24	6	
Tube's volum	e	10–15 ml 50 ml		ml	
		4,200 rpm 4,200 rpm 4,200 rpm		) rpm	
Max. RCF:	LMC-3000	1,610 × g	Not applicable	1,61	0 × g
IVIAX. RCF.	LMC-4200R	3,160 × g	2,860 × g	3,160 × g	
Cat. number:		BS-010208-BK	BS-010212-JK	BS-010208-DK	BS-010208-XK

# Plastic conical bottom centrifuge tube

Manufacturers: Falcon, Greiner Bio-one, Sarstead, Corning, Nunc, TPP, etc.





# HOW TO CHOSE ROTOR

# Interchangeable Rotors for LMC-3000 and LMC-4200R









Rotor <b>R-2</b>	Rotor <b>R-24GC</b>	Rotor <b>R-12/15</b>	Rotor RMT-24
Swing	g-out	Angled Swing-out	Fixed-angle
128 × 85.6 mm	53 × 74 mm	Ø17 × 120 mm	Ø11.1
2	24	12	24
up to 45 mm	-	15 ml	1.5-2 ml
2,000 rpm	1,500 rpm	4,200 rpm	4,200 rpm
560 × g	280 × g	1,610 × g	Not applicable
560 × g	280 × g	3,160 × g	1,837 × g
BS-010208-AK	BS-010208-VK	BS-010208-EK	BS-010221-BK

Standard 96-well microtitre plates, skirted PCR plates and deepwell plates up to 45 mm

Manufacturers: Nunc, Greiner, Greiner Bio-one, etc.




2



96-well semi-/ unskirted PCR plate

Manufacturers: Nunc, Greiner, Greiner Bio-one, etc.









R-2

Material: Ertacetal® POM-C and is autoclavable

### **Gel Cards**

Manufacturers: Grifols®, DiaMed®, Bio-Rad® etc.



**R-24GC,** Rotor for Gel Cards for blood group serology testing (Forward Group, Reverse Group, RhD Type and 3 cell antibody screen). Recommended centrifugation time – 9 minutes



R-24GC

Cat. number

### ORDERING INFORMATION: optional accesories for rotors

Adapter\* for R-2:

AP-96 2 adapters for 96-well semi-skirted and non-skirted PCR plates BS-010219-DK

Adapters\*\* for R-12/10, R-24/10: Vacutainers dimensions (Ø × length)

Adapters** for R-12/10, R-24/10:		vacutainers dimensions (Ø×iength)	
BN-13/75	for vacutainers® 2–5 ml	13 × 80 mm	BS-010208-PK
BN-13/100	for vacutainers® 4–8 ml	13 × 105 mm	BS-010208-QK
BN-16/100	for vacutainers® 8–10 ml	16 × 105 mm	BS-010208-RK

Rack for rotors

**RR-U** BS-010208-UK

<sup>\* –</sup> Set of 2 adapters, made of **Ertacetal® POM-C and** is autoclavable

<sup>\*\* -</sup> Set of 12 adapters, made of **Kocetal® POM.** Max. temperature +100°C

100 6 000 rpm

Speed regulation range

# LMC-56, Laboratory Centrifuge

**LMC-56** is a modern benchtop, low-speed laboratory centrifuge designed for convenient sedimentation, centrifugation and collection of necessary samples. It provides operation with tubes, blood collection systems, gel cards, microtest and ELISA plates. Our centrifuges are designed for safe work (metal protecting housing), easy maintenance and wide application range in medical, biochemical, chemical, industrial and other type of laboratories.

for centrifuge tubes	100–6,000 rpm (3,750 × g)		
Speed regulation range for microtitre plates	100–2,000 rpm (560×g)		
Setting resolution	100 rpm		
Digital time setting	1-90 min (inc. 1 min)		
Timer sound signal	yes		
Rotor imbalance diagnostics (automatic stop, "IMBALANCE" warning)			
Rotator detection	LCD, 2×16 signs		
Display	LCD, 2×16 signs		
Chamber diameter	340 mm		
Overall dimensions (W×D×H)	560 × 480 × 315 mm		
Weight	16.5 kg		
Nominal operating voltage	230 V, 50/60 Hz or 120 V, 50/60 Hz		



Power consumption (230/120 V)





Cat. number

320 W / 330 W

LMC-56 without rotors

BS-010221-A01




### Features:

- User-friendly centrifugation parameter input and simultaneous display of the set and actual parameter values.
- Safe assay performance: metal protective housing and metal lid, automatic imbalance switch-off, lid lock during the centrifuge operation provide safe operation at all speeds.
- Rotor imbalance automatic diagnostics (emergency stop, imbalance indication).
- Automatic rotor detection with active rotational speed limit.
- Relatively high speed maximum 6,000 rpm or 3,750 g.
- Wide choice of accessory rotors and adapters.
- A wide selection of rotors up to 13.
- Improved chamber to reduce sample heating during centrifugation.
- Different modes of acceleration and deceleration, including deceleration mode with switched off forced braking.
- Possibility to set the speed both in revolutions per minute and by relative centrifugal force.



# Interchangeable bucket adapters for **LMC-56**



\* - Included in BA-4/VT

### ORDERING INFORMATION:

Rotor	Description	Capacity	Cat. number
BR-4U	Bucket rotor without bucket adapters	4	BS-010221-AK
Bucket adapters	Description	Capacity	Cat. number
<b>1</b> BA-14/2U	Bucket adapters for 14 microtest tubes of 2 ml (set of 4 pcs.)	14 × 4	BS-010221-JK
2 BA-8/5	Bucket adapters for 8 tubes of 5 ml (set of 4 pcs.)	8 × 4	BS-010221-FK
<b>3</b> BA-4/VT	Bucket adapters for 4 vacutainer adapters (set of 4 pcs.). Including 8-9 ml [16×100] adapters, BN-16/100-16	4 × 4	BS-010221-4K
4 BA-4/RT	Bucket adapters for 4 tubes of 10 ml (set of 4 pcs.)	4 × 4	BS-010221-3K
<b>⑤</b> BA-4/FT	Bucket adapters for 4 tubes of 15 ml (set of 4 pcs.)	4 × 4	BS-010221-YK
<b>6</b> BA-2/50	Bucket adapters for 2 tubes of 50 ml (set of 4 pcs.)	2 × 4	BS-010221-CK

Adapters for <b>BA-4/VT</b>	Description	Cat. number
3.1 BN-13/75-16	For bucket adapters BA-4/VT optional adapter set (16 pcs) for vacutainers 2-5 ml	BS-010221-OK
3.2 BN-13/100-16	For bucket adapters BA-4/VT optional adapter set (16 pcs) for vacutainers 4-8 ml	BS-010221-PK
3.3 BN-16/100-16	For bucket adapters BA-4/VT optional adapter set (16 pcs) for vacutainers 8-9 ml	BS-010221-QK

Adapters for <b>BA-2/50</b>	Description	Cat. number
<b>6.1</b> BI-25-8	For bucket adapters BA-2/50 optional adapter set (8 pcs) for 25 ml conical tubes	BS-010221-UK

# HOW TO CHOSE ROTOR

# Interchangeable Rotors for LMC-56










	Rotor <b>Ri-12/10</b>	Rotor <b>Ri-24/10</b>	Rotor <b>Ri-6</b>	Rotor Ri-6P
Rotor type	Swing		g-out	
Dimensions	Ø16 × 105 mm Ø29 × 115 mm		15 mm	
Capacity	12	24	6	
Tube's volume	10–15 ml 50 ml		ml	
Max. speed	4,200 rpm	4,000 rpm	4,200 rpm	
Max. RCF:	3,160 × g	2,860 × g	3,160 × <i>g</i>	
Cat. number:	BS-010221-MK	BS-010221-LK	BS-010221-HK	BS-010221-IK


## Plastic conical bottom centrifuge tube

Manufacturers: Falcon, Greiner Bio-one, Sarstead, Corning, Nunc, TPP, etc.



# Plastic round bottom centrifuge tube, Vacutainers

Manufacturers: Nunc, Greiner, Greiner Bio-one, TPP, etc.



# ■ HOW TO CHOSE ROTOR

# Interchangeable Rotors for **LMC-56**












Rotor Ri-2	Rotor Ri-24GC	Rotor <b>Ri-12/15</b>	Rotor RMT-24	Bucket rotor <b>BR-4U</b>
Swin	Swing-out		Fixed-angle	Swing-out
128 × 85.6 mm	53 × 74 mm	Ø17 × 120 mm	Ø11.1	-
2	24	12	24	4 buckets
up to 45 mm	-	15 ml	1.5-2 ml	2-50 ml
2,000 rpm	1,500 rpm	4,200 rpm	6,000 rpm	4,200 rpm
560 × g	280 × g	3,160 × g	3,750 × g	2070-3260 × g
BS-010221-GK	BS-010221-NK	BS-010221-KK	BS-010221-BK	BS-010221-AK

Standard 96-well microtitre plates, skirted PCR plates and deepwell plates up to 45  $\,\mathrm{mm}$ 

Manufacturers: Nunc, Greiner, Greiner Bio-one, etc.






Ri-2

### 96-well semi-/ unskirted PCR plate

Manufacturers: Nunc, Greiner, Greiner Bio-one, etc.







Ri-2

Material: Ertacetal® POM-C and is autoclavable

# **Gel Cards**

Manufacturers: Grifols®, DiaMed®, Bio-Rad® etc.




**Ri-24GC**, Rotor for Gel Cards for blood group serology testing (Forward Group, Reverse Group, RhD Type and 3 cell antibody screen). Recommended centrifugation time – 9 minutes



Ri-24GC





# Interchangeable Rotors for **LMC-56**

# ORDERING INFORMATION:

Adapters	Description	Cat. number
BI-25-6	For rotors R-6, R-6P, <b>Ri-6, Ri-6P</b> optional insert set (6 pcs) for 25 ml conical tubes	BS-010221-VK
BN-13/75	For rotors R-12/10, <b>Ri-12/10</b> optional adapter set (12 pcs) for vacutainers 2-5 ml	BS-010208-PK
BN-13/100	For rotors R-12/10, <b>Ri-12/10</b> optional adapter set (12 pcs) for vacutainers 4-8 ml	BS-010208-QK
BN-16/100	For rotors R-12/10, <b>Ri-12/10</b> optional adapter set (12 pcs) for vacutainers 8-9 ml	BS-010208-RK
BN-13/75-24	For rotors R-24/10, <b>Ri-24/10</b> optional adapter set (24 pcs) for vacutainers 2-5 ml	BS-010221-RK
BN-13/100-24	For rotors R-24/10, <b>Ri-24/10</b> optional adapter set (24 pcs) for vacutainers 4-8 ml	BS-010221-SK
BN-16/100-24	For rotors R-24/10, <b>Ri-24/10</b> optional adapter set (24 pcs) for vacutainers 8-9 ml	BS-010221-TK
AP-96	2 adapters for 96-well semiskirted and non-skirted PCR plates - made of Ertacetal® C. Autoclavable.	BS-010219-DK
AP-384	2 adapters for 384-well plates made of Ertacetal® C.Autoclavable.	BS-010219-EK
RR-U	Rack for rotors	BS-010208-UK

### Rotor RMT-24



Bucket rotor BR-4U with BA-2/50 adapter



# THERMOSTATED EQUIPMENT:

# THERMOSTATS – DRY BLOCK, **HEATING/COOLING SYSTEMS**



Dry block thermostat

# Bio TDB-100 and TDB-120, Dry Block Thermostats

**Bio TDB-100 / TDB-120 –** compact, easy-to-use thermostat for Eppendorf type micro tubes. It is specially designed for long incubation at different temperatures. The thermostat has an undeniable advantage in working with microquantities of reagents in microtubes. The thermostat possesses unprecedentedly high precision and uniformity of temperature distribution over the block.

With the help of the software-enabled temperature calibration function, the user can calibrate the unit in the range of several percent of the selected temperature to compensate for differences in the thermal behaviour of tubes from different manufacturers.







### Heat up times for **Bio TDB-100**:



# **BLOCKS (BUILT IN) SPECIFICATIONS:**

### Bio TDB-100

TDB-120

2 Block A-53 21 × 0.5 ml + 32 × 1.5 ml

microtubes

3 Block A-103 21 × 0.5 ml +32 × 1.5 ml + 52 × 0.2 ml

microtubes

Block for Bio TDB-100









### Heat up times for TDB-120:



# Bio TDB-100 and TDB-120, Dry Block Thermostats

	Bio TDB-100	TDB-120	
Temperature setting range	+25°C +100°C	+25°C +120°C	
Temperature control range	5°C above ambient +100°C	5°C above ambient +120°C	
Temperature setting resolution	0.1	°C	
Temperature stability	±0.	1°C	
Temperature uniformity @ +37°C	±0.	1°C	
Temperature calibration coefficient range	0.936-1.063 (± 0.063)	0.968-1.031 (± 0.031)	
Digital time setting	1 min – 96 h/non-stop (increment 1 min)		
Timer sound signal	yes		
Display	LCD, 2×16 signs		
Block capacity	24 × 2/1.5 ml +15 × 0.5 ml +	<b>A-53</b> 21 × 0.5 ml +32 × 1.5 ml	
	10 × 0.2 ml microtubes	microtubes	
		<b>A-103</b> 21 × 0.5 ml +32 × 1.5 ml	
		+ 52 × 0.2 ml microtubes	
Overall dimensions (W×D×H)	Overall dimensions (W×D×H) 210 × 230 × 115 mm 230 × 21		
Weight	2.8 kg		
Nominal operating voltage	230 V, 50/60 Hz or 120 V, 50/60 Hz		
Power consumption	200 W (870 mA)		

### **ORDERING INFORMATION:**

Cat. number

Bio TDB-100 with built-in block

BS-010412-AAA

TDB-120 with built-in block A-103

BS-010401-QAA

TDB-120 with built-in block A-53

BS-010401-PAA















DESCRIPTION

# CH-100, Heating/Cooling Dry Block

**CH-100** is the result of combining two popular Biosan instruments:

- 1. Heating Dry block and
- 2. Cooling Dry block thermostat

The combined construction of the aluminium block and Peltier element module cooled with the forced ventilation radiator provides fast switching of the cooling and heating modes.

**CH-100** is a very effective instrument for sample preparation during enzyme reactions, hybridization reactions, DNA analysis.

Microprocessor controlled time and temperature. Simultaneous indication of the and actual temperature and time.

Temperature setting range	−10 °C +100 °C
Temperature control range	30°C below ambient+100°C
Temperature setting resolution	n 0.1°C
Temperature stability	±0.1°C
Temperature uniformity @ +3	7 °C ±0.25°C
Temperature calibration coefficient range	0.936-1.063 (±0.063)
Digital time setting	1 min–96 h/non-stop (increment 1 min)
Timer sound signal	yes
Display	LCD, 2 × 16 signs
Overall dimensions (W × D × H)	240 × 260 × 165 mm
Weight	3.2 kg
Input current/power consump	otion 12 V, 4.4 A / 55 W
External power supply	Input AC 100–240 V 50/60 Hz; Output DC 12 V

### **BLOCKS (BUILT IN) CAPACITY:**

Block CH-1	20 × 0.5 ml +12 × 1.5 ml microtubes
Block CH-2	20 × 1.5 ml microtubes
Block CH-3	20 × 2 ml microtubes



 CH-100 with block CH-1
 BS-010410-BAI

 CH-100 with block CH-2
 BS-010410-CAI


 CH-100 with block CH-3
 BS-010410-UAI

Ice on block CH-2





### Heat up and cool down times for CH-100:















### Heat up and cool down times for CH3-150:

from 25°C	+	12 min	to 100°C	
from 25°C	+	60 r	nin	to -3°C

**B2-50** Ø48 mm × 2 sockets, depth 58 mm

### INTERCHANGEABLE THERMOBLOCKS:

•		
0	B10-16	Ø16 mm × 10 sockets, depth 56 mm
<b>3</b>	B6-25	Ø25 mm × 6 sockets, depth 40 mm
4	B23-1.5	23 sockets for 1.5 ml microtubes, depth 35 mm
6	B10-13	Ø13 mm × 10 sockets, flat bottom, depth 30 mm
0	B5-29	Ø29 mm × 5 sockets, flat bottom, depth 40 mm
0	B18-12	18 sockets for Ø12 mm round bottom tubes, depth 58 mm
<b>3</b>	B-1L	Mini water or beads bath, 55×100×73 mm
		and the second s

### Different block types can be provided on request

ORDERING INFORMATION:	Cat. number
CH 3-150 without blocks	BS-010418-AAA
Optional blocks:	
B2-50	BS-010418-AK
B10-16	BS-010418-BK
B6-25	BS-010418-CK
B23-1.5	BS-010418-DK
B10-13	BS-010418-LK
B5-29	BS-010418-KK
B18-12	BS-010418-EK
B-1L	BS-010418-WK
<b>1</b> B2-50 <b>2</b> B10-16	B B6-25 4 B23-1.5

# CH 3-150, Combitherm-2

Combitherm-2 **CH3-150** is specially designed to thermostabilise materials at temperatures from –3 °C to +150 °C according to analysis methods. To obtain useful functionality and decrease foot-print of instruments Combitherm-2 thermoblocks combined in a common electronic circuit board as well as inside a common external body. The front keyboard's left part is responsible for setting parameters for cooling plug-in blocks and the right part – for heating plug-in blocks. Both of them are regulated independently and can realise up to 16 programs, including temperature and time in each program. Peltier technology is used for cooling below room temperature; PCB is used for heating till +150°C.

Separation of cooling and heating parts from each other increases durability of the instrument and speed of temperature changing after setting a new program.

Heating Block Specifications:			
Temperature setting range	+25°C +150°C		
Temperature control range	5°C above ambient+150°C		
Setting resolution	1°C		
Stability	±0.1°C		
Temperature calibration coefficient range	0.9361.063 (±0.063)		

Cooling Block Specifications:			
Temperature setting range	−3°C +20°C		
Temperature control range	23°C below ambient 5°C below ambient		
Setting resolution	0.1°C		
Stability	±0.1°C		

General Specifications	
Digital time setting	1 min–99 h 59 min (increment 1 min)
Timer sound signal	yes
User adjustable programs (temperature and time)	16 (heating) +16 (cooling)
Display	LCD
Overall dimensions (W $\times$ D $\times$ H)	295 × 285 × 220 mm
Weight (without block)	5.6 kg
Nominal operating voltage	230 V, 50/60 Hz
Power consumption	430 W (1.8 A)

**7** B18-12











**6** B10-13



**6** B5-29



❸ B-1L

# **QB Series,** Dry Block Heating Systems with Interchangeable Blocks

Equipment presented on pages 68–69 is produced by Grant Instruments (Cambridge) Ltd. Biosan is the sole distributor of Grant Instruments products in CIS and the Baltic States (Latvia, Lithuania, Estonia) and the official distributor for a number of other regions.

A market-leading range of versatile, high-quality dry block heating systems with excellent temperature control, providing a source of precision heating for many sensitive analytical procedures.

### A premium product range at an affordable price:

- Accurate, reproducible and safe heating of your samples advanced temperature control combined with high quality, precision-engineered blocks providing excellent thermal contact;
- Versatile range of interchangeable heating blocks to fit any tube or plate you are using for your samples;
- Full range of models and options to cater for basic through to more sophisticated applications;
- · Wide range of accessories.





Product video







Model (Cat. Num.)	QBD1/QBD2/QBD4	QBH2	
Туре	Digital	Digital	
Number of blocks	1/2/4	2	
Temperature range	amb. +5°C to 130°C	amb. +5°C to 200°C	
Temperature setting range	+15°C to 130°C	+15°C to 200°C	
Temperature stability @ 37°C	±0.1	±0.1	
Temperature uniformity within the block @ 37°C	±0.1	±0.1	
Display / Resolution	LED / 0.1°C	LED / 0.1°C	
Safety: Overtemperature	Therm	nal fuse	
Timer with a sound alarm	1 min up to 72 h		
Heat up time from 25°C to 100°C	15 min		
Power consumption	150/300/600 W	300 W	
Power supply	120 V c	or 230 V	

SPECIFICATIONS

# **QB Series,** Dry Block Heating Systems with Interchangeable Blocks: Accessories

Interchangeabl	e blocks (Cat. Num.)	QBD1	QBD2	QBD4	QBH2	QBA1	QBA2
No. of blocks		1	2	4	2	1	2
QB-0 Plain bloo	ck without holes	+	+	+	+	+	+
QB-10 24 × 10	mm Ø holes, 50 mm hole depth	+	+	+	+	+	+
QB-12 24 × 12	mm Ø holes, 50 mm hole depth	+	+	+	+	+	+
QB-13 12 × 13	mm Ø holes, 50 mm hole depth	+	+	+	+	+	+
QB-16 12 × 16	mm Ø holes, 50 mm hole depth	+	+	+	+	+	+
QB-17H for 10 17 mm diam, 7	× Falcon tubes tall '5 mm deep	+	+	+	+	+	+
QB-18 12 × 18	mm Ø holes, 50 mm hole depth	+	+	+	+	+	+
QB-24 5 × 24 m bottles, 50 mm	nm Ø holes and universal I hole depth	+	+	+	+	+	+
	nl centrifuge tubes, s, 50 mm hole depth	+	+	+	+	+	+
QB-H 56 × 0.2 i	ml microtube, 14 mm hole depth	+	+	+	+	+	+
QB-E0 24 × 0.5	ml microtube, 30 mm hole depth	+	+	+	+	+	+
QB-E1 24 × 1.5	ml microtube, 35 mm hole depth	+	+	+	+	+	+
QB-E2 24 × 2.0	ml microtube, 35 mm hole depth	+	+	+	+	+	+
QB-E5 12 x 5.0 16.7 mm diam	ml microtube, 53.5 mm hole depth, eter	+	+	+	+	+	+
QB-DN Dolphir 24 × Ø 11.13 m	n nose tube Im to Ø 6.1 mm	+	+	+	+	+	+
External Pt100	0 temperature probe						
QBEP	Standard probe. For in-sample or in-block temperature control; encased in stainless steel sheath, Ø 3 mm × 30 mm long, with 350 mm of cable	+	+	+	+	-	-
QBEP-WM	Short-form probe. For in-sample or in-block temperature control; encased in stainless steel sheath, Ø 3 mm × 14 mm long, with 350 mm of cable	+	+	+	+	-	-
	cks of molecular biology and biotechr ocks 140 × 100 × 75 mm supplied with			ool			
QDP-H	96 holes in microplate configuration for 0.2 ml microplates, strips or individual tubes. Uniformity ± 0.3°C within tubes across the block; 6.2 mm Ø holes, 14 mm hole depth	-	+	-	+	-	+
ODP-FI	Universal block for standard 96-well plates (u-well, v-well, flat bottom, high temperature) Uniformity ± 0.5°C between wells; supplied with hinged, double layer lid to create an insulated incubation chamber	-	+	-	+	-	+
Safety covers (	not required with QDP-FL Microtiter b	locks)					
	Made from tough clear acrylic for maxi- mum visibility whilst preventing acci- dental touching of a hot block or con- tamination of samples from splashes. Clearance height 85 mm	QBL1	QBL2	QBL4	QBL2	QBL1	QBL2



# CRFT, Controlled rate freezer and thawer

Liquid nitrogen and cryogen free controlled rate freeze and thaw system for cryopreservation – **CRFT.** 

The **CRFT** brings accuracy, precision and reproducibility to biological cryopreservation. This unit is ideal for research into the cryopreservation of a wide range of material including embryos, stem cells, mammalian cells, spermatozoa, antibodies, tissue sections and rodent organs.

The freezing and thawing rate is precisely controlled, ensuring accuracy and reproducibility throughout the freezing and thawing profiles, especially for the important nucleation/seeding phase. This provides optimal recovery of cells on thawing. Operation is simple and can be carried out with or without a PC.

Max plate temp	30°C
Min plate temp	-100°C
Temp accuracy	± 0.5°C
Temp stability	± 0.1°C
Temp uniformity	± 1.0°C Ref
Max cooling rate	10°C/min
Max thawing rate	10°C/min from -100°C
Max temp stabilisation time	5 min
Overal dimensions	345 × 290 × 380 (lid off)
(W×D×H)	345 × 290 × 460 (lid on)
Weight	5 kg
Power supply	100-240V, 50-60Hz (24V DC)

### **APPLICATIONS**

The **CRFT** is highly versatile and can be used for the cryopreservation research of a wide range of samples in cryovials, straws, bags, microplates and Matrix-96-well block plates in the following areas:

- · Transgenic embryology research;
- · Stem cell research;
- Clinical and research samples, e.g. lymphocytes and tissue cell lines in conventional cryovials;
- Various mammalian cells including cardiomyocytes, adipose, liver and muscle;
- · Cord blood derived stem cells;
- · Adherent cells and stem cells in microplates;
- · Cell suspensions in numbered/barcoded arrays;
- · System integrations;
- · Suitable for applications in veterinarian IVF.



### **ACCESSORIES**

H00 - Head for 16×0.5ml straws;

H01 - Head for 18×0.3ml IMV straws;

H02 - Head for 55×1.8ml cryovials (0.5ml max fill);

H03 - Head for Flat universal plate;

H05 - Head for 1×SBS microplate;

H06 - Head for 55×1.8ml deep cryovials (1ml max fill);

**CRFTCRYOPEN** – Cryopen nucleating tool;

**CRFTCARTRIDGE** - Replacement cartridge for cryopen.

### **KEY FEATURES AND BENEFITS**

- Accurate and reproducible control of cooling rates and sample temperatures;
- · A focus on maximising cell survival rates;
- · Controlled and customisable freezing and thawing;
- · Linear and non-linear freezing profiles;
- Liquid nitrogen, alcohol and cryogen free no need for 'topping up';
- · Interchangeable heads for flexibility (ordered separately);
- Free of charge software allows you to design, track, export and report temperature profiles;
- Easy to use and samples can be nucleated/seeded in-situ;
- Low running costs: estimated at 1% of liquid nitrogencontrolled rate freezing and runs on a 24V supply;
- · Command protocols available.

# THERMOSTATIC EQUIPMENT:

# WATER BATHS, ORBITAL/LINEAR SHAKING BATHS, UNSTIRRED WATER BATHS, HEATING/COOLING CIRCULATORS



Optima<sup>™</sup> Series
Stirred Thermostatic Baths
and Heating Circulators


**Ecocool**Energy Efficient Refrigerated /
Heating Circulating Baths



**TDB-120** Dry block thermostat

# How to choose thermostat?

#### How to choose thermostat?



SPECIFICATIONS

#### WB-4MS, Stirred water bath

Stirred water bath **WB-4MS** is designed for chemical, pharmaceutical, medical and biological laboratory research, for processes requiring constant temperature ranging from ambient temperature to 100°C.

**WB-4MS** provides increased temperature stabilization (up to 0.1°C) due to a built-in magnetic stirrer (speed control range 250–1000 rpm).

Easy setup, high-temperature maintenance accuracy, compact size and attractive modern design make this water bath widely used.

Tank capacity	41
Temperature setting range	+25°C +100°C
Temperature control range	5 °C above ambient +100°C
Temperature setting resolution	n 0.1°C
Temperature stability	±0.1°C
Temperature uniformity @ +3	7°C ±0.1°C
Stirring speed control range	250-1000 rpm
Digital time setting	1 min–96 h/non-stop (increment 1 min)
Timer sound signal	yes
Display	LCD, 2 × 16 signs
Digital setting of temperature	, time and mixing speed
Plastic lid with stainless steel i	nterior included
Quiet operation	
Working volume	235 × 135 × 110 mm
Overall dimensions (W×D×H)	340 × 270 × 250 mm
Weight	3.4 kg
Nominal operating voltage	230 V, 50/60 Hz or 120 V, 50/60 Hz
Power consumption	230 V, 50 Hz/600 W (2.6 A) 120 V, 60 Hz/670 W (5.6 A) 100 V, 50/60 Hz/600 W (6.0 A)





WB-4MS with two adapters





#### ORDERING INFORMATION:

Cat. number

#### WB-4MS with base BP-1 and lid

BS-010406-AAA

#### Optional accessories:

Racks:	Diameter/tube volume	Capacity	Cat. number
1 TR-5/30	Ø 30 mm	5 tubes	BS-010406-KK
<b>2</b> TR-16/19	Ø 16–19 mm	16 tubes	BS-010406-FK
<b>3</b> TR-30/13	Ø 10–13 mm	30 tubes	BS-010406-IK
<b>4</b> TR-44/11	2/1,5 ml	44 tubes	BS-010406-JK

OLS26

#### Combined Orbital/Linear Shaking Bath OLS26



Patented, combined orbital and linear shaking mechanism of the **OLS26** allows optimisation of aeration and shear forces mixing for reproducible results.

- · Precision digital temperature control;
- 0°C to 99°C operating range;
- Stability ±0.1°C;
- · Easy changeover from linear to orbital shaking;
- · Adjustable shaking speed and stroke length;
- Polycarbonate lid included as standard;
- Drain tap for convenient emptying;
- · 3 year warranty;

Tank size

TU26 included, other trays sold separately.





Tarik Size		201
Minimum working dep	th	70 mm
Temperature control range		ambient +5 to 99°C. 0 to 99°C with accessory cooling
Temp. uniformity (DIN	12876-3) @ 70	D°C ±0.1°C
Temp. stability (DIN 12	876-3) @ 70°C	±0.1°C
Display		vidual displays and controls for emperature and shaking speed)
Orbital and Linear sha	king speed	20 to 200 rpm (depending on load)
Orbital shaking radius		9 mm
Shaking speed display	resolution	1 rpm
Linear shaking stroke l	ength	18, 28, and 36 mm
Shaking tray area		380 × 235 mm
Timer		1 to 999 min
Dimensions (W×D×H)		335 × 565 × 325 mm
Heater power 120 V/23	30 V	1.05/1.4 kW
Drain tap		yes
Safety		over temperature protection/ low liquid level cut-out
Supply voltage		110–120 V or 220–230 V
ORDERING INFORMA	TION:	Cat. number

Equipment presented on pages 64–65, 75–87 is produced by Grant Instruments (Cambridge)
Ltd. Biosan is the sole distributor of Grant Instruments products in CIS and the Baltic States (Latvia, Lithuania, Estonia) and the official distributor for a number of other regions.

OLS26 with TU26 tray

DESCRIPTION

#### Linear shaking bath - LSB Aqua Pro range

World-renowned shaking water baths. High quality, robust design with unique magnetically coupled shaking mechanism for maximum reliability, consistency and quiet operation. Extensive range of accessories to provide the right solution for your application. Varied vessels types can be securely held using high-quality springs, clamps or racks.

#### **FEATURES**

- Ambient +5°C to 99°C operation;
- Stability ±0.1°C;
- · Choice of two models 12 and 18 litres;
- · Drain tap for convenient emptying;
- 3 year warranty;
- · Polycarbonate lid included;
- Extensive choice of accessory shaking trays. Tray sold separately.



LSB18

	9.2 kg W: 360 mm D: 385 mm H: 425 mm		11.2 kg W: 335 mm D: 565 mm H: 425 mm			
Tank size		12		18 I		
Minimum working depth		60 r	mm			
Temperature range	ambient +5 to 99°C					
Uniformity (DIN 12876-3) @ 70°C	±0.1°C					
Stability (DIN 12876-3) @ 70°C	±0.1°C					
Display	LED					
Linear shaking speed		20 to 200 strokes/mir	(depending on l	oad)		
Shaking speed display resolution		1 strok	es/min			
Linear shaking stroke length		20 r	mm			
Shaking tray area	240 >	× 235 mm	420	) × 235 mm		
Timer		1 to 99	99 min			
Heater power 120/230V	0.8	/0.8 kW	1.	05/1.4 kW		
Drain tap		y€	es			
Safety		over-temperature prote	ection/low liquid o	cut-out		
Supply voltage	110-120 V or 220-230 V					

LSB12



#### ORDERING INFORMATION:

Cat. number

LSB12, Linear shaking bath 12 I with TU12 tray

LSB12

LSB18, Linear shaking bath 18 I with TU18 tray

LSB18

#### Accessories for Shaking Baths: LSB 12, LSB 18 & OLS 26


Accessories <b>LSB</b> and <b>OL</b>	OLS26	LSB12	LSB18	
		Cat	alogue numb	er
A FA	Universal tray with adjustable springs. Highly versatile for a variety of vessel types.	TU26	TU12	TU18
	Flask/plate tray – with threaded holes to accept flask clamps or holder for deep well plates (≥2ml). See option below.	TF26	TF12	TF18
	Test tube tray – compatible with SR racks or can be used alone to accommodate bags and miscellaneous vessels. See rack option below.	TS26 (holds up to 5 SR racks)	TS12 (holds up to 3 SR racks)	TS18 (holds up to 5 SR racks)
	Base tray – perforated stainless steel, allows bath to be used as an unstirred bath.	SBT26	SBT12	SBT26
	Cooling coil – the source of constant cooling to enable bath to be operated at or below ambient, down to 0 °C. LS200 lid (with an access hole for cooling coil) recommended.	CC26R	-	
	Heat exchange coil – attach to a cold water supply or refrigerated circulator. Can be used down to 2°C above the temperature of the coolant. LS200 lid (with an access hole for cooling coil) recommended.	CW26		-
	Stainless steel sloping lid, gabled.	LS200	LU14	LU28
	Replacement polycarbonate lid, clear, gabled.	AQL26	AQL12	AQL26

Cat. number	Description	<b>OLS26</b> Capacity	LSB12 Capacity	LSB18 Capacity
SC-25	for 25 ml flask	28	18	33
SC-50	for 50 ml flask	24	14	26
SC-100	for 100 ml flask	15	9	17
SC-250	for 250 ml flask	8	5	14
SC-500	for 500 ml flask	6	4	6
SC-1000	for 1,000 ml flask	3	2	4
SH-DWP	1 × deep well plate	4	1	4

Test tube racks / microtube racks for TS tray					
Cat. number	Tube diameter (mm)	Rack capacity			
SR-10	10	48			
SR-13	13	44			
SR-16	16	24			
SR-19	19	21			
SR-25	25	12			
SR-30	30	10			
Cat. number	Microtube size (ml)	Rack capacity			
SR-SE	0.5	119			
SR-LE	1.5	48			

#### **Unstirred Water Bath**





- The reliability, quality and consistent performance of Grant products have made Grant a leading manufacturer of water baths for decades;
- A new era for Grant water baths now all models from basic to advanced with digital controls;
- Proven performance technology to deliver temperature control you can rely on;
- Set and Forget<sup>™</sup> technology minimal bath setup, maximum time for your work.

#### **SUB Aqua Pro** Digital Unstirred Water Bath



Built to the highest standard and specifications and incorporating the latest technology, the SUB Aqua Pro advanced water bath range supports even the most demanding applications, that require accurate temperature control. Choose from eight models with

- base tray and lid included as standard.Ambient +5°C to 99°C operation;
- Set and Forget<sup>™</sup> technology fast heat-up, accurate temperature control;
- Stability ±0.2°;
- Adjustable over temperature alarm protect samples from overheating;
- · Advanced dry start and run-dry protection;
- · Three programmable temperature presets;
- 3 year warranty.

	SUB Aqua Pro digital unstirred water bath range – summary of specifications							
	SAP2	SAP2S	SAP5	SAP12	SAP18	SAP26	SAP34	SAPD
	2.5 kg W: 185 mm D: 200 mm H: 200 mm	3 kg W: 335 mm D: 215 mm H: 150 mm	3 kg W: 335 mm D: 215 mm H: 200 mm	6 kg W: 360 mm D: 380 mm H: 225 mm	9.5 kg W: 335 mm D: 590 mm H: 275 mm	9 kg W: 335 mm D: 590 mm H: 275 mm	14.5 kg W: 335 mm D: 590 mm H: 275 mm	9 kg W: 545 mm D: 380 mm H: 225 mm
Tank capacity	21	2 l (shallow)	51	12 I	18 I	261	34	5   & 12
Temperature range				ambient	t°C + 5 to 99			
Temp. display and setting resolution		0.1℃						
Temp stability (DIN 12876) @ 70 °C				±	0.2°C			
Temperature setting/energy regulation				c	ligital			
User adjustable over temp. alarm					+			
Fixed thermal cut-out					+			
Dry start/boil dry protection					+			
Programmable temp. presets					3			
Countdown timer with audible alarm				1 to	999 min			
Working area D×W (mm)	117 × 131	139 × 289	131 × 281	281 × 306	485 × 281	481 × 278	635 × 281	131 × 281 & 281 × 306
Minimum fill level	50 mm	32 mm	50 mm	50 mm	50 mm	70 mm	70 mm	50 mm
Maximum fill level			2	25 mm below	the top of the	e tank		
Drain tap included	_	_	-	+	+	+	+	+
Heater power 120 V/230 V kW	0.25/0.25							
Supply voltage V				120	or 230			

#### SUB Aqua Pro Digital Unstirred Water Bath

OPTIONS A	OPTIONS AND ACCESSORIES							
SAP2	SAP2S	SAP5	SAP12	SAP18	SAP26	SAP34	SAPD	
21	2 l (shallow)	51	12	18	26 I	34 I	5 &12	
Replaceme	ent polycarbo	nate transpare	nt lids*					
AQL2	AQL5	AQL5	AQL12	AQL26	AQL26	-	AQL5, AQL12	
Directs cor	ndensation av	vay from imme	rsed vessels, avoid	s contamination,	, reduces evapo	ration and saves	s energy	
Stainless s	teel sloping li	ds*						
-	LU6	LU6	LU14	LU28	LU28	LU36	LU6 & LU14	
Flat lids*								
-	-	LF6 (2 ring sets)	LF14 (4 ring sets)	LF28 (6 ring sets)	LF28 (6 ring sets)	LF36 (8 ring sets)	LF6 / LF14	
With ring s	ets of variabl	e hole diamete	r to accommodate	tall vessels whils	t reducing evap	oration		
Polypropyl	ene spheres*	(packs per bat	h)					
1 × PS20	1 × PS20	1 × PS20	1 × PS20	2 × PS20	2 × PS20	3 × PS20	2 × PS20	
Useful alterna	ative to a lid, min	imises evaporation	and heat loss whilst allo	wing easy access to	vessels in the bath; p	particularly useful fo	r tall vessels	
Raised she	lves – reversi	ble, allows two	shelf depths. h = sl	helf height above	e tank base (mm	1)		
-	-	_	RS14H (h 40 or 78) shelf covers half area of SAP12	RS18H (h 40 or 135) shelf covers half area of SAP18	RS28H (h 45 or 135) shelf covers half area of SAP26	RS36H (h 45 or 135) shelf covers half area of SAP34	RS14H (h 40 or 78) shelf covers half area of SAPD	
Racks (no.	per bath)							
_	_	1 × J2	2 × J2	4 × J2	4 × J2	6 × J2	1 + 2 × J2	
Choice of e	eight variants	to accommoda	te different tube d	iameters and mi	crotubes (see be	elow)		
Replaceme	ent base trays							
AQBT2	AQBT5	AQBT5	AQBT12	AQBT26	AQBT26	SBT36	AQBT5 & AQBT12	

Required if flat-bottomed flasks are to be placed directly on the base of the bath and to promote thermal convection in the bath

<sup>\* –</sup> Lid or spheres recommended for use above  $60^{\circ}\text{C}$ 

Unstirred Bath R	acks				
J2 Racks	Tube size Ø	Capacity	J2 Racks	Tube size Ø	Capacity
J2-10	10 mm	84	J2-25	25 mm	18
J2-13	13 mm	55	J2-30	30 mm	12
J2-16	16 mm	36	J2-SE	0.5 ml	105
J2-19	19 mm	32	J2-LE	1.5 ml	65

# **Optima™ Series**, Stirred Thermostatic Baths and Heating Circulators







Heating Circulators Specifications on page 82 and all available accessories on page 84 A cost-effective range of multi-purpose systems combining Grant's legendary quality and reliability. Precise temperature control for a wide range of laboratory applications.

- Accurate and safe temperature control for samples and users;
- Intuitive programming and thoughtful design features makes working with Grant heated baths and circulators easy;
- Robust, durable construction for longevity, reliability and long-term low cost of ownership;
- A complete range 32 models to cover basic through to sophisticated needs, each model represents excellent value for money.

#### **APPLICATIONS**

Grant stirred baths and circulators provide a source of precision heating and cooling for many routines and sensitive analytical procedures including sample incubation, calibration and quality control testing. All models from the TC120 upwards are suitable for unnecessary both open and closed-loop circulators (i.e. remote vessel open or closed).

For more powerful heating requirements, i.e. above 200 °C, contact marketing@biosan.lv for advice.

#### Model selection (see next page):

Any of the four **Grant Optima**<sup>™</sup> digital thermostats can be combined with any of eight Grant tanks (five stainless steel and three plastic) to provide a choice of 32 models.

#### **Optima™ Series,** Heating Circulators Specifications









Grant Optima™ Heating Specifications	Circulators	General pur	pose Digital	Digita	al High Performance
Specifications		T100	TC120	TX150	TXF200
Stability (DIN 12876) @ 70°C	°C	±0.05	±0.05	±0.01	±0.01
Uniformity (DIN 12876) @ 70°C	°C	±0.1	±0.1	±0.05	±0.05
Setting resolution	°C	0.1	0.1	0.1 (	(0.01 with Labwise™)
Display		4 digi	it LED	fu	ll colour QVGA TFT
Timer function		-	1 to 6,000 min	11	min to 99 h 59 min
No. preset temperatures		3	3	3	3
Re-calibration points		2	2	5	5
Offset adjustment		-	_	+	+
Socket for external probe (	TXPEP, TXSEP)	-	_	+	+
Communication interface		-	_	USB & RS232	USB & RS232
Programmable		-	_	remote via PC/laptop 1 program/ 30 segments	direct via user interface or remote via PC/laptop 10 programs / 100 segments
Relays		-	_	1	1
Safety	overtemperature	fixed		adjustable	e cut-out
Safety	fluid level – float switch	+	+	+	+
Alarms (can be configured	to switch a relay)	-	high, without relay	high and low	high and low
Heater power 230 V	kW	1.3	1.3	1.9	1.9
Electrical power 230 V	kW	1.4 (50-60 Hz)	1.4 (50 Hz)	2.0 (50 Hz)	2.0 (50-60 Hz)
Height above tank rim	mm	200	200	200	200
Depth below tank rim	mm	135	135	135	135
Grant Optima™ thermost	at pumps (integral)				
Maximum pressure	water, mbar	-	210	310	530
Maximum flow	water, l/min	-	16	18	23 (adjusted flow rate)
Pipe bore	inlet/outlet, mm	-	6/11	6/11	6/11
Dimensions (H×D×W)	mm			315 × 145 × 115	

7			
\	ODDEDING	INFORMATION	
	ORDERING	INFORMATION	10

Cat number	T100 EURO	TC120 EURO	TX150 EURO	TXF200 EURO	
Cat. number:	1 100 EURO	IC120 EURO	1X150 EURO	TXF200 EURO	

# **Optima™ Series,** Water Bath Combinations and Accessories

Capacity (l)	Outer tank dimensions  1. Dimensions (H×D×W) Weight (kg)  2. Working area (D×W)  3. Min/max fluid depths  4. Inner tank dimensions (H×D×W)	<b>T100</b> Temperature setting range	TC120 Temperature setting range	<b>TX150</b> Temperature setting range	<b>TXF200</b> Temperature setting range
ST5 – 5 l Stainless steel	1. 215 × 335 × 187 mm, 2.9 kg 2. 150 × 260 mm 3. 85/140 mm 4. 150 × 300 ×150 mm	T100-ST5 amb.+15 to 100 °C	TC120-ST5 0 to 120 °C	TX150-ST5 0 to 150 °C	TXF200-ST5 0 to 200 °C
ST12 – 12 l Stainless steel	1. 215 × 332 × 360 mm, 4.5 kg 2. 205 × 300 mm 3. 85/140 mm 4. 150 × 325 × 300 mm	T100-ST12 0 to 100 °C	TC120-ST12 0 to 120 °C	TX150-ST12 0 to 150 °C	TXF200-ST12 0 to 200 °C
ST18 – 18 l Stainless steel	1. 215 × 545 × 340 mm, 7.3 kg 2. 385 × 300 mm 3. 75/130 mm 4. 150 × 505 × 300 mm	T100-ST18 0 to 100 °C	TC120-ST18 0 to 120 °C	TX150-ST18 0 to 150 °C	TXF200-ST18 0 to 200 °C
ST26 – 26 l Stainless steel	1. 270 × 535 × 340 mm, 7.7 kg 2. 385 × 300 mm 3. 125/180 mm 4. 200 × 505 × 300 mm	T100-ST26 0 to 100 °C	TC120-ST26 -15 to 120 °C	TX150-ST26 -15 to 150 °C	TXF200-ST26 -15 to 200 °C
ST38 – 38 l Stainless steel	1. 260 × 733 × 338 mm, 11.9 kg 2. 575 × 300 mm 3. 125/180 mm 4. 200 × 690 × 300 mm	T100-S38 0 to 100 °C	TC120-S38 -15 to 120 °C	TX150-S38 -15 to 150 °C	TXF200-S38 -15 to 200 °C
P5 – 5 l Plastic	1. 180 × 323 × 220 mm, 2.2 kg 2. 120 × 150 mm 3. 85/140 mm 4. 155 × 240 × 160 mm	T100-P5 amb.+15 to 99 °C	TC120-P5 amb.+15 to 99 °C	TX150-P5 amb.+15 to 99 °C	TXF200-P5 amb.+15 to 99 °C
P12 – 12 l Plastic	1. 180 × 412 × 340 mm, 3.4 kg 2. 210 × 280 mm 3. 85/140 mm 4. 155 × 325 × 280 mm	T100-P12 amb.+5 to 99 °C	TC120-P12 amb.+5 to 99 °C	TX150-P12 amb.+5 to 99 °C	TXF200-P12 amb.+5 to 99 °C
1. 180 × 589 × 340 mm, 5.1 kg P18 – 18 l 2. 375 × 280 mm Plastic 3. 85/140 mm 4. 155 × 510 × 290 mm		T100-P18 amb.+5 to 99°C	TC120-P18 amb.+5 to 99°C	TX150-P18 amb.+5 to 99°C	TXF200-P18 amb.+5 to 99°C
OPTIONS AND AC	CESSORIES				
Labwise™ PC soft	ware (optional)				
Allows two-way co programming and	ommunication for status display, d data capture	_	_	+	+
External probes (optional)					
TXPEP flexible plastic probe, 3 m cable TXSEP stainless steel probe, 3 m cable		- -		+ +	++
Remote switching device (optional)					
For switching app	liances on and off (up to max. 8 Amps)	-	-	1	2
Vertical turbine pumps (optional)					
Low noise, compa	act design. Supplied with pipe connections	and special lid	for fitting to tar	nk, pipe bore 12	.7 mm
VTP 1 max. pressure 1,000 mbar Required only what max. flow 9 l/min + Required only what a higher pressure 1,650 mbar.					elivered by the

#### **Optima™ Series,** Water Bath Accessories

ACCESSORII	ΞS						
	Lids to help reduce evaporation/ heat loss and avoid sample	Polypropylene spheres (no. of packs required, 300	Rack systems to optimise use of available bath capacity (no. of racks	Raised shelves to allow shallow vessels to be accommodated	to allow sys room temp coil dipped	y cooling systems ystems to operate at or below iperature by means of cooling d into the bath; designed for mpact on working area	
	contamination	spheres in one pack)	accommodated)		the bath co	coolers cooling ted to a n unit by a e. Extract uously, with	Heat exchange coil Designed to be attached to a supply of cooling tap water or a refrigerated circulator
					C1GR (0 to 40°C)	C2GR (-15 to 40°C)	CW5 (2°C above coolant temperature)
ST5 – 5 L stainless steel	STL5 flat stainless steel	1 × PS20	1×QR	=	7	-	7
ST12 – 12 L stainless steel	<b>STL12</b> gabled, hinged (removable) stainless steel	1 × PS20	2×VR	RS14	7	-	7
ST18 – 18 L stainless steel	STL26 gabled, hinged (remo steel	2 × PS20	4×VR	RS22	7	_	7
ST26 – 26 L stainless steel	STL26 gabled, hinged (remo steel	2 x PS20	4×VR	RS28	7	7	3
ST38 – 38 L stainless steel	STL38 gabled, hinged (remo steel	3 × PS20	6×VR	RS28 or RS38	-7	7	3
P5 – 5 L plastic	PL5 flat, stainless steel	1 × PS20	1×QR	-	-	-	_
P12 – 12 L plastic	PL12 curved plastic	1 × PS20	2×VR	RS14	_	_	_
P18 – 18 L plastic	PL18 curved plastic	2×PS20	4×VR	RS22		_	-

#### Ecocool

#### Energy Efficient Refrigerated / Heating Circulating Baths



- Choice of two models, temperature range -30 °C to +150 °C (model dependent);
- 3 year warranty with renowned service and support, no registration required;
- · Active cooling through the whole temperature range;
- True energy saving of up to 80% against standard compressor units.

A new range of innovative, eco-friendly, refrigerated heating circulating baths offering significant running cost savings whilst delivering powerful cooling.

All products in the Ecocool range are supplied assembled as ready to use kits, complete with accessory hosing, clips and connectors as standard.

H: 640 mm D: 430 mm W: 245 mm		ecocool 100R	ecocool 150R
Temperature range	°C	-20 to 100	-30 to 150
Temperature stability	±°C	0.05	0.02
Refrigerant		R2	90
Flow rate (max)	L/min	17	14 – 22 (adjustable)
Pump pressure (max)	mbar	250	530
Tank volume	L	4.5	5.5
Working Area (D × W)	mm	280 × 140	280 × 150
Min/Max liquid level	mm	105/140	160/190
Cooling power	@ 20°C W @ 0°C W @ -10°C W @ -20°C W	250 200 100 50	450 350 300 200
Programs		-	1 × 30 segments via Labwise™
Communication interface		-	USB
Temperature probe socket		-	6 pin mini DIN
Display		4 digit LED	Full colour QVGA TFT
Languages		-	5 (EN, FR, DE, IT, ES)
Timer		1 min to 99	hrs 59 mins
Temperature presets			3
Alarms		High	High and low
Electrical power (max) kW	120V/230V	2.16/2.07 (50-60 Hz)	2.28/2.76 (50-60 Hz)
Safety		Adjustable over temperature cut-out	
Ready to use kits		Assembled and supplied with standard	tubing, insulation, clips and connectors

#### **Ecocool**

#### Energy Efficient Refrigerated / Heating Circulating Baths

#### **APPLICATIONS**

- PHARMACEUTICAL Mini pilot plant reactors
- EDUCATION Rotary evaporator cooling, replacement of running tap water cooling, immersing small samples, photometry, chromatography systems
- INDUSTRIAL QC testing, sample preparation, general cooling, reaction chemistry, temperature control, semi-conductor manufacturing, rheometry
- FOOD Refractometry
- LIFE-SCIENCE Electrophoresis cooling
- HIGH TEMPERATURE COOLING Active up to 150°C



Options and accessories	ecocool 100R	ecocool 150R	
Labwise™ PC software (optional)			
Allows two-way communication for status display, programming and data capture + USB cable provided	_		
External probes (optional)			
PEP plastic probe	-	+	
SEP stainless steel probe	-	+	
Vertical turbine pumps (optional) when pump is fi	tted, available wo	orking area is reduced.	
Low noise, compact design. Supplied with pipe connspecial lid for fitting to tank, pipe bore 12.7 mm	Required only where application demands higher pressure than that delivered by the		
VTP1-LT max. pressure 1,000 mbar; max. flow 9 l/min		internal to maintain flow.  Note: The optional VTP pumps will transfer additional heat to the baths and reduce the	
VTP2-LT max. pressure 1,650 mbar; max. flow 12 l/min	*	net cooling power of the refrigeration unit. The above figures must be taken into consideration when choosing the refrigeration unit, when ordering a VTP pump, please specify which refrigeration base unit it is to be used with.  Note: Other sizes of heat exchange coil can be made to your specification. Contact us for further information.	

#### **Hose Kits**

HOSE100 General purpose hose kit: –40 to 100  $^{\circ}\text{C}$  HOSE200 High temperature hose kit: –50 to 200  $^{\circ}\text{C}$ 



Hose kit  $2 \times 2$  m, assembled with Optima<sup>™</sup> pump outlet plate and simple hose clips, no tools required.

#### Optima™ R series, Refrigerated Thermostatic Bath and Circulators



We recommend using the following liquids with refrigerated thermostatic baths and circulators:

- -50 to 50°C: Silicone oil low viscosity (Bayer silicone M3);
- -30 to 30°C: 50% water 50% antifreeze (inhibited ethylene glycol);
- 0 to 30°C: 80% water 20% antifreeze (inhibited ethylene glycol);
- 5 to 99.9°C: Water.

Cost-effective and efficient multi-purpose systems for low-temperature applications.

- Powerful precision cooling, whether used in openloop or closed-loop format
- Combining legendary quality, reliability and design for everyday usage – useful features, straightforward maintenance, compact design
- Robust, durable construction for longevity, reliability and long-term low cost of ownership
- · Up to 4 years warranty

Grant low-temperature circulators provide a source of precision cooling for many sensitive analytical procedures, including spectrophotometry, viscometry, refractometry and electrophoresis. They are suitable for use in both open and closed-loop circulation (i.e. remote vessel open or closed).

Alternatively, Grant RC series of recirculating chillers (closed circulators) can be used. These are generally needed for more powerful cooling requirements, e.g. the removal of mechanical or electrical heat produced in apparatus or machinery. Please contact marketing@biosan.ly for advice.

#### Model selection:

The R4R refrigeration range consist of one refrigeration unit which can be combined with four heating circulators to offer a temperature range of  $-30^{\circ}$ C to  $100^{\circ}$ C.

Capacity (l) Outer tank dimensions	Working area (L×W)     Min/max liquid depths     Weight	<b>T100</b> H: 333 mm D: 172 mm W: 120 mm	TC120 H: 333 mm D: 172 mm W: 141 mm	<b>TX150</b> H: 342 mm D: 172 mm W: 141 mm	<b>TXF200</b> H: 342 mm D: 172 mm W: 141 mm
R4R – 20 l stainless steel H: 550 mm D: 515 mm W: 393 mm; <i>Cat.num.: R4</i>	• 230 × 305 mm • 80/140 mm • 40.6 kg	T100-R4R 0°C to 100°C	TC120-R4R -25°C to 100°C	TX150-R4R -30°C to 100°C	TXF200-R4R -30°C to 100°C
Options and accessories					
Labwise™ PC software (optional)					
Allows two-way communication fo programming and data capture +		-	-	+	+
External probes (optional)					
TXPEP flexible plastic probe, 3 m c TXSEP stainless steel probe, 3 m c	- -	- -	++	+ +	
Remote switching device (optional)					
For switching mains power appliances on and off (up to max. 8 Amps)		-	-	1	1
Vertical turbine pumps (optional)					
Low noise, compact design. Suppl lid for fitting to tank, pipe bore 12	and special		y where applica		
VTP 1 max. pressure 1,000 mbar; i	+	a higher pressure than that delivered by the internal pump to maintain flow			
VTP 2 max. pressure 1,650 mbar; i	max. flow 12 l/min	+			



# MAGNETIC STIRRERS, OVERHEAD STIRRER



MM-1000 Overhead Stirrer Multi Mixer

#### MS-3000 and MMS-3000, Magnetic Stirrers

**MS-3000** and **MMS-3000** are compact magnetic stirrers with stainless steel working surface. Units provide stirring of liquids with the rotation speed of magnetic element up to 3,000 rpm. Up to date, it is the highest value of the maximal speed for magnetic stirrers of global producers.

Strong magnets hold the driven magnetic element firmly in the magnetic clutch. Stirring is performed without undesirable heating and noise.

Enclosures of stirrer **MS-3000** are made of strong steel and painted with powder enamel, which is chemically resistant to acids and alkali.


The stirrers are supplied with a cylinder-shaped magnetic stirring bar  $(6 \times 25 \text{ mm})$  encapsulated in PTFE for universal use.

**MMS-3000** is equipped with a detachable stand for supporting various sensor elements (temperature, pH and others) inside the stirred liquid.

Magnetic stirrer is an ideal laboratory instrument for PH-metering, extraction and dialysing with small quantities of substances.

Operation temperature range +4°C to +40°C (from cold rooms to incubators) at maximal relative humidity 80%.





#### MS-3000 and MMS-3000, Magnetic Stirrers

	MS-3000	MMS-3000
Speed control range	0–3,00	0 rpm
Stirring volume up to (H <sub>2</sub> O)	5 L	20 L
Working surface material	Stainles	ss steel
SR-1, attachable stand size	-	Ø 8 × 320 mm
Max. length of magnetic stirring element (bar)	50 mm	70 mm
Stirring liquid viscosity	up to 1,170 mPa×s	
Maximum continuous operation time	24 h	
Operation in closed laboratory rooms	at ambient temperatu	re from +4°C to +40°C
Working plate size	110 × 110 mm	Ø 160 mm
Overall dimensions (W×D×H)	120 × 150 × 65 mm	185 × 230 × 75 mm
Weight	0.8 kg	1.5 kg
Input current/power consumption	12 V, 220 mA / 2.6 W 12 V, 250 mA / 3 W	
External power supply	Input AC 100-240 V, 50	/60 Hz; Output DC 12 V

MMS-3000



MMS-3000



ORDERING INFORMATION:

MS-3000 blue (on request)

MS-3000 white

MMS-3000

Cat. number

BS-010301-AAF

BS-010301-ABF

BS-010305-AAF

Optional accessories for MMS-3000:

HTP-1, Holder for temperature probe (see page 93)

BS-010309-FK

#### MSH-300 and Intelli-Stirrer MSH-300i, Magnetic Stirrers with hot plate

MSH-300 and Intelli-Stirrer MSH-300i are magnetic stirrers of the new generation. Enclosures of stirrers are made of metal painted with powder enamel chemically resistant to acids and alkali. The stirrers are equipped with a detachable stand for supporting various sensor elements (temperature, pH and others) inside the stirred liquid.

The stirrers are supplied with a cylinder-shaped magnetic stirring bar (6 × 25 mm) for universal use covered with Teflon.

Units are equipped with overheat protection providing an automatic switch-off of the device when overheating for the set temperature difference occurs.

Magnetic stirrers with heating can be used for laboratory operations such as organic synthesis, extraction, analysis of oil products, pH-measurements, dialysis, soil suspending, preparing buffer solutions, etc.

Additional protection disables the heating if the temperature of the plate exceeds the set temperature for 30°C.

Operation temperature range +4°C to +40°C (from cold rooms to incubators) at maximal relative humidity 80%.



Intelli-Stirrer MSH-300i is a digital version of magnetic stirrer with heating; it is designed for laboratories with higher requirements. It offers digital setting and control of temperature and rotation speed.

A powerful magnet allows mixing solutions with glycerine viscosity level. Maximum volume of stirred liquid (water) is 20 litres.

An external probe provides direct control of the stirred liquids temperature.

External temperature probe:			
Probe type	thermocouple		
Connection	type K		
The cable is covered with Teflon, me	echanically strong,		
elastic and chemically stable against oils, acids,			
aggressive reagents and liquids			
Cable length	1 m		
Operation temperature range	-50°C to +250°C		





#### MSH-300 and Intelli-Stirrer MSH-300i, Magnetic Stirrers with hot plate

	MSH-300	Intelli-Stirrer MSH-300i
Speed control range	250-1,250 rpm	100–1,250 rpm (10 rpm increment)
Max. stirring volume (H <sub>2</sub> O)	15 L	20 L
Plate temperature regulation range	+30°C +330°C	+30°C +330°C (1°C increment)
Temperature control range with external probe	-	20 °C +150°C
Display	-	LCD
Temperature uniformity on the plate	±3	°C
Working plate heating time till 330°C	15 min	11 min
Diameter of working plate	160 mm	
Plate material	Aluminium alloy	
SR-1, attachable stand size	Ø 8×3	20 mm
Length of magnetic stirring element	10–50 mm 20–70 mm	
Max. stirring liquid viscosity	up to 1,1	70 mPa.s
Maximum continuous operation time	16	8 h
Fault indication	Outputs sound signal and turns off the heating	Outputs an error code on display, turns off the heating
Overall dimensions (W×D×H)	190 × 270 × 100 mm	
Weight	2.9 kg 3.2 kg	
Nominal operating voltage	230 V; 50/60 Hz c	or 120 V; 50/60 Hz
Power consumption (Stirring)	8.5 W	
Power consumption (Heating)	550	) W

Connecting external probe to the Intelli-Stirrer MSH-300i



#### Plate heat up time for MSH-300:

from 25°C — 15 min to 330°C

#### Plate heat up time for Intelli-Stirrer MSH-300i:

from 25°C 11 min to 330°C

#### **EXTERNAL SENSOR INSTALLATION:**



ORDERING INFORMATION:	Cat. number	<u></u>
MSH-300 with stand	BS-010302-OAA	
Intelli-Stirrer MSH-300i with stand	BS-010309-AAA	
Optional accessories:		
External temperature probe	BS-010309-BK	
HTP-1, holder for temperature probe	BS-010309-FK	

DESCRIPTION

1 Rotation:

#### MM-1000, Overhead Stirrer Multi Mixer

Overhead Stirrer Multi Mixer **MM-1000** is designed for stirring liquids up to 20 litres. Quiet and reliable mixer can provide stable continuous mixing up to seven days. It can realise three types of motion:

- 1 Rotational
- 2 Reciprocal
- 3 Vibration.

**MM-1000** performs separate (mono–) (1; 2; 3), consecutive binary cycles (c) (1–2) × c; (1–3) × c and (2–3) × c and complex tri-cycles (1–2–3) × c.

Speed, angle and time of stirrer rotation are under microprocessor control.

Multi Mixer can be used for stirring solutions up to the "medium viscosity" range (from 1,000 to 10,000 mPa.s). It is an ideal instrument for biotechnology, organic synthesis, analytical laboratories.


The innovative combination of three motion types provides a high level of homogeneity due to consecutive combination of laminar and turbulent flows that cause substances to dissolve faster.

Electrically safe and energy efficient – powered by 12 V external power supply.

#### Specifications of movement types:

Speed regulation range	40-1,000 rpm
Time	0-250 s
2 Reciprocal motion:	
Turning angle	0°-360° (increment 30°)
Time	0-250 s
3 Vibro motion:	
Turning angle	0°-5° (increment 1°)
Timer	0-5 s
Timer sound signal	yes
Stirring volume up to (H <sub>2</sub> O)	20 L
Digital time setting	1 min–96 h/non-stop (increment 1 min)
Overall dimensions (W×D×H)	140 × 135 × 250 mm
Weight	2.4 kg
Input current/power consump	tion 12 V, 700 mA / 8.4 W
External power supply	Input AC 100–240 V 50/60 Hz, Output DC 12 V






Multi mixing



#### Accessories for MM-1000

### ORDERING INFORMATION: Cat. number MM-1000 without stirrers BS-010306-AAH

Optional accessories: Dimensions Cat. number Type **⋒** MP-1 378 × (70 × 70) × 8 mm BS-010306-AK Paddle stirrer **B** MP-2 Propeller stirrer 2 folding blades (326 × 55 × 8 mm) BS-010306-BK **MP-3** Propeller stirrer 3 folding blades (325  $\times$  50  $\times$  8 mm) BS-010306-CK MA-1 Anchor stirrer 332 × 90 × 8 mm BS-010306-DK MC-1 BS-010306-EK Centrifugal stirrer 358 × 60 (110) × 8 mm Double clamp For device mounting BS-010306-LK





#### **BIOSAFETY EQUIPMENT:**

# BIOSAFETY AIR, BIOSAFETY SURFACE, WATER PURIFICATION SYSTEMS



#### UVR-M and UVR-Mi, UV Cleaner-Recirculators

#### How does UV-Air Flow Cleaner-Recirculator work?

Operation principle is based on a constant, forced air circulation through recirculator's chamber in close vicinity to UV lamps, thus ensuring maximal efficiency of disinfection. The inner mirror surface of the recirculator chamber reflects ultraviolet rays thereby increasing the UV radiation density and enhancing the disinfection effect.

#### What does UV Air Flow Cleaner-Recirculator consist of?

UV Air Flow Cleaner-Recirculator consists of a germicidal UV lamp, a fan unit equipped with dust filters and a control unit confined in a flow-through chamber.

#### What are the Benefits of UVR-M and UVR-Mi recirculators?

- UV Air Recirculators are ideal for air disinfection in hospitals (especially in outpatient departments, operating rooms, emergency rooms, delivery rooms etc.), kindergartens, research laboratories, veterinary clinics
- Recirculators are effective against common airborne diseases by disinfecting the air and efficiently destroying disease-causing agents (viruses, microorganisms) by UV radiation
- Provide complete protection from UV radiation
- Easy to install, operate and maintain. Very low noise level
- Built-in timer allows controlling the UV lamp operating time (UVR-Mi model)
- Digital control unit allows tracking overall UV lamp operating time (UVR-Mi model)


# Product Class Product Class UVR-M UVR-Mi



UVR-Mi

#### Recirculator fixation:

Convenient fixation on walls (standard)
 Mounting on a movable stand (optional)




# DESCRIPTIC

#### **UVR-M** and **UVR-Mi,** UV Cleaner–Recirculators

	UVR-M	UVR-Mi
UV radiation source bactericidal UV-C, TUV 25W 1SL/25	1 lamp	2 lamps
UV radiation level	18 mW/cm <sup>2</sup> /s	36 mW/cm <sup>2</sup> /s
Air-flow productivity	21 m³/h	29 m³/h
Full user protection from direct UV light	Yes	
Display	-	LCD
UV lamp operation indicator	Yes	Yes
UV lamp lifetime counter	-	Yes
Timer	-	1 min-24 h/non-stop
Clock and ON/OFFscheduler	-	Yes
Lamp fault detection	-	Yes
Overall dimensions (W×D×H)	110 × 135 × 660 mm	110 × 135 × 660 mm
Weight	2.7 kg	2.4 kg
Nominal operating voltage	230 V, 50 Hz or 120 V, 60 Hz	230 V, 50 Hz
Power consumption (230/120 V)	125 VA (540 mA)/160 VA (1.3 A)	110 W (0.5 A)

#### **OPERATION PRINCIPLE**



 $T \ \ - \ \text{thymine formations}$ 



ORDERING INFORMATION:

Cat. number

BS-040105-AAA

UVR-Mi BS-040110-AAA

Optional adapters:

UVR-M

**UVR-S** (stand) BS-040105-AK

UVR-M and UVR-Mi, UV Air Recirculators Test Report. biosan.luvr-test





DNA/RNA UV-cleaner cabinets (UVC/T-AR, UVC/T-M-AR, UVT-B-AR, UVT-S-AR) are designed for dependable, contamination-free operations with DNA samples, supporting stable results in any laboratory environment.

All models are bench-top units, constructed with a robust metal framework, glass (or plexiglass) walls, and a working surface finished in powder enamel or stainless steel (see specifications on page 102). This solid construction ensures durability and a consistently safe workspace as well as makes the transportation easy and reliable.

To minimize contamination risks, PCR cabinets are equipped with an open UV lamp in the upper hood, providing reliable disinfection of the surface by inactivating DNA/RNA fragments within 15–30 minutes. A digital timer ensures the exact needed time for effective and direct UV irradiation.

A new addition is the touchscreen control panel, which allows straightforward cabinet operations, programming of UV exposure cycles, access to the User Manuals, Service access and real-time monitoring of lamp life—making daily operation easy, precise, consistent, and reproducible. For workspace illumination, an energy-efficient 12W LED daylight lamp is now standard. This delivers bright, stable light for your tasks, while reducing energy consumption by around 50% compared to traditional fluorescent lighting—supporting long-term, cost-effective operation.

Each PCR-cabinet includes a flow-type bactericidal UV recirculation unit for continuous decontamination during use. The recirculator, comprising a UV lamp, fan, and dust filters in a dedicated housing, protects the user from direct UV exposure and increases the effectiveness of DNA/RNA inactivation. It processes 100 cabinet volumes per hour, maintaining steady aseptic conditions inside the working area.

Biosan UV-cleaner cabinets are recommended for operations with DNA/RNA, providing a stable and trustworthy environment so you don't have to worry about direct and cross-contamination.



Development and evaluation of DNA amplicon quantification video is available on the website

#### Advantages of Biosan UV-cleaner boxes:

- · Ozone free high-density UV decontamination;
- · Long living UV lamps (9,000 hours average);
- Automatic switching off of UV-lamps when the protective screen is opened;
- Bactericidal flow-type recirculator providing permanent decontamination inside UV-cleaner box during operation;
- · Shockproof glass walls;
- · Low noise, low energy consumption;
- Tables for installation of UV-cleaner boxes;
- UV-cleaner boxes with the bactericidal UV cleanerrecirculator AR is the invented by Biosan solution.














#### Germicidal, shortwave (254 nm) ultraviolet energy is used for complete destruction of various biological agents




#### Average dosage for different surfaces

Surface	Dosage after 15 min	Dosage after 30 min
Working surface (40–50 cm)	570-680 mW/cm <sup>2</sup>	1140-1360 mW/cm <sup>2</sup>
Side walls (10–50 cm)	570-2500 mW/cm <sup>2</sup>	1140-5000 mW/cm <sup>2</sup>
Front window (10–50 cm)	570-2500 mW/cm <sup>2</sup>	1140-5000 mW/cm <sup>2</sup>


ORDERING INFORMATION:	Cat. number
UVC/T-AR with inlet	BS-040102-AAA
UVT-B-AR with internal socket and inlet	BS-040109-A06
UVC/T-M-AR with internal socket and inlet	BS-040104-A06
UVT-S-AR with internal sockets and inlet	BS-040107-AAA



Model	UVT-S-AR (double size)
Wall materials	Rear: stainless steel. Sides and front: glass (EUROGLASS,Germany)
Working surface material	Stainless steel
Open UV-lamp	2×30 W built-in bactericidal UV-C, TUV 30W 1SL/25
Recirculator UV radiation level	18 mW/cm²/s
Radiation type	UV ( $\lambda$ = 253.7 nm), ozone-free
Digital display & Lamp life timer	yes
Display	3.5", medical glove compatible
Digital time setting of direct UV exposure	1 min–24 h/non-stop (increment 1 min)
UV–recirculator	1 × 30W (efficiency >99% per 1 h)
Daylight lamp (for working area illumination)	1 × LED-12W
Thickness of side panels	4 mm
Thickness of upper front panel	8 mm
Thickness of the front protective screen	5 mm
Optical transmission	95%
UV protection	>96% UV-protection film
Working area dimensions	1,210 × 500 mm
Opening size (W×H, fully raised protective screen)	1,210× 190 mm
Safety features	Automatic open UV-lamp switching off when screen is open
Power outlets inside the unit (230/120 V)	3 built-in sockets max. 1000 W/600 W, Inlet for power cords
Nominal operating voltage	100-240 V, 50/60 Hz
Power consumption	113 W
Overall dimensions(W×D×H)	1250 × 600 × 590 mm
Weight (net/gross)	51/62 kg
Optional table	<b>T-4L</b> (W×D×H: 1290×600×770 mm)









#### UVC/T-AR(compact)

Plexiglass: Polymethyl methacrylate (ALTUGLAS EX)

UVC/T-M-AR (compact)

Rear: stainless steel. Sides and front: glass (EUROGLASS, Germany) **UVT-B-AR** (compact)

Rear: stainless steel. Sides: steel with chemical resistant powder coating. Front: glass (EUROGLASS,Germany)

Steel with chemical resistant powder coating

Stainless steel

1 × 25 W built-in bactericidal UV-C, TUV 25W 1SL/25

18 mW/cm<sup>2</sup>/s

UV ( $\lambda$  = 253.7 nm), ozone-free

yes

3.5", medical glove compatible

1 min-24 h/non-stop (increment 1 min)

1 × 25W (efficiency >99% per 1 h)

	1 × LED-8W	
4 mm	4 mm	2 mm
	8 mm	
8 mm	4 mm	4 mm
92%	95	5%
>99.90% Polymethyl methacrylate ALTUGLAS EX	>96% UV-protection film	

methacrylate ALTUGLAS EX

•

645 × 490 mm

645 × 165 mm 645 × 190 mm

Automatic open UV-lamp switching off when screen is open


Inlet for power cords


Inlet for power cords and 1 built-in socket, max. 1,000 W/600 W

100-240 V, 50/60 Hz

64 W

720 × 535 × 555 mm	700 × 580	× 555 mm
21.4/30 kg	27,4/36 kg	28,8/37.5 kg
<b>T-4</b> (W×D×H: 800×600×745 mm)		





#### P-5, F-1, Shelves for DNA/RNA UV-cleaner boxes

Two types of shelves have been developed for DNA/RNA UV-cleaner boxes to increase the effective area of the box: **P-5** – shelf-holder for five pipettes and **F-1** flat shelf. On the **F-1** shelf, you can place laboratory glassware, reagents and other items that are convenient to keep in close proximity.

P-5, shelf for pipettes	BS-040104-DK
ORDERING INFORMATION:	Cat. number
A B	
Dimension (W×D)	400 × 140 mm
<b>B</b> F-1, shelf:	
Capacity	5 pipettes
Dimension (W×D)	230 × 140 mm
<b>A</b> P-5, shelf for pipettes:	





#### PDS-250 and PDS-10L, DNA/RNA decontamination solution

BS-040104-CK

Contamination is especially problematic in the highly sensitive PCR technique. Originating from aerosolized fragments, contaminant DNA can lead to cross-contamination, thus resulting in inaccurate data and, as a result, misinterpreted analysis.

F-1. shelf

**PDS** is a ready-to-use solution for eliminating DNA and RNA from the surface prior PCR reaction preparation. DNA/RNA is removed within seconds after use. The solution contains a non-alkaline and non-carcinogenic agent. **PDS** is intended for use at PCR cabinets and laminars (e.g. UVT-S-AR), lab devices – BioMagPure 12, TS-100, pipettors – Assist series pipettes etc.

**PDS** is effective against amplicon, plasmid, or genomic DNA and RNA from most surfaces except light or non-ferrous metals (e.g. aluminium, copper, lead, nickel, tin, titanium, zinc etc.)

The use of **PDS** both before and after PCR analysis is fast, easy and ideal to maintain a clean work area, thereby saving time and expenses.

PDS is heat resistant and stable for several years.

The decontamination solution is also available in 10 l containers – **PDS-10L.** 





ORDERING INFORMATION:

Cat. number

#### Laboratory furniture







**G** LF-1, laboratory chest of drawers

ORDERING INFORMATION:	Cat. number
<b>T-4</b> , table	BS-040101-BK
<b>T-4L</b> , table	BS-040107-BK
<b>LF-1</b> , laboratory chest of drawers	BS-050101-BK

Modular design of laboratory furniture provides flexibility and ease of use.

A T-4, table for – UVC/T-AR, UVC/T-M-AR, UVT-B-AR	
Maximum load	50 kg
Drawers	1
Mobility	Wheels with brakes
Material	Laminated particle board
Overall dimensions (W×D×H)	800 × 600 × 745 mm
Weight	23 kg

<b>B</b> T-4L, table for – UVT-S-AR	
Maximum load	75 kg
Drawers	1
Mobility	Wheels with brakes
Material	Laminated particle board
Overall dimensions (W×D×H)	1290 × 600 × 770 mm
Weight	36 kg

<b>C</b> LF-1, laboratory chest of drawers	
Drawers	5
Mobility	Wheels with brakes
Material	Laminated particle board
Overall dimensions (W×D×H)	300 × 450 × 705 mm
Weight	28 kg

UVT-S-AR on T-4L table with two LF-1 laboratory chests



#### LT-120, LT-150 and LT-180, Height adjustable desks

LT-120, LT-150 or LT-180 laboratory table is the perfect solution for any research or laboratory facility. It features a chemically resistant high pressure laminate worktop, making it suitable for use in a variety of laboratory settings. The worktop is easy to clean and maintain, ensuring a hygienic work environment. The table is also adjustable in height, which allows you to customize the table to your preferred working height for added comfort and convenience. In addition to its standard features, this laboratory table also offers an optional drawer with an anti-slippery surface and a retractable block of electrical sockets. This added feature allows for convenient storage of lab equipment and tools, and the anti-slippery surface ensures that items stay securely in place. The retractable block of electrical sockets allows for easy access to power outlets, eliminating the need for extension cords and reducing clutter on the work surface. This table is an all-in-one solution for any laboratory or research facility, providing a sturdy, reliable and practical work surface that can stand up to the demands of daily use.

- · 2 legs with 2 motors;
- Controller with memory 4 memory modes, digital height indicator, button with up/down movement;
- · Adjustable height: 625-1280 mm;
- · Lifting capacity: 89/86/83 kg;
- Work surface from chemically resistant HPL, width 1200/1500/1800 mm, depth 600 mm, thickness of 25 mm, blue edge.

Drawer and retractable socket block comes as an accessories.



S	A LT-120, Height adjustable desk				
IONS	Maximum load	89 kg			
FICAT	Overall dimensions (W×D×H)	1200 × 600 × 625–1280 mm			
PECI	Weight	45 kg			
S					

LT-150, Height adjustable desk				
Maximum load	86 kg			
Overall dimensions (W×D×H)	1500 × 600 × 625–1280 mm			
Weight	49 kg			

LT-180, Height adjustable desk				
Maximum load	83 kg			
Overall dimensions (W×D×H)	1800 × 600 × 625–1280 mm			
Weight	52 kg			

**LTD-1**, drawer with non slip silicone mat and soft close function. Can be mounted afterwards.

<b>D</b> LTD-1, Drawer	
Drawer internal dimensions (W×D×H)	310 × 350 × 40 mm
Overall dimensions (W×D×H)	430 × 410 × 70 mm

SBR-1, socket block, retractable, with blue LED.

- 3× power sockets
- 2× USB-A ports, 5V, 2.1A

If needed, should be indicated before placing an order.

<b>SBR-1</b> , Socket block	
Overall dimensions (W×D×H)	120 × 120 × 230 mm

ORDERING INFORMATION:	Cat. number	$\overline{}$
LT-120, Height adjustable desk	BS-040107-EK	
LT-150, Height adjustable desk	BS-040107-JK	
LT-180, Height adjustable desk	BS-040107-HK	
Optional		
LTD-1, Drawer	BS-040107-AAA	
SBR-1, Socket block	BS-040107-NK	

#### Ultrapure water systems: Labaqua series





#### The Labaqua systems include:

- · Boost pump
- Pre-filter set
- · Reverse osmosis module
- · Deionization module
- · Final stage polishing module
- 30 L storage tank with an integrated Grade 2 dispensing valve
- · Recirculation system

#### Model specific modules:

- Labaqua Trace Point-of-use microfilter
- Labaqua HPLC Point-of-use microfilter, TOC monitor
- Labaqua Bio Point-of-use ultrafilter,
   UV sterilization module, TOC monitor

Labaqua ultrapure systems are multi-purpose water purification systems. The Labaqua systems produce ultrapure and pure water directly from tap water.

Ultrapure (Grade 1) water is dispensed through the point-of-use filter on the front panel. Pure (Grade 2) water is dispensed directly from the storage tank.

Labaqua ultrapure water can be used for the most demanding applications, including, but not limited to: Inorganic trace analysis, Liquid chromatography, Cell culture, Molecular biology. With resistivity of 18.2 Mega – Ohm × cm (0.055  $\mu$ S/cm), ultrapure water produced by a Labaqua system exceeds requirements of all relevant standards (ISO 3696 Grade 1, ASTM Type I, CLSI Type I). Purified water is collected in a storage tank. An integrated recirculation system ensures consistent quality of water and reduces total organic carbon (TOC) to very low levels: <5ppb. Pure water produced by the Labaqua systems complies with ISO 3696 Grade 2 water requirements and can be used for labware washing, wet chemistry methods, flame spectrophotometers, etc

All cartridges and filters are easily accessible, and no tools are required to replace them. The Labaqua system can be installed on a laboratory bench or mounted on a wall.

#### FEATURES:

- Volumetric dispense enables the user to set accurate dispensing volume for each dispense cycle. The dispense volume can be set either from the keyboard or by using the "teaching" mode.
- Water quality embedded recirculation loop ensures stable premium water quality and enables practical elimination of Total Organic Carbon (TOC).
- Low running costs performance of the deionization and polishing modules is constantly monitored. Monitoring algorithm enables cutting running costs, as replacement of the modules is requested only when service life is close to the end.
- Total organic carbon (TOC) monitor organic contaminants may not affect the conductivity of water, so conductivity sensors cannot be used for TOC monitoring. Therefore, a special TOC monitoring module is needed to measure TOC level.
- Color graphic LCD display system component status is reflected on the display in an intuitive colour pattern (Green/ Yellow/Red).
- System flowchart shows all component status and water quality parameters at a glance.

## Ultrapure water systems: Labaqua series

Purified water specifications	Labagua Trace	Labaqua HPLC	Labagua Bio
Ultrapure (Grade 1) water resistivity	18.2 MΩ × cm		
Ultrapure (Grade 1) water conductivity	0.055 μS/cm		
Pure (Grade 2) water resistivity		>10 MΩ × cm	
Pure (Grade 2) water conductivity		<0.1 µS/cm	
TOC	<10 ppb	<5;	opb
RNase	i i		<0.01 ng/ml
DNase	-	-	<4 pg/ml
Bacteria	<0.01 CFU/ml		
Endotoxins	<0.15 EU/ml <0.001 EU/ml		
Particles >0.22 μm		<1/ml	
Deionization module life (standard module)		1 m³	
Dimensions (W×D×H)	320 × 560 × 620 mm		
Storage tank	30 I		
Feed water pressure	0.5–5 bar		
Feed water conductivity	<1,300 µS/cm		
Weight	24 kg 25 kg 26 kg		26 kg
Nominal operating voltage	100-240 V, 50/60 Hz		
Power consumption	130 W		

	Application	Labaqua Trace	Labaqua HPLC	Labaqua Bio
	Glassware rinsing	+	+	+
	Laboratory washers	+	+	+
	Autoclaves	+	+	+
General	Electrochemistry	+	+	+
laboratory applications	Wet chemistry	+	+	+
аррисацонз	Spectrophotometry	+	+	+
	Buffer and media preparation	+	+	+
	Reagent preparation	+	+	+
	Flame atomic absorption spectrophotometry	+	+	+
Inorganic	Graphite atomizer atomic absorption spectrophotometry	+	+	+
analysis	Plasma mass-spectrometry (ICPMS)	+	+	+
methods	Plasma spectrophotometry (ICPOES)	+	+	+
	Ion chromatography	+	+	+
Organic	Liquid chromatography (HPLC/ UHPLC)		+	+
analysis	Gas chromatography		+	+
methodes	Total organic carbon measurements		+	+
Mala da	Flow cytometry			+
Molecular Biology	Cell and tissue culture			+
Diology	Molecular biology			+


Cat. number

410222

## Ultrapure water systems: Labaqua series










ORDERING INFORMATION

Filter set for BS-070104-KK (Carbon/PP, PP  $1\mu m$ )



Labaqua Trace include 30 l tank, power cord	BS-070105-A02
Labaqua HPLC include 30 l tank, power cord	BS-070104-A02
Labaqua Bio include 30 l tank, power cord	BS-070106-A02
Water purification system sets including different storage tanks are available	please inquire
Optional accessories:	
External pre-filter set (polyphosphate/carbon/1 μm) with manometer	BS-070104-LK
External pre-filter set (carbon/1µm) with manometer	BS-070104-KK
Storage tank "Comfort" with base, tap and multipoint level switch, 60 l	BS-070102-SK
Storage tank "Comfort" with base, tap and multipoint level switch, 100 l	BS-070102-FK
Universal remote dispenser set with 3 m supply hose and water distribution module	BS-070104-JK
External input pressure reducer	10175
Replacement parts	
Internal prefilter set	BS-070104-AK
Deionization module	BS-070104-IK
Polishing module	BS-070104-BK
RO membrane	BS-070104-NK
Point-of-use microfilter – 0.22 μm non sterile	BS-070104-EK
Point-of-use microfilter – 0.22 μm sterile	BS-070104-FK
Point-of-use ultrafilter	BS-070104-GK
UV bulb 254 nm	BS-070104-RK
UV bulb 185 nm	BS-070104-DK
0.22 μm air vent filter for the storage tank	BS-070102-AK
Air filter for storage tank	BS-070104-PK
Filter set for BS-070104-LK (Polyphosphate, Carbon/PP, PP 1µm )	410223



# DENSITOMETERS, PHOTOMETER



### **DEN-1** and **DEN-1B**, McFarland Densitometers

Densitometers **DEN-1** and **DEN-1B** are designed for measurement of cell suspension's turbidity in the range:

0.0-6.0 McFarland units

 $(0-180 \times 10^7 \text{ cells/ml})$ :

Densitometers provide the opportunity to measure solution turbidity in a wider range (up to 15.0 McFarland units), however, it is necessary to remember that, in this case, the standard deviation values increase.

A densitometer is used for measurement of cell concentration (bacterial, yeast cells) during the fermentation process, determination of microorganism sensitivity to antibiotics, microorganism identification using various test-systems, for measurement of absorption at the definite wavelength, as well as for quantitative estimation of the colour solution concentration, absorbing green light.

The operation principle is based on the measurement of optical density with digital presentation of results in McFarland units. The unit is calibrated at the factory (for operation with 16 mm diameter glass tubes) and keeps calibration without power supply. However, if necessary, it is possible to calibrate the unit by 2–6 points in 0.0–6.0 McFarland unit range. Both commercial standards offered by Biosan and the cell suspensions prepared in a laboratory can be used for calibration.

## Following polymer microparticles calibration kits and glass tubes are available on request:

- **CKG16** for glass tubes with diameter 16 mm, set of 0.5; 1.0; 2.0; 3.0; 4.0 McFarland Turbidity Standards (latex particles)
- Calibration kit for glass tubes with diameter 18 mm, set of 0.5; 1.0; 2.0; 3.0; 4.0; 5.0 McFarland Turbidity Standards (BaSO4)
- Calibration kit for glass tubes with diameter 12 mm, set of 0.0 (blank); 0.5; 2.0; 3.0 McFarland Turbidity Standards (latex particles)
- Glass Test Tube Ø16 or Ø18 × 100 mm, high borosilicate, PP Cap with silicone pad.
   Packing – 100 pcs/box

Up to date information on calibration kits can be found on the website: www.biosan.lv

#### Two versions of the product are available:

- 1. DEN-1 powered from external energy supply;
- DEN-1B powered both from external energy supply and batteries (AA).







DEN-1B rear side with calibration controls




## **DEN-1** and **DEN-1B**, McFarland Densitometers

	DEN-1	DEN-1B
Light source	LED	
Wavelength	λ = 565 ±1	15 nm
Measurement range	0.00-15.0	0 McF
Display resolution	0.01 M	lcF
Accuracy	(0.5–4.0 Mc	cF) ±3%
Measurement time	1 s	
Sample volume	not less tha	an 2 ml
Tube external diameter	18 mm (withou 16 mm (using include 12 mm (using option	ed A-16 adapter)

Possibility to restore factory calibration settings

Display	LCD	
Overall dimensions (W $\times$ D $\times$ H)	165×115×75 mm	
Weight	0.7	′ kg
Independent power supply	_	3 × AA batteries
Input current/power consumption	12 V, 7 n	nA/0.1 W
External power supply	Input AC 100–240 V, 50	V60 Hz, Output DC 12 V
Standard set	Adapter A-16 External power supply	Adapter A-16 External power supply







ORDERING INFORMATION:	Cat. number	$\overline{}$
DEN-1 with A-16 adapter	BS-050102-AAF	
DEN-1B with A-16 adapter	BS-050104-AAF	
Optional accessories:		
A-12 adapter for 12 mm tubes	BS-050102-IK	
<b>CKG16</b> for glass tubes with diameter 16 mm, set of 0.5; 1.0; 2.0; 3.0; 4.0 McF	BS-050102-BK	
Calibration kit for glass tubes with diameter 18 mm, set of 0.5; 1.0; 2.0; 3.0; 4.0; 5.0 McF	70900	
Calibration kit for glass tubes with diameter 12 mm, set of 0.0; 0.5; 2.0; 3.0 McF	21255	
<b>Glass Test Tube</b> Ø16 × 100 mm, high borosilicate, PP Cap with silicone pad. Packing – 100 pcs/box	BS-050102-MK	
$\textbf{Glass Test Tube} \ \emptyset 18 \times 100 \ \text{mm, high borosilicate, PP Cap with silicone pad. Packing - 100 \ \text{pcs/box}$	BS-050102-NK	

Light source

## **DEN-600** Photometers

**DEN-600** is a compact, portable, rechargeable battery-powered photometer. It comprises of 600 nm wavelength optical system, which enables to apply:

- 1. OD600 method estimates the total number of cells.
- 2. McFarland (McF) turbidity measurement method.
- Bradford protein assay method for protein concentration measurement.
- 4. Other methods that can be adjusted or optimized using 600 nm wavelength.

The device serves as an inexpensive alternative to a spectrophotometer commonly used for these applications. Because **DEN-600** is battery powered and compact, it can be comfortably located in a biosafety cabinet, anaerobic chamber or quickly moved to another lab room. Additionally, the vessel holding mechanism allows accommodating round bottom, conical vials or falcon tubes, therefore enabling to measure the absorbance (Abs) and turbidity in Abs, OD and McFarland units.

USB connectivity and DEN software allow for data transfer, data processing and calculation, software calibration for Bradford protein assay method or a custom calibration for a specifically applicable vessel.


Light 30dice	LLD
Photodetector	Silicone photodiode
Measurement wavelengt	th (λ) 600 nm ±10 nm
Vessel type	Cuvettes, round bottom tubes, falcon tubes
Battery type	LiPo
PC system requirements:	Intel/AMD Processor, 1 GB RAM, Windows 10/11, USB
Dimensions (W×D×H)	120 × 145 × 65 mm
Weight	0.5 kg
External power supply	Input AC 100–240 V 50/60 Hz, Output DC 12 V
Input current / power co	nsumption 12 V, 0,2 A / 2,5 W



#### **COMMON APPLICATIONS:**

- · Cell concentration measurement
- · Cell growth data estimation
- · Log phase estimation for microbial cells induction
- · Competent cell preparation
- · Bradford protein assay method
- · Antibiotic susceptibility testing
- · Inhibitory tests

LED





Measurement modes	Absorbance	McFarland
Measurement range	0-3.0 Abs	0-16.00 McF
Resolution	0.001 Abs	0.01 McF
Accuracy	±0.006 @ 1 Abs	±0.1 @ 0-8 McF
Repeatability	±0.003 @ 1 Abs	±0.05 @ 0-8 McF

#### ORDERING INFORMATION

Cat. number

#### Optional accessories:

**DEN-600** 

**Verification set for Abs.** Certified reference material, neutral density glass filter set of 4 Abs reference points – 0.3532, 1,0512, 2,0425, 2,927 (the values may vary slightly from batch to batch)

BS-050109-AK

CKG16 for glass tubes with diameter 16 mm, set of 0.5; 1.0; 2.0; 3.0; 4.0 McF

BS-050102-BK

Calibration kit for glass tubes with diameter 12 mm, set of 0.0; 0.5; 2.0; 3.0 McF

21255

Glass Test Tube 16 × 100 mm, high borosilicate, PP Cap with silicone pad. Packing – 100 pcs/box

BS-050102-MK

## ASPIRATORS, PIPETTES



**FTA-2I**Aspirator with Trap Flask

DESCRIPTION

## FTA-2i, Aspirator with Trap Flask

Aspirator with trap flask **FTA-2i** is designed for aspiration or removal of alcohol, buffer and liquid from reaction vessels (e.g. during DNA/RNA purification or other macromolecule reprecipitation techniques).

The device can be applied for routine operations of cells washing from culture medium and resuspension in a buffer. Aspirator operation principle is based on creating negative pressure in trapping flask using built-in microcompressor. The collecting tip is connected with polyethylene tube to the trapping flask. Liquid is removed from the reaction vessel when the collecting tip is in contact with the solution. A tube holder-organizer is conveniently located at **FTA-2i** right-hand side; it accommodates two 1.5–2 ml tube slots (e.g. for hydrochloric acid solution and distillate) necessary for collecting tip washing and storing, so that a tip can be re-used.

**FTA-2i** is equipped with a level sensor that detects excess liquid with consequential prevention of the overflow by automatically switching off the pump with a sounding alarm indication.

The devices come, as standard, with a vacuum regulation control knob that allows to select a preferable aspiration speed smoothly.

Additionally, a hand operator can be purchased for a more comfortable usability of the new accessories. The autoclavable hand operator features a pressure-sensitive button that can control the aspiration speed.

#### **Common applications:**

Aspiration speed

Removal and disposal of liquid from various reaction vessels

up to 10 l/min (air)

, ispiration speed	ap to 10 11 11 (all)
Vacuum regulation	-200 to -800 mbar (adjustable)
Trap flask	2 l, polypropylene (autoclavable)
Liquid level sensor type	Invasive
Timer sound signal	yes
Overflow protection	Motor stops, light and sound signal
,	microbiologic filter 2200/02 mination from the trap flask infected particles
Filter pore diameter	0.027 micron
Input current/power consumption	12 V, 1 A/10.8 W
External power supply	Input AC 100–240V 50/60 Hz; Output DC 12 V
Dimensions (W×D×H)	185 × 290 × 390 mm
Weight*	1.85 kg

<sup>\* -</sup> Accurate within ±10%.




## FTA-2i, Aspirator with Trap Flask

#### **Optional accessories:**

HAS-1, hand operator set

- Handheld vacuum controller;
- 1-channel adapter (with ejector) for 200 μL tips;
- 3 1-channel adapter with 125 mm stainless steel pin;
- 4 1-channel adapter with 40 mm stainless steel pin;
- **3** 8-channel adapter (with ejector) for 200 μL tips;
- 6 8-channel adapter with 35 mm stainless steel pin;
- **7** Stand SH-6.



ORDERING INFORMATION:	Cat. number
<b>FTA-2i</b> , with 2l trap flask, universal adapter <b>MA-U</b> (for 200/1000 $\mu$ L single use tips)	BS-040120-A02
Optional accessories:	
HAS-1, hand operator set	BS-040118-PK
MA-8, 8-channel adapter manifol	BS-040108-BK
Extended tubing 2 m long, with fittings and MA-U adapter	BS-040120-DK
Replacement parts:	
Suction microbiologic hydrophobic filter	BS-040120-S10
<b>MA-U</b> , universal adapter for 200/1000 μL single use tips	BS-040118-AK

DESCRIPTION

## FTA-U, Universal vacuum aspirator NEW

**FTA-U** is designed for the precise aspiration and removal of alcohol, buffer solutions, and other liquids from reaction vessels, including applications such as DNA/RNA extraction, ELISA plate washing, and cell culture media exchange. The device is suitable for use with microtubes, plates, and flasks, and can be configured with either a 2 L or 4 L polypropylene trapping flask.

Suction power is regulated continuously between –200 mbar and –950 mbar, with a maximum aspiration speed of 12 l/min (air) and a suction flow rate up to 47 ml/s (liquid). Aspiration intensity is controlled by a linear regulator on the front panel, and operation status is clearly indicated with a light ring around the control knob. Overflow protection is provided by an invasive level sensor, which stops the pump and triggers sound and light alarms to prevent accidental spills.


**FTA-U** is equipped with a hydrophobic microbiological air filter (pore size 0.027 µm) that removes up to 99.99% of bacteria and viral particles from the airstream, ensuring contamination-free operation and laboratory safety. The standard set includes the MA-U universal adapter for singleuse 200 µL or 1000 µL tips, with optional accessories such as an 8-channel aspiration tip (MA-8), HAS-1 hand operator for ergonomic handling, and a foot switch for hands-free operation.

All parts in contact with liquids are chemically resistant and suitable for use with common laboratory reagents. The trapping flask, lid, and fittings are autoclavable at 121°C for 15 minutes. Biosan recommends replacing the hydrophobic filter every 30 days and silicone tubing every six months for optimal performance.


optimal performance.	
Aspiration speed	up to 12 l/min (air)
Vacuum regulation	-200 to -950 mbar (-95 kPa)
Trap flask	2–4 L, polypropylene (PP)
Liquid level sensor type	Invasive
Timer sound signal	yes
Noise level , dBA	100%/62% vacuum ~ 56/ 46 db
Microbiological filter	hydrophobic air filter, pore diameter 0.027 μm
Pump type	Durable mini diaphragm vacuum pump
Quick-coupling system	+
Overflow protection	Pump stop, light and sound signal
Filtration: Hydrophobic mi eliminates risk of contamir by bacteria, viruses and inf	nation from the trap flask
Filter pore diameter	0.027 micron
Input current/power consumption	12 V, 2 A/24 W
External power supply	Input AC 100–240V 50/60 Hz; Output DC 12 V
Dimensions (W×D×H)	$220 \times 320 \times 380$ mm (2L flask) $220 \times 320 \times 460$ mm (4L flask)

2.40 kg / 2.82 kg

Weight with 2L flask / 4L flask



## **FTA-U,** Universal Aspirator with Trap Flask



ORDERING INFORMATION:	Cat. number
FTA-U, aspirator with 2 L trap flask	BS-040122-A03
FTA-U, aspirator with 4 L trap flask	BS-040122-A02
Optional accessories:	
HAS-1, hand operator set	BS-040118-PK
MA-8, 8-channel adapter manifol	BS-040108-BK
FS-1, Foot switch	BS-010177-AK
Trap flask set, 2l	On request
Trap flask set, 4l	On request
Replacement parts:	
Suction microbiologic hydrophobic filter	BS-040120-S10
MA-U, universal adapter for 200/1000 $\mu L$ single use tips	BS-040118-AK

External power

supply

## FTA-1, Aspirator with Trap Flask

Aspirator with trap flask **FTA-1** is designed for aspiration/ removal of alcohol/buffer remaining quantities from microtest tube walls during DNA, RNA purification and other macromolecule reprecipitation techniques.

The device can also be used for routine operations of cells washing from culture medium and resuspension in a buffer. Aspirator operation principle is based on creating negative pressure in trapping flask using built-in microcompressor. The collecting tip is connected with polyethylene tube to the trapping flask. Liquid is removed from the microtest tube when the collecting tip touches the solution surface. A tube holder-organizer is conveniently located at FTA-1 right-hand side; it accommodates two tubes (e.g. for hydrochloric acid solution and distillate) necessary for collecting tip washing and storing, so that a tip can be reused.

● Suction microbiological hydrophobic filter type 2200/02: Suction microbiologic filter eliminates the risk of contamination with bacteria, viruses and infected particle from patient to suction pump or central vacuum distribution. Suction microbiological filter is hydrophobic with very high bacterial blocking efficiencies, up to 99.99999% particles bigger than 0.027 µm (which is smaller than Hepatitis A, B and C).

Vacuum	–500 mbar
Trap flask volume	1 L
Dimensions with trap flask (W×D×H)	186 × 213 × 327 mm
Weight with trap flask	1.8 kg
Input current/power consumption	12 V, 300 mA / 3.6 W

Input AC 100-240 V; 50/60 Hz;

Output DC 12 V

BS-040108-S25





Suction microbiologic hydrophobic filter









## Assist, pipette series



The Assist series pipettes are single, 8 or 12 channel variable volume pipettes designed to measure and transfer volumes.

Single-channel pipettes are produced in ten ranges of volumes from 0.1  $\mu$ l to10,000  $\mu$ l depending on the model.

Multichannel pipettes are produced in four ranges of volumes: 0.5– $10 \mu l$ , 5– $50 \mu l$ , 20– $200 \mu l$ , 50– $300 \mu l$ .

The pipettes are equipped with an analogue counter which shows the pipetting volume. The volume setting is done by turning the pipetting pushbutton knob or the black adjustment knob in the right direction. The volume range is shown on the pipetting pushbutton. New versions with volume lock function also available.

Common pipettes usage depending on the volume		
AP2, AP10, AP8-10, AP12-10	Measurement and transfer of micro-volumes, DNA sequencing and enzyme-assay applications.	
AP20, AP50, AP100, AP200, AP250, AP1000, AP8-50, AP12-50, AP8-200, AP12-10, AP8-300, AP12-300	Measurement and transfer of general aqueous solution, acids and bases.	
AP5000, AP10000	Measurement and transfer of large volumes.	


Pipette		Volume (µl)	Colour code	Fit to tips	Cat. nu	umber
No lock	With lock				No lock	With lock
Single chan	inel:					
AP2	APL2	0.1-2.0	•	10	BS-010501	BS-010526
AP10	APL10	0.5–10.0	•	10 µl	BS-010502	BS-010527
AP20	APL20	2–20			BS-010503	BS-010528
AP50	-	5–50		200 µl	BS-010504	-
AP100	APL100	10–100		200 μι	BS-010505	BS-010530
AP200	APL200	20-200			BS-010506	BS-010531
AP250	-	50-250		300 μl	BS-010507	-
AP1000	APL1000	100–1,000		1,000 µl	BS-010508	BS-010532
AP5000	APL5000	500-5,000	$\circ$	5,000 µl	BS-010509	BS-010533
AP10000	APL10000	1,000–10,000		10,000 µl	BS-010510	BS-010534
Multichann	el:					
AP8-10 AP12-10	APL8-10 APL12-10	0.5–10	-	10 μΙ	BS-010511 BS-010512	BS-010535 BS-010539
AP8-50 AP12-50	APL8-50 APL12-50	5–50	-	200 1	BS-010513 BS-010514	BS-010536 BS-010540
AP8-200 AP12-200	APL8-200 APL12-200	20–200	-	200 µl	BS-010515 BS-010516	BS-010537 BS-010541
AP8-300 AP12-300	APL8-300 APL12-300	50–300	-	300 µl	BS-010517 BS-010518	BS-010538 BS-010542
Sets:						
AP10, AP20, AP200, AP1000, 4 position stand, demo tips				BS-010519	_	
AP10, AP100, AP1000, AP5000, 4 position stand, demo tips				BS-010520	-	
APL Starter Kit 1: APL10, APL20, APL200, APL1000, 4 place stand, demo tips 3× boxes – 10 μl, 200 μl, 1000 μl				-	BS-010543	
APL Starter Kit 2: APL10, APL100, APL1000, APL5000, 4 place stand, demo tips 3× boxes – 10 μl, 200 μl, 1000 μl, 5× tips 5000 μl				-	BS-010544	

## Assist, pipette series

#### **Features:**

- · Contoured shape of the handle and light weight;
- · Proven accuracy and precision;
- · UV resistant & fully autoclavable;
- · 5 & 10 ml shaft protected by filter;
- · Available in 8 and 12-channel version;
- Colour coded for easy volume identification;
- The adjustable ejector height system to accommodate virtually all brands of tips;
- Dual volume setting using the pushbutton or the thumbwheel:
- · Soft spring system for smooth, effortless pipetting.
- · Versions with volume lock function also available.

Colour coded for easy volume identification





#### ○ ORDERING INFORMATION:

Pipette stands:	Cat. number		
Carousel stand (rotating) for 6 pipettes	BS-010522		
Multiple stand (fixed) for 8 pipettes	BS-010523		
3 1-position stand	BS-010524		
4 4-position stand	BS-010525		

#### Pipette tips features:

- The tip is made of imported high-quality medical grade polypropylene, and the filter element is made of ultra-high molecular weight polyethylene, which has good hydrophobicity.
- The filter element effectively protects the pipette from the sample to ensure the safety of sample aspiration.
- The filter element can reduce the hazard of the gas residue in the tip body to the operator during the pipetting process and improve the repeatability of the experiment.
- · No inner surface coating, will not pollute the sample
- · Have precise tick marks
- · DNase, RNase and Pyrogen free
- · Can be autoclaved
- Suitable for conventional pipettes such as Biosan (Assist series) Eppendorf and Gilson.
- Provided in sterile racks.



10 μl Filtered/racks



2 200 µl Filtered/racks



3 1000 μl Filtered/racks

#### ₩ OF

#### ORDERING INFORMATION:

ORDERING INFORMATION.					
Pipette tips:	Size (Hר)	Packaging	Color	Compatibility	Cat. number
10 μl Filtered/racks	32 × 0.4 mm	96 tips/rack	Natural	Universal	BHZ01R1W-FS
2 200 μl Filtered/racks	59 × 0.5 mm	96 tips/rack	Natural	Universal	BHZ03R1W-FS
3 1000 µl Filtered/racks	86 × 0.8 mm	96 tips/rack	Natural	Universal	BHZ05R1W-FS



#### Charging stand



## Assistboy, pipette controller

**Assistboy** pipette controller is a device intended for pipetting liquids with the use of measuring pipettes. It can work with all types of glass or plastic serological pipettes in the volume range from 0.5 ml to 100 ml.

Controller is equipped with an exchangeable filter membrane which protects shaft mechanism from aggressive liquid fumes.

Two dispense modes permit selection of dispensing intensity depending on the user's needs. The selected setting of the pipette controller mode is shown on display.

#### Safe and efficient work

- Protected by a PTFE filter blocking any liquid from entering the unit
- Autoclavable filter, the pipette holder and the nosepiece
- · UV resistant body for safe sterilization
- Powerful, environmentally friendly 3 Ni-MH batteries enable many hours of continuous work
- LCD display showing battery charge level

#### Speed and working mode adjustment

- Function buttons for SPEED and working MODE control in a reach of a thumb
- Additional speed adjustment by the pressure applied to the trigger buttons

#### **Working comfort**

- Suitable for glass & plastic volumetric pipettes 0.5–100 ml
- · Ergonomically shaped handle
- · Well located function buttons
- · Convenient charging stand

#### ORDERING INFORMATION:

Cat. number



**Assistboy** with charging stand

BS-010521







## **BIOPROCESSING:**

## CO<sub>2</sub> INCUBATOR, SHAKER-INCUBATORS, PERSONAL BIOREACTORS



RTS-1 and RTS-1C
Personal Bioreactor







Reverse-Spin® Innovative Mixing Technology



USB connection



Product video

Users articles: biosan.lv/report

	RTS-1	RTS-1C
E.coli BL21 Factory calibration measurement range, in $OD_{850}$ : at 10–20 ml volume at 20–30 ml volume	0-10 OD (0-19 0-8 OD (0-15.2	OD <sub>600</sub> equivalent) 2 OD <sub>600</sub> equivalent)
Factory calibration measurement precision	±0.3	OD <sub>850</sub>
Mass transfer coefficient k <sub>L</sub> a (h <sup>-1</sup> )	Up to 350 ±2	26 h <sup>.1</sup> at 5 ml
Measurement Wavelength (λ)	850 ±	15 nm
Light source	LE	ED
Real time measurement	1–60	) min
Temperature setting range	+25°C +70°C (increment 0.1°C)	+4°C +70°C (increment 0.1°C)
Bottom control range point	5°C above ambient	15°C below ambient
Top control range point	70	0°C
Stability	±0.°C	
Sample temperature accuracy: 20–45°C <20°C >45°C	±	.1 2 3
Sample temperature heating/cooling rate	0.7°C/min	
Sample volume	10-30 ml	
Speed control range	50–2,000 rpm (in	crement 10 rpm)
Speed control precision	±15	rpm
Reverse-Spin Time	1–60 s (inci	rement 1 s)
Display	LC	CD .
Minimum PC requirements		Processor, s 10/11, 2.0 USB port
Optimal PC requirements	Intel/AMD Processor, 3 GB RAM, Windows 10/11, 2.0 USB port	
Overall dimensions (W×D×H)	130 × 212 × 200 mm	
Weight	1.7 kg	2.2 kg
Input current/power consumption	12 V DC, 3.3 A/40 W	12 V DC, 5 A/60 W
External power supply	Input AC 100–240 V 50/60 Hz; Output DC 12 V	

<sup>\* –</sup> Highest k<sub>L</sub>a (h<sup>-1</sup>) is achieved at 5 ml working volume which is optimal for aerobic cultivation

<sup>\*\* –</sup> Conversion coefficients from OD<sub>850</sub> to OD<sub>600</sub> vary between strains and phases of growth



Reverse-Spin® Technology – Innovative Principle of Microbial Cultivation on web page biosan.lv/rts-tech






Figure 1. 3D graph of E.coli BL21 growth kinetics showing the effect of different temperatures in 7 parallel RTS bioreactors.

#### TYPICAL APPLICATIONS

- · Fermentation real-time growth kinetics
- · Clone candidate screening
- · Protein expression
- · Temperature stress and fluctuation experiments
- · Media screening and optimization
- · Growth characterization
- · Inhibition and toxicity tests

· Strain quality control

RTS-1 and RTS-1C are personal bioreactors that utilize unique Reverse-Spin® technology that applies non-invasive, mechanically driven, low energy consumption, innovative type of agitation where cell suspension is mixed by the single-use falcon bioreactor tube rotation around its axis with a change of direction of rotation motion resulting in highly efficient mixing and oxygenation for aerobic cultivation. Combined with a near-infrared optical system, it is possible to register cell growth kinetics non-invasively in real-time.

#### **FEATURES**

- Reverse-Spin® mixing principle in 50 ml falcon tubes allows to achieve high k<sub>L</sub>a (h-1) up to 450, which is essential for efficient aerobic cultivation;
- Individually controlled bioreactor accelerates optimization process;
- Possibility to cultivate microaerophilic and obligate anaerobic microorganisms (not strict anaerobic
- Reverse-Spin® mixing principle enables non-invasive biomass measurement in real-time;
- Near-infrared optical system makes it possible to register cell growth kinetics;
- Free of charge software for storage, demonstration and analysis of data in real-time;
- Compact design with a low profile and small footprint for personal application;
- Temperature control for bioprocess applications;
- Active cooling for rapid temperature control, e.g. for temperature fluctuation experiments;
- Task profiling for process automatization;
- Cloud data storage possibility to remotely monitor the process of cultivation while at home or using a mobile phone.

#### **SOFTWARE FEATURES**

- Real-Time cell growth logging;
- 3D graphical representation of OD or growth rate over time over unit;
- Pause option;
- Save/Load option;
- Report option: PDF and Excel;
- Connect up to 10 units simultaneously to 1 computer;
- Remote monitoring option (requires internet connection);
- Cycling/Profiling options;
- · User manual calibration possibility for most cells.

ORDERING INFORMATION	Cat. number
RTS-1C including TubeSpin® Bioreactor 50, TPP®, 20 pcs.	BS-010160-A04
RTS-1 including TubeSpin® Bioreactor 50, TPP®, 20 pcs.	BS-010158-A04
Optional accessories:	
TubeSpin® Bioreactor 50, TPP®, 20 pcs.	BS-010158-AK
TubeSpin® Bioreactor 50, TPP®, 180 pcs.	BS-010158-CK
USB 2.0 Hub 10 × ports	BS-010158-BK

# Recommendations for creating personal settings for cultivation of microorganisms. Points that should be considered:

- The growth rate directly depends on the tube's rotation speed, since it is directly proportional (in the range from 1,500 to 2,500 rpm) with the rate of saturation of the medium with oxygen.
- Naturally, with aerobic metabolism, the change in OD over time will also proportionally increase depending on same as above.
- 3. This will also affect the specific growth rate  $\Delta OD/\Delta t$ .
- As well as the time for the growth curve to reach the stationary growth phase during aerobic fermentation (the higher the tube rotation speed, the faster the culture's exit to the stationary phase)
- The saturation of the medium with oxygen will depend on the frequency of switching the tube rotation to the opposite (RST) (the more often the direction of rotation of the tube is reversed, the higher the oxygen mass transfer)
- 6. OD λ=850 this wavelength is used to measure microorganism cell concentration because nutrient media and microorganism cells have colour. This must be taken into account when monitoring the specific dynamics of microorganism growth. In order to go into the "shadow" region (independent of the colour of the medium and microorganism), we offer the near infrared (not visible to the human eye) light scattering measurement range of 850 nm. Since we are still in the sensitive range and, at the same time, are independent of the natural colouration of microorganism colonies. The conversion factor OD<sub>850</sub>/OD<sub>500</sub> is about 2.

It is known that the aerobic bacterial growth is influenced by efficient gas exchange. Figure 2 a-c, serves as an example of growth optimization and illustrates the relationship between RST and gas exchange. As RST decreased the specific growth rate, and biomass yield increased, the highest aeration and optimal growth conditions for *E.coli BL21* optimized at 2000 RPM 1 s RST.

#### Cell growth depending on rotation intensity




Figure 2, a-c. Influence of Reverse Spin Time (RST) on growth kinetics of *E.coli BL21* in  $OD_{600}$ . (a-c) Biomass growth; (b-d) Specific growth rate; throughout cultures were grown in 50 ml TPP Bioreactor tubes, 30% filling volume, 2,000 RPM, RST 1, 2, 4, 8, 16, 30 seconds, LB medium and 37°C temperature, to convert  $OD_{850}$  to  $OD_{600}$  simply multiply  $OD_{850}$  by 1.9.

#### k<sub>1</sub>a (h-1) results in RTS-1/C




Figure 3. Determination of  $k_L a$  in 50 ml TPP Bioreactor tubes. The bioreactor vessels were filled with 5, 10, 20, 30 ml deionized water, and measurements were made by non-invasive  $O_2$  sensors and optics (PreSens, Regensburg, Germany) at 37°C using the gassing-out method. Mean and standard deviation of at least five independent experiments are shown.

The  $k_La$  was measured in 5, 10, 20, 30 ml of deionized water in 50 ml TPP Bioreactor tubes at agitation rate of 2,000 rpm and 1 s RST, this agitation rate was found optimal for Reverse-Spin® mixing principle during initial optimization studies. Over the working volume range, the  $k_La$  increased with the decrease of liquid volume (Figure 3). At the smallest working volume of 5 ml, the highest  $k_La$  of 350 ±26 h<sup>-1</sup> was reached.

#### Cells successfully cultivated

Saccharomyces cerevisiae, Pichia pastoris, Yarrowia lipolytica, Bacillus subtilis, Escherichia coli, Lactobacillus acidophilus, Bifidobacterium bifidum, Pseudomonas aeruginosa, Hybridoma, Jurkat and CHO cells.

#### Types of recommended tubes

For aerobic microorganisms, it is recommended to use tubes that are supplied by TPP – TubeSpin® Bioreactor 50ml. For obtaining optimal results growing aerotolerant anaerobes, it is required to seal the screw cap of TPP TubeSpin® Bioreactor 50ml by tape or purchase TPP TubeSpin® 50 ml falcon tubes without the membrane filter. It is also possible to use other manufacturer tubes of the same type, e.g. Corning® 50 ml Mini Bioreactor, but the device rotor must be modified. It is possible to request this specific modification.

#### Factory calibration particle size and calibration coefficients 600nm/850nm

Factory calibration of the instrument is designed for rod-shaped bacteria size of *E.coli BL21*. In case of exceeding this size, the measurement system will not work correctly. Optical density  $OD_{850}$  to  $OD_{600}$  conversion coefficient of the factory calibration is equal to 1.9.

#### Factory calibration growth phase influence on measurement accuracy

During the growth transition of Escherichia coli culture from exponential growth to the stationary phase, many morphological and physiological changes occur, including cell volume decrease and cell shape change. Therefore, if cells were taken for referent measurement using a spectrophotometer at different stages from the stationary phase, then the correctness of measurement will be worse than specified.

#### **Conversion rate coefficient of user calibration**

Optical density  $OD_{850}$  to  $OD_{600}$  nm conversion rate coefficient depends on the cell size and volume. Therefore, the coefficient will be different for other cell sizes. The device can be calibrated at desired reference wavelength to meet user's needs, e.g.  $OD_{600}$ .

#### Do you want to test this system?

We can provide demo units for 50% of the price for testing or creating an application note. For such, inquiries please contact our R&D department directly at igor@biosan.lv.

## RTS-8 and RTS-8 Plus, Multi-channel Bioreactors

RTS-8 and RTS-8 plus are multi-channel bioreactors that utilise unique Reverse-Spin® technology that applies non-invasive, mechanically driven, low energy consumption, innovative type of agitation where cell suspension is mixed by the single-use falcon bioreactor tube rotation around its axis with a change of direction of rotation motion resulting in highly efficient mixing and oxygenation for aerobic cultivation.

Combined with a near-infrared measurement system, it is possible to register cell growth kinetics and additionally on RTS-8 plus fluorescence and luminescence measurement systems used to register pH and  $O_2$  non-invasively in real-time. For pH and  $O_2$ , innovative single-use sensor spots are used inside the tubes. Although  $O_2$  supply is one of the major issues in the cultivation of aerobic organisms, especially in oxygen-limited conditions, adequate methods for real monitoring of dissolved oxygen were missing, and sufficient  $O_2$  supply was usually assumed. Innovative non-invasive oxygen sensors integrated into falcon tubes now enable online oxygen monitoring and give new insights into metabolic activities.

The pH is one of the major issues in the cultivation of cells, yeast or bacteria. Cultivation vessels, that are sensor-limited, are widely applied in academic and industrial bioprocess development. As adequate methods for real monitoring of pH were not available, cumbersome at-line sampling was used, lacking high data density and interfering with growth. Non-invasive real-time pH measurement provides new insights into metabolic activity and changes in metabolic pathways.

#### SOFTWARE FEATURES

- · Real-Time cell growth logging;
- 3D graphical representation of OD or growth rate over time over unit;
- · Pause option;
- · Save/Load option;
- · Report option: PDF and Excel;
- Remote monitoring option (requires internet connection);
- · Cycling/Profiling options;
- · User manual calibration possibility for most cells.

only RTS-8 Plus

Real-Time pH and O<sub>2</sub> measurement and logging

Tube for RTS-8 Plus with sensor









#### **FEATURES**

- Parallel cultivation of 8 tube bioreactors enables to save time and resources for bioprocess optimization;
- Individually controlled bioreactor accelerates optimization process;
- Possibility to cultivate microaerophilic and obligate anaerobic microorganisms (not strict anaerobic conditions);
- Reverse-Spin® mixing principle enables non-invasive biomass measurement in real time;
- Near-infrared optical system makes it possible to register cell growth kinetics;
- Free of charge software for storage, demonstration and analysis of data in real time;
- Compact design with low profile and small footprint for personal application;
- Individual temperature control for bioprocess applications;
- Active cooling for rapid temperature control, e.g. for temperature fluctuation experiments;
- Task profiling for process automatization;
- Cloud data storage to remotely monitor the process of cultivation while at home or using a mobile phone.

only RTS-8 Plus

 Non-invasive O<sub>2</sub> and pH measurement allows for accurate monitoring of metabolic activitiesf cultivation while at home or using a mobile phone

#### Advantages of the sensor spots:

- They are small;
- Their signal does not depend on the flow rate of the sample;
- They can be physically divided from the measuring system which allows a non-invasive measurement;
- · They can be used in disposables;
- Therefore, they are ideally suited for the examination of small sample volumes, highly parallelized measurements in disposables, and biotechnological applications.

## RTS-8 and RTS-8 Plus, Multi-channel Bioreactors

#### **TYPICAL APPLICATIONS**

- · Fermentation real time growth kinetics;
- · Clone candidate screening;
- · Protein expression;
- · Temperature stress and fluctuation experiments;
- · Media screening and optimization;
- · Growth characterization;
- · Inhibition and toxicity tests;
- · Strain quality control;
- · Initial bioprocess optimization studies.



Product video

	RTS-8	RTS-8 Plus
Light source	Laser	
Measurement wavelength (λ)	850 ±1	5 nm
Measurement range	0-100	OD <sub>600</sub>
E. coli factory calibration measurement range	0-50 C	DD <sub>600</sub>
S. cerevisiae factory calibration measurement range	0-75 C	DD <sub>600</sub>
Achievable user calibration measurement error (range 0.1–6 OD <sub>600</sub> )	±0.	3
Achievable user calibration measurement error (range 6–50 OD <sub>600</sub> )	≤59	6
Achievable user calibration measurement error (range 50–75 OD <sub>600</sub> )	≤10	%
Measurement periodicity per hour	1-60 (increm	ent 1 min)
Temperature setting range +15°C +60°C		+60°C
Temperature control range	+15°C below ambient +60°C (increment 0	
Temperature stability	±0.3°C	
Sample temperature accuracy (20–37°C)	±1°	С
Tube sockets	8	
Sample working volume range	3–50	ml
Speed control range	150-2,700 rpm (in	crement 1 rpm)
Display	LCI	)
Dimensions (W×D×H)	350 × 690 × 300 mm	
Weight	20 kg	
Nominal operating voltage	AC 230 V, 50 Hz 3.15 A / 500 W	
Power consumption		
O <sub>2</sub> sensor* / pH sensor**	-/-	+/+

#### \*O<sub>2</sub> sensor

Range	0–100%
Accuracy	$\pm 0.05\%$ O <sub>2</sub> at 0.2%, $\pm 0.4\%$ O <sub>2</sub> at 20.9%
Drift	<0.03% O <sub>2</sub> within 30 days
Temperature range	up to 40°C
Response time (t90)	<6 s
Storage stability	18 months

#### \*\*pH sensor

Range	4.0-8.5 pH
Accuracy	±0.10 pH at pH 7
Drift	<0.005 pH per day
Temperature range	up to 40°C
Response time (t90)	<120 s
Storage stability	18 months

#### ORDERING INFORMATION

Cat. number

Cat. number

RTS-8

Calibration *E.coli* BS-010168-A01
Optional calibration

Calibration *S.Cerevisiae* BS-010168-A09
Calibration *E.coli and S.Cerevisiae* BS-010168-A10

Including TPP TubeSpin® Bioreactor vessels 50 ml, 20 pcs

#### RTS-8 Plus

Calibration *E.coli* BS-010170-A01 Optional calibration

Calibration *S.Cerevisiae*Calibration *E.coli and S.Cerevisiae* 

BS-010170-A08 BS-010170-A11

Supports pH and pO $_2$  measurement, including TPP TubeSpin® Bioreactor vessels 50 ml, 20pcs and sterile TPP TubeSpin® Bioreactor vessels 50 ml with pH and O $_2$  sensors, 10 pcs

#### Optional accessories:

TubeSpin® Bioreactor 50, TPP®, 20 pcs.	BS-010158-AK
TubeSpin® Bioreactor 50, TPP®, 180 pcs.	BS-010158-CK
USB 2.0 Hub 10 × ports	BS-010158-BK
Sterile TPP TubeSpin® Bioreactor vessel, 50 ml, with pH and O₂ sensors, 1 pce.	BS-010170-AK

DESCRIPTION

## ES-20, Orbital Shaker-Incubator

The **ES-20** is a compact bench-top Shaker-Incubator used for mixing of biological liquids and incubation and cultivation of biological liquids according to the operator set program.

Built-in microprocessor thermocontroller provides constant temperature control in the incubator chamber. Forced heated air circulation inside the transparent plexiglass chamber guarantees even temperature distribution. Dismountable construction makes transportation easy.

Orbital shaking is controlled by the digital tachometer (rpm) and Digital time setting regardless of the temperature. The unit is equipped with a direct-drive system, ensuring the most reliable, stable long-time operation (up to 30 day nights).

The **ES-20** is extremely easy to operate, with a very straightforward setup of temperature, speed and time, using the two-line set-up and status display, which clearly indicates both set and actual values for each of the three parameters.

## DIFFERENT INTERCHANGEABLE PLATFORMS ALLOW USING **ES-20** FOR:

- Growing cell cultures in flasks and other laboratory glassware
- Extracting tissue samples at physiological temperatures
- Other sample preparation processes

Temperature setting range	+25°C +42°C
Speed control range	50-250 rpm
Temperature control range	5°C above ambient +42°C
Setting resolution	0.1°C; 10 rpm
Temperature stability	±0.5°C
Temperature accuracy	±0.5°C
Temperature uniformity	±0.5°C
Orbit	10 mm
Display	LCD, 2 × 16 signs
Digital time setting	1 min–96 h/non-stop (1 min increment)
Timer sound signal	yes
Plexiglas walls thickness	7 mm
Maximum load	2.5 kg
Overall dimensions (W×D×H)	340 × 340 × 435 mm
Dimensions of the inner chambe	er 305 × 260 × 250 mm
Weight	13.2 kg
Nominal operating voltage	230 V, 50/60 Hz or 120 V, 50/60 Hz
Power consumption (230/120 V)	160 W (0.7 A)/ 170 W (1.6 A)





Product video



#### Heat up time for ES-20:

from 25°C 16 min to 42°C

ORDERING INFORMATION:		Cat. number		
ES-20 without platform		BS-010111-AAA		
Optional acc	essories:			
Platforms:	UP-12	BS-010108-AK		
	PP-4	BS-010108-BK		
	P-12/100	BS-010108-EK		
	P-6/250	BS-010108-DK		
	P-16/88	BS-010116-BK		

Description of all platforms for ES-20 on page 20



## ES-20/60, Orbital Shaker Incubator









#### Heat up time for ES-20/60





Orbital Shaker-Incubator ES-20/60 for biotechnological and pharmaceutical laboratories is a professional category equipment designed to cultivate microorganisms and eukaryotic cells, including animal, plant and insect cells. It is also possible to cultivate thermophilic bacteria in ES-20/60 shaker-incubator.

Shaker is equipped with a direct-drive mechanism for platform motion. It provides a reliable and stable operation for the long term experiments needed for cell growth.

Shaker-Incubator ES-20/60 provides smooth or intensive mixing in flasks installed on the platform.

Built-in noiseless thermoresistant brushless fan provides precise temperature distribution inside the chamber (adjustable for up to +80°C). The inner chamber is made of stainless steel. State-of-the-art motor, newest thermal insulation materials, soft-start of the platform motion and temperature PID-control decrease the energy consumption and make the Shaker-Incubator highly energy efficient despite its relatively large size.

Temperature setting range	+25°C +80°C
Speed control range	50-250 rpm
Temperature control range	10°C above ambient +80°C
Setting resolution	0.1°C; 10 rpm
Temperature stability	±0.5°C
Temperature accuracy	±2°C
Temperature uniformity	±2°C
Orbit	20 mm
Display	LCD, 2×16 signs
Digital time setting	1 min–96 h/non-stop (1 min increment)
Timer sound signal	yes
Maximum load	8 kg
Overall dimensions (W×D×H)	590 × 525 × 510 mm
Dimensions of the inner cham	ber 460 × 400 × 310 mm
Weight	41.1 kg
Nominal operating voltage	230 V, 50/60 Hz or 120 V, 50/60 Hz
Power consumption	450 W (2 A)/450 W (4.5 A)
ORDERING INFORMATION:	Cat. number

ES-20/60 without platform

BS-010135-AAA

### ES-20/80 and ES-20/80C Orbital Shaker-Incubator

**ES-20/80, ES-20/80C** shakers-incubators for biotechnological and pharmaceutical laboratories is a professional category equipment. The typical applications include microbial and cell culture cultivation, protein expression, solubility studies, general mixing, as well as other various applications in the fields of biology and chemistry. The unit is equipped with a newly developed triple eccentric mechanism for platform motion that provides supreme balancing characteristics, superior reliability and quiet operation. The achieved stability of the unit during vigorous mixing allows for stacking installation of up to 3 units which enables to save space.

The new display and easy to use user interface provide a clear and intuitive control of parameters and also allow data logging, storage and display over time. Additional features like out-of-balance sensor and automatic thermostat failure detection make this shaker-incubator an advanced and safe product. Bluetooth connectivity to PC allows for data management, data logging, parameter control and profiling in the included software.

A combination of an electric heating element and a built-in heat-resistant brushless fan provides precise temperature distribution inside the chamber up to +80 °C. In **ES-20/80C** shaker-incubator cooling of the chamber is provided by additional Peltier modules allowing to cool until 12.5 °C below the ambient temperature. Additionally, excellent sample temperature uniformity is achieved. The inner chamber is made of stainless steel. State-of-the-art motor, thermal insulation materials and parameter PID-control decrease the energy consumption and make the shaker-incubator highly energy efficient despite its relatively large size.

Stacking kit for 3× ES-20/80, ES-20/80C





## ES-20/80 and ES-20/80C, Orbital Shakers Incubators

	ES-20/80	ES-20/80C	
Temperature setting range	+25°C +80°C +4°C +80°C		
Speed control range	50–400 rpm (increment 10 rpm)		
Temperature control range	5°C above ambient +80°C 12.5°C below ambient +80°C		
Setting resolution	0.1°C; 10 rpm		
Temperature stability	-	±0.1°C at 37°C	
Temperature accuracy	±0.1°C at 37°C		
Temperature uniformity	±0.3°C at 37°C	±0.2°C at 37°C	
Orbit	20 mm		
Display	TFT, 5 inches		
Digital time setting	1 min-96 h/non-stop (1 min increment)		
Timer sound signal	yes		
Maximum load	10.6 kg		
Data transfer	Bluetooth®		
Stacking	up to 3*		
PC system requirements	Intel/AMD Processor, 1 GB RAM, Windows 10/11, USB, Bluetooth		
Overall dimensions (W×D×H)	620 × 530 × 510 mm 620 × 620 × 510 mm		
Inner chamber dimensions (W×D×H)	460 × 400 × 310** mm 460 × 400 × 325** mm		
Weight	48 kg	50 kg	
Nominal operating voltage	230 V, 50/60 Hz		
Power consumption	500 W (2.2 A)		

<sup>\*</sup> Additional stacking kit required

<sup>\*\*</sup> Height is measured from platform surface

ORDERING INFORMATION	Cat. number
ES-20/80 with software, without platform	BS-010167-A05
ES-20/80C with software, without platform	BS-010173-A01
Optional accessories:	
USB Bluetooth® adapter	BS-010425-FK
Stacking kit for 2× <b>ES-20/80</b> , <b>ES-20/80C</b>	BS-010167-OK
Stacking kit for 3× ES-20/80. ES-20/80C	BS-010167-PK

Platforms cat. numbers for ES-20/80, ES-20/80C can be found on next page





Platforms cat. numbers for ES-20/80, ES-20/80C can be found on next page

## Platforms for ES-20/80 and ES-20/80C

Platform		Description	Dimensions (Working Area)	Cat. number
HSP-30/100		Platform with 30 tight fit clamps for 100–150 ml flasks	360 × 400 mm (360 × 400 mm)	BS-010167-KK
HSP-16/250		Platform with 16 tight fit clamps for 250–300 ml flasks	360 × 400 mm (360 × 400 mm)	BS-010167-MK
HSP-9/500		Platform with 9 tight fit clamps for 500 ml flasks	360 × 400 mm (360 × 400 mm)	BS-010167-NK
HSP-6/1000		Platform with 6 tight fit clamps for 1,000 ml flasks	360 × 400 mm (360 × 400 mm)	BS-010167-LK
PP-400	3	Flat platform with non-slip silicone mat	360 × 400 mm (360 × 400 mm)	BS-010135-FK
UP-168		Universal platform for different flasks (Clamps ordered separately)	360 × 400 mm (360 × 400 mm)	BS-010135-JK
HSC-50 HSC-100 HSC-250 HSC-500 HSC-1000 FC-2000 (m	aax 200 rpm)	Tight fit clamp for 50, 100, 250, 500, 1,000, 2,000 ml flask (for UP-168)	Ø 50 mm Ø 65 mm Ø 85 mm Ø 105 mm Ø 130 mm Ø 160 mm	BS-010167-DK BS-010167-EK BS-010167-FK BS-010167-JK BS-010167-JK BS-010126-NK
SPML		Set of 3 double-sided adhesive strips as an alternative for regular flask clamps (for UP-168)	390 × 80 × 3 mm	BS-010135-MK
TR-21/50		Test tube rack for 50 ml with 21 drillings (for UP-168)	340 × 124 mm (2 per platform)	BS-010135-KK
TR-44/15		Test tube rack for 15 ml with 44 drillings (for UP-168)	340 × 124 mm (2 per platform)	BS-010135-LK
P-EX For ES-20/60 ES-20/80, ES-20/80C		Platform Exchange Set (Slider + Trato install 4 × 2l flasks.Accommodate P-6/1000, P-9/500, P-15/250, P-30/1000, HSP-9/500, HSP-15/250, HSP-30/100	e: UP-168, 00, HSP-6/1000,	BS-010173-CK



## S-Bt Smart Biotherm, Compact CO<sub>2</sub> Incubator

**S-Bt Smart Biotherm** is designed for work in the areas of cell biology (operations with animal cell cultures and tissues), molecular biology (DNA/RNA reaction analysis, hybridization reactions), biotechnology (synthesis of target proteins and other molecules), immunology (synthesis of antibodies and other proteins of the immune system). Unit provides six-sided heating: the heating elements are located on the walls and on the door, thus providing excellent uniform temperature distribution, regardless of external factors, such as ambient temperature and positioning of the device.

Built-in infrared  $\rm CO_2$ -sensor allows accurate control of the  $\rm CO_2$  level. The sensor makes measurement non-sensitive to changes in temperature and humidity inside the incubator.

The chamber is made of stainless steel with smoothed seams to minimize contamination and to facilitate cleaning.

**S-Bt** is equipped with a UV air recirculation system – 1 UV lamp and a fan are mounted behind the rear wall, providing decontamination of the working volume.

A convenient access port is built in the wall of the incubator for easy output of wire sensors or devices' installed inside. The access port is heated independently to prevent the formation of condensate.

Unit is equipped with error tracing and alarm systems, which significantly lower potential risks during operation.

Unit is equipped with a "black box" system that records temperature, humidity and CO<sub>2</sub> levels, as well as statuses for door opening, UV lamp, fan and errors, to the internal memory. Bluetooth® connection to PC is available.

Chamber Material	Stainless steel (1 mm)
Temperature setting range	+25°C +60°C
Temperature stability	±0.1°C
Temperature uniformity @37	7°C ±0.3°C
Timer sound signal	yes
Working volume	46 I
Number of shelves	3 (max 6)
Inner door	Glass
Relative humidity	>90% @ 37°C
Humidity delivery	Water bath
CO <sub>2</sub> control range*	0-20%
CO <sub>2</sub> sensor	Infrared sensor
Temperature and CO <sub>2</sub> level in	nput Digital
UV lamp	1 × 6 W, TUV G6T5
Data transfer	Wireless
Access port	1 (Ø 26 mm)
Working voltage	230V, 50/60 Hz; 115 V, 50/60 Hz
Power consumption	600 W
Weight	37.7 kg
Dimensions (W×D×H)	500 × 560 × 550 mm
Inner chamber dimensions (W	/×D×H) 350 × 330 × 390 mm

<sup>\* –</sup> At set temperature from ambient to 50 °C



#### APPLICATION AREAS

- Cell biology: operations with animal cell cultures and tissues;
- Molecular biology: DNA/RNA reaction analysis, hybridization reactions;
- Biotechnology: synthesis of target proteins and other molecules;
- Immunology: synthesis of antibodies and other proteins of immune system.

#### **FEATURES**

- Six-sided heating provides uniform distribution of the temperature inside the chamber;
- Infrared CO<sub>2</sub> sensor, non-sensitive to temperature and humidity changes;
- UV recirculation system for decontamination cycles;
- · Bluetooth data transfer to PC;
- «Black box» parameter logging system;
- · Error tracing and alarm system;
- Separately heated lockable port for chamber access for cables.

## **S-Bt Smart Biotherm,** Compact CO<sub>2</sub> Incubator

Simple CO<sub>2</sub> tank connection



Air UV recirculation system in the chamber



Gas purification filter



PC software









ORDERING INFORMATION	Cat. number
S-Bt Smart Biotherm, PC software included + RS6, rack with 3 shelves	BS-010425-A01
S-Bt Smart Biotherm, PC software included + RS2, rack for CPS-20/CTR-6 installation	BS-010425-A10
Optional accessories:	
CPS-20, CO <sub>2</sub> Shaker	BS-010172-A01
CTR-6, CO <sub>2</sub> Tube Roller	BS-010174-A01
Shelf	BS-010425-AK
USB Bluetooth® adapter	BS-010425-FK
Incubator stacking device	BS-010425-CK

## CPS-20, CO<sub>2</sub> Shaker

 ${\rm CO_2}$  Shaker **CPS-20** provides regulated orbital motion of the platform and is designed for use specifically in  ${\rm CO_2}$  incubators. **CPS-20** is specifically designed for unnecessary harsh environments such as  ${\rm CO_2}$  and humidity and provides reproducible results for cell culture growth. A choice of five interchangeable platforms provides the possibility of performing various procedures and techniques in various cultivation vessels.

**CPS-20** incorporates a brushless motor with a guaranteed service life of up to 35,000 hours and a direct-drive system, ensuring the most reliable, stable long–time operation (up to 30 day nights). The unit is equipped with a triple eccentric mechanism for platform motion that provides supreme balancing characteristics, superior reliability and quiet operation. The specially designed remote controller allows for the protection of electronics from a  ${\rm CO}_2$  incubator environment, as well as, the remote control minimizes interference with the incubator environment and the ongoing experiment.

Speed control range	50-250 rpm
	(increment 10 rpm)
	max. speed depends
	on the load and vessels' shape

Digital time setting	1 min-96 h/non-stop
	(increment 1 min)

	(increment i iiiii)
Digital speed control	+
Maximum continuous operation tir	me 168 h
Orbit	20 mm
Maximum load	3 kg
Overall dimensions (W×D×H)	255 × 255 × 100 mm
Weight	3.4 kg
Input current/power consumption	470 mA / 5.7 W

External power supply Input AC 100–240 V; 50/60 Hz Output DC 12 V











## Platforms for **CPS-20**

Platform	Description	Dimensions (mm)	Working Area (mm)	Cat. number
UP-12	Universal platform with adjustable bars for different types of flasks, bottles and beakers with silicone mat	285 × 220 × 40	270 × 185 × 40	BS-010108-AK
Bio PP-4	Flat platform with silicone mat for Petri dishes, culture flasks, agglutination cards	255×255	230 × 230	BS-010116-AK
P-12/100	Platform with clamps for flasks, 100–150 ml (12 places)	250×190	250 × 190	BS-010108-EK
P-6/250	Platform with clamps for flasks, 250–300 ml (6 places)	250×190	250×190	BS-010108-DK
P-16/88	Platform with spring holders for up to 88 tubes up to 30 mm diameter (e. g. 10 ml, 15 ml, 50 ml tubes)	275 × 205 × 75	275 × 205 × 75	BS-010116-BK

S-Bt Smart Biotherm, Compact CO<sub>2</sub> Incubator



CPS-20, CO<sub>2</sub> Shaker BS-010172-A01

RS2, rack for CPS-20/CTR-6 installation

BS-010425-HK

DESCRIPTION

## CTR-6, CO<sub>2</sub> Tube Roller

CO2 Tube Roller CTR-6 provides regulated rocking and rolling of maximum up to 6 rollers and is designed for use specifically in CO2 incubators. CTR-6 is specifically designed for use in harsh environments such as CO2 and humidity and provides reproducible results for cell culture growth. Possibility to remove rollers makes the unit flexible and allows for performing various procedures and techniques in various cultivation vessels. The specially designed remote controller allows for protection of electronics from CO2 incubator environment, as well as does not interfere with the experiment.

Tube roller CTR-6 incorporates a stepper motor with a guaranteed service life up to 10,000 hours. It is possible to stack up to 2 units, saving valuable bench space. Typical applications include cells cultivation (eukaryotic, microbial) and general mixing (resuspension, viscous and liquid-solid suspensions).

Speed control range*	5–80 rpm
Increment	1 rpm
Digital time setting	1 min – 96 h or non-stop
Increment	1 min
Maximum load	3 kg
Tilt angle	4°
Dimensions (W×D×H)	310 × 262 × 80 mm
Input current/power consumpt	tion 12 V, 415 mA/ 5 W
Weight, accurate within ±10%	3 kg
External power supply	Input AC 100–240 V; 50/60 Hz Output DC 12 V

<sup>\*</sup> Maximum speed depends on the load on the platform and the shape of the vessels













#### ORDERING INFORMATION

Cat. number

CTR-6, CO<sub>2</sub> Tube Roller

BS-010174-A01

Optional accessories:

Stacking kit for 2 × CTR-6

BS-010174-BK

Stacking kit for 3 × CTR-6

BS-010174-CK

## LAB DIAGNOSTICS:

## DNA/RNA PURIFICATION, IMMUNODIAGNOSTICS



**HiPo MPP-96**Microplate Photometer

### **Automatic DNA/RNA extraction**





Incubation TS-100



TS-100C



Sample resuspendation and droplet spin down



**Automatic extraction** 

BioMagPure 12 Plus and reagents



### **Manual DNA/RNA extraction** using magnetic beads technology

**Sample resuspendation** and droplet spin down FVL-2400N



**Sample preparation** in UV-Cabinet for PCR

UVC/T-M-AR, or similar, see UV-Cabinets for PCR

**Mixing and** resuspendation





Multi Bio RS-24



**Capture** of magnetic beads





Centrifugation •





Vacuum aspiration •



FTA-2i



Incubation TDB-120



TS-100C




# **BioMagPure 12 Plus,** Compact Bench-Top Robotic Workstation For Automated Nucleic Acid Purification

The **BioMagPure 12 Plus** consists of a compact bench-top robotic workstation for automated nucleic acid purification. Usage of pre-filled reagent cartridges and disposable consumables enable a true walk-away automation and high-quality nucleic acid extraction solution. Proven magnetic separation technology makes purification efficient, easy to use, reliable, safe and cost-effective.

**BioMagPure 12 Plus** has an ingeniously designed polygonal reaction chamber with patented parts that ensure high efficiencies of lysis and elution through a large contact area of magnet and heating element allowing maximisation of magnetic bead recovery, minimisation of the residues of magnetic beads and alcohols in the final elute product. Specific formation of reaction chamber ensures unrivalled mixing ability and exclude conventional mixing by tip or pipetting thus eliminates cross-contamination possibility.

Reagent kit contains everything for extraction procedure performance, including all necessary plastics, pre-filled reagent cartridges, incubation buffers and solutions for sample pre-treatment (if needed),

With the flexibility of processing 1–12 samples per run, the **BioMagPure 12 Plus** is tailor-made to fit small clinics and early-stage laboratories. By occupying minimal counter space and greatly reducing technician man-hours, this series allows organizations to operate facilities in a much more cost-effective fashion.



	Separation
Processing time	45–60 min
Processing capability	1–12 samples per run
Extraction technology	magnetic particle separation technology
Protocol	programmed by scanning a barcode
Protocol input	barcode scanner
Sample volume	10–2,000 µl (depending on the kit)
Elution volume	50–400 μl
Connection to PC	not required
Display	LCD (20 × 4)
Certification	CE IVD
Nominal operating voltage	110-240 V, 50/60 Hz
Dimensions (W×D×H)	560 × 590 × 510 mm
Weight	55 kg





#### 3 easy steps







#### **Features:**

- Advanced magnetic bead technology;
- Reaction chamber with patented parts;
- Piercing-pin system for elimination of cross-contamination;
- · Walk-away automation;
- · Reliable quality;
- · No PC required;
- · Ready-to-use reagent cartridges;
- 3 easy steps: LOAD-RUN-OBTAIN.

## Reagents for **BioMagPure 12 Plus**

#### ORDERING INFORMATION:

Name	Description	Cat. number
BioMagPure 12 Plus	Compact Bench-Top Robotic Workstation For Automated Nucleic Acid Purification.	BS-060202-AAA
Blood DNA Extraction Kit 200	Blood DNA Extraction Kit is used with the BioMagPure 12 Plus instrument to extract DNA from 10-400µl mammalian whole blood, suspension of mammalian blood cells.	BS-060201-AK
Blood DNA Extraction Kit 1200*	Blood DNA Extraction Kit is used with the BioMagPure 12 Plus instruments to extract gDNA from 400-1000 $\mu$ l mammalian blood, suspension of mammalian blood cells.	BS-060201-BK
Viral Nucleic Acid Extraction Kit*	Viral Nucleic Acid Extraction Kit is used with the BioMagPure 12 Plus instrument to extract Viral DNA or RNA from human biological specimens such as serum, plasma, and other cell-free fluids.	BS-060201-CK
Tissue DNA Extraction Kit*	BioMagPure 12 Plus Tissue DNA Extraction Kit is used with the BioMagPure 12 Plus instrument to extract genomic DNA from a variety of animal tissues, swab samples and bloodstain.	BS-060201-DK
Cultured Cell DNA Extraction Kit*	Cultured Cell DNA Extraction Kit is used with the BioMagPure 12 Plus instrument to extract genomic DNA from culture cells and buffy coat.	BS-060201-EK
Bacterial DNA Extraction Kit	Bacterial DNA Extraction Kit is used with the BioMagPure 12 Plus instrument to extract genomic DNA from both Gram-positive and Gram-negative bacteria.	BS-060201-FK
HPV DNA Extraction Kit for Swab*	HPV DNA Extraction Kit is used with the BioMagPure 12 Plus instrument for DNA extraction of the Human Papillomavirus (HPV) from cervical cell samples which collected by cervical brush or genital swab in liquid-based Medium (e.g. Hologic Thinprep PreservCyt®, BD SurepathTM, etc.) or other STM (sample transport media) preservation solutions(e.g. QIAGEN DNA PAP Cervical sampler, Roche Cobas® PCR Cell Collection Media, Hybribio cell preservation solution, etc.).	BS-060201-GK
TB DNA Extraction Kit*	TB DNA Extraction Kit is used with the BioMagPure 12 Plus instrument for extraction of genomic DNA of <i>Mycobaceteria spp.</i> (e.g. <i>Mycobacterium tuberculosis</i> ) from different specimen.	BS-060201-IK
FFPE DNA Extraction Kit*	FFPE DNA Extraction Kit is used with the BioMagPure 12 Plus instrument for extraction of genomic DNA from FFPE (Formalin-Fixed, Paraffin-Embedded) tissue samples. Providing good quality, high integrity DNA for Molecular diagnosis and research works	BS-060201-JK
Forensic DNA Extraction Kit*	Forensic DNA extraction kit is used to extract and isolate genomic DNA from forensic samples.	BS-060201-KK
Viral/Pathogen Nucleic Acids Extraction Kit A*	Viral/Pathogen Nucleic Acids Extraction Kit A is used with the BioMagPure 12 Plus instrument to extract Viral and bacterial DNA/RNA from cell-free samples, such as serum, plasma, and other cell-free body fluids.	BS-060201-LK
Viral/Pathogen Nucleic Acids Extraction Kit B	Viral/Pathogen Nucleic Acids Extraction Kit B is used with the BioMagPure 12 Plus instrument to extract viral and bacterial DNA/RNA from swab samples (cell-rich samples).	BS-060201-MK
Viral RNA Extraction Kit*	Viral Nucleic RNA Extraction Kit is used with the BioMagPure 12 Plus instrument to extract Viral RNA from human biological specimens such as serum, plasma, and other cell-free fluids.	BS-060201-NK
Plant DNA Extraction Kit*	Plant DNA Extraction Kit is used with the BioMagPure 12 Plus instrument to extract genomic DNA from plants (leaf, seeds and spores) and fungal tissues. Up to 100 mg of tissue can be used for purification.	BS-060201-OK
Total RNA E xtraction Kit*	Total RNA Extraction Kit is used with the BioMagPure 12 Plus instrument to extract total RNA from whole blood, blood cells, animal tissue, plant tissue, yeast or cultured cells.	BS-060201-PK
Viral Nucleic Acid Large Volume Extraction Kit*	Viral Nucleic Acid Large Volume Extraction Kit is used with the BioMagPure 12 Plus instrument to extract Viral DNA or RNA from human biological specimens such as serum, plasma, and other cell-free fluids.	BS-060201-QK
CFC DNA Extraction Kit Large Volume*	CFC DNA Extraction Kit Large Volume - is used with the BioMagPure 12 Plus instrument to extract circulating DNA from plasma serum or cell-free body fluids sample volume ranging up to 5 ml.	BS-060201-RK

<sup>\* -</sup> MOQ is 8 kits (384 reactions)

### MagSorb-16, NEW

### Magnetic rack for manual nucleic acid extraction

MagSorb-16 is a magnetic rack that easily accommodates up to 16 single use tubes (1.5-2 ml).

Magnetic beads based NA extraction methods offer quick and efficient separation of genetic material from cellular leftovers. Wide variety of commercially available extraction kits manufactured by companies like ThermoFisher Scientific, Promega, Qiagen etc. are available in the market.

Our MagSorb-16 magnetic rack can be used together with any manufacturers magnetic beads based NA extraction reagents that allow working in 1.5-2.0 ml tube format.

Number of places in stand	16
Tube's volume	1.5–2 ml
Type of tubes	microtubes
Dimensions (W×D×H)	160 × 50 × 50 mm
Weight	0.24 kg



MagSorb-16, magnetic rack	BS-010601	
ORDERING INFORMATION	Cat. number	$\overline{}$

## **PCR Analysis**

Sample resuspendation and droplet spin down



FVL-2400N (tubes)



MSC-2P (PCR-plates)



**Real time PCR** detection



### BioQuant-96, Real-time PCR detection system



#### DESCRIPTION

**BioQuant-96** is the newest product of Biosan Molecular diagnostic product family.

It has adopted innovative thermoelectric refrigeration technology, brand-new light source and light path design. Detection from the top allows using different consumables – 0,2 tubes, 8-tube strips and semi-skirted and non-skirted 96-well PCR plates. The unique constant current power and 6-zone independent temperature control method ensure more rapid, correct and stable fluorescence quantitative analysis, while maintaining its excellent performance in lowest possible energy consumption.

Device is available in 5-channel configuration. Meanwhile, it has been added with functions including independent temperature control, low temperature storage of sample at 4°C and FAST mode for more faster cycling (confirm reagent compatibility with fast mode).

**BioQuant-96** is comprehensively realizing automatic gain setting and improving user experience. It will fully meet the demand of scientific research laboratories.



Program Max. 20 segments for each program, max. 99 cycles  Continuous operation +  Scan period 5.5 s		
Excitation wavelength  Emission wavelength  Channels  Channels  Temperature setting range  Heating/coolling rate  Heating/coolling rate  Heating/coolling rate  Temperature accuracy  Temperature fluctuation  Temperature uniformity  Sample volume  Gradient temperature range  Hot-lid temperature range  Continuous operation  Continuous operation  Max. 20 segments for each program, max. 99 cycles  Software  BioQuant-96 Software included on the program windows 10/11  Minimum PC requirements  Nominal operating voltage  Dimensions (W×D×H)  Number of sockets  Max. 20 segments for each program, e		cap), 96 × 0.2 (0.1) ml Plate transparent cap), 8-strip tubes
Emission wavelength Channels 5 Temperature setting range Heating/coolling rate Heating/coolling rate Femperature accuracy Temperature fluctuation Temperature uniformity  Emperature control mode Sample volume  Fadient temperature range Hot-lid temperature range Hot-lid temperature range  Frogram  Max. 20 segments for each program, max. 99 cycles  Continuous operation  Continuous operation  Minimum PC requirements  Nominal operating voltage  Power consumption  Dimensions (W×D×H)  Number of sockets  Pose Scapper (JSB adapter, Bluetooth adapter)  Sun Sockets  Solvare Solvare Bloeved (JSB adapter, Bluetooth adapter)  Bush At-99.9°C (increment: 0.1°C)  \$4-99.9°C (increment: 0.1°C)  \$4-99.9°C (increment: 0.1°C)  \$40.1°C  \$40.1°	Dynamic range	1~10 <sup>10</sup> copies
Channels       5         Temperature setting range       4~99.9°C (increment: 0.1°C) range         Heating/coolling rate       6.0°C/s / 5.5°C/s         Temperature accuracy       ≤±0.1°C         Temperature fluctuation       ≤±0.1°C         Temperature uniformity       ≤±0.3°C         Temperature control mode       BLOCK/Tube simulation mode (automatic control based on sample volume)         Sample volume       1~100 μl         Gradient temperature range       30~110°C (adjustable default: 105°C, Automatic Hot-lid)         Repeatability       5%         Scan mode       Entire plate         Program       Max. 20 segments for each program, max. 99 cycles         Continuous operation       +         Scan period       5.5 s         Software       BioQuant-96 Software included         Operating system       Windows 10/11         Minimum PC requirements       Memory: 512 MB, HD: 10 GB, CPU: Pentium* 4, Virtual Memory: ≥1,000 MB         Nominal operating voltage       100-240 V; ~50/60 Hz         Power consumption       600 W         Dimensions (W×D×H)       490 × 290 × 391 mm         Number of sockets       USB adapter, Bluetooth adapter	Excitation wavelength	300-800 nm
Temperature setting range  Heating/coolling rate  6.0°C/s / 5.5°C/s  Temperature accuracy  ≤±0.1°C  Temperature fluctuation  Temperature control mode  BLOCK/Tube simulation mode (automatic control based on sample volume)  Sample volume  Foradient temperature range  Hot-lid temperature range  Hot-lid temperature range  Hot-lid temperature range  Temperature on trol mode  BLOCK/Tube simulation mode (automatic control based on sample volume)  Sample volume  1~100 μl  Gradient temperature range  Hot-lid temperature range  30~110°C (adjustable default: 105°C, Automatic Hot-lid)  Repeatability  5%  Scan mode  Entire plate  Program  Max. 20 segments for each program, max. 99 cycles  Continuous operation  +  Scan period  5.5 s  Software  BioQuant-96 Software included  Operating system  Windows 10/11  Minimum PC requirements  Memory: 512 MB, HD: 10 GB, CPU: Pentium® 4, Virtual Memory: ≥1,000 MB  Nominal operating voltage  Power consumption  600 W  Dimensions (W×D×H)  490 × 290 × 391 mm  Number of sockets  USB adapter, Bluetooth adapter	Emission wavelength	500–800 nm
range  Heating/coolling rate  6.0°C/s / 5.5°C/s  Temperature accuracy  5±0.1°C  Temperature fluctuation  5±0.1°C  Temperature uniformity  5±0.3°C  Temperature control mode  BLOCK/Tube simulation mode (automatic control based on sample volume)  Sample volume  1~100 µl  Gradient temperature range  Hot-lid temperature range  Hot-lid temperature range  Togram  Max. 20 segments for each program, max. 99 cycles  Continuous operation  +  Scan period  Software  BioQuant-96 Software included to Garding system  Windows 10/11  Minimum PC requirements  HD: 10 GB, CPU: Pentium® 4, Virtual Memory: 512 MB, HD: 10 GB, CPU: Pentium® 4, Virtual Memory: ≥1,000 MB  Nominal operating voltage  Power consumption  600 W  Dimensions (W×D×H)  490 × 290 × 391 mm  Number of sockets  USB adapter, Bluetooth adapter	Channels	5
Temperature accuracy  Temperature fluctuation  Temperature uniformity  ≤±0.1°C  Temperature uniformity  ≤±0.1°C  Temperature control mode  BLOCK/Tube simulation mode (automatic control based on sample volume)  Sample volume  1~100 μl  Gradient temperature range  Hot-lid temperature range  Hot-lid temperature range  30~110°C (adjustable default: 105°C, Automatic Hot-lid)  Repeatability  5%  Scan mode  Entire plate  Program  Max. 20 segments for each program, max. 99 cycles  Continuous operation  +  Scan period  5.5 s  Software  BioQuant-96 Software included  Operating system  Windows 10/11  Minimum PC requirements  Memory: 512 MB, HD: 10 GB, CPU: Pentium® 4, Virtual Memory: ≥1,000 MB  Nominal operating voltage  100~240 V; ~50/60 Hz  Power consumption  600 W  Dimensions (W×D×H)  490 × 290 × 391 mm  Number of sockets  USB adapter, Bluetooth adapter		4~99.9°C (increment: 0.1°C)
Temperature fluctuation  Temperature uniformity  ≤±0.3°C  Temperature control mode  BLOCK/Tube simulation mode (automatic control based on sample volume)  Sample volume  1~100 µl  Gradient temperature range  Hot-lid temperature range  Hot-lid temperature range  30~110°C (adjustable default: 105°C, Automatic Hot-lid)  Repeatability  5%  Scan mode  Entire plate  Program  Max. 20 segments for each program, max. 99 cycles  Continuous operation  +  Scan period  5.5 s  Software  BioQuant-96 Software included  Operating system  Windows 10/11  Minimum PC requirements  Memory: 512 MB, HD: 10 GB, CPU: Pentium® 4, Virtual Memory: ≥1,000 MB  Nominal operating voltage  Power consumption  600 W  Dimensions (W×D×H)  490 × 290 × 391 mm  Number of sockets  USB adapter, Bluetooth adapter	Heating/coolling rate	6.0°C/s / 5.5°C/s
Temperature uniformity  Temperature control mode alto mode (automatic control based on sample volume)  Sample volume  Temperature range  Tempera	Temperature accuracy	≤±0.1°C
Temperature control mode mode (automatic control based on sample volume)  Sample volume 1~100 μl  Gradient temperature range 1~36°C  Hot-lid temperature range 30~110°C (adjustable default: 105°C, Automatic Hot-lid)  Repeatability 5%  Scan mode Entire plate  Program Max. 20 segments for each program, max. 99 cycles  Continuous operation +  Scan period 5.5 s  Software BioQuant-96 Software included of the control of the con	Temperature fluctuation	≤±0.1°C
Monde (automatic control based on sample volume)         Sample volume       1~100 μl         Gradient temperature range       1~36°C         Hot-lid temperature range       30~110°C (adjustable default: 105°C, Automatic Hot-lid)         Repeatability       5%         Scan mode       Entire plate         Program       Max. 20 segments for each program, max. 99 cycles         Continuous operation       +         Scan period       5.5 s         Software       BioQuant-96 Software included         Operating system       Windows 10/11         Minimum PC requirements       Memory: 512 MB, HD: 10 GB, CPU: Pentium® 4, Virtual Memory: ≥1,000 MB         Nominal operating voltage       100-240 V; ~50/60 Hz         Power consumption       600 W         Dimensions (W×D×H)       490 × 290 × 391 mm         Number of sockets       USB adapter, Bluetooth adapter	Temperature uniformity	≤±0.3°C
Gradient temperature range  Hot-lid temperature range  Hot-lid temperature range  30~110°C (adjustable default: 105°C, Automatic Hot-lid)  Repeatability  5%  Scan mode  Program  Max. 20 segments for each program, max. 99 cycles  Continuous operation  +  Scan period  5.5 s  Software  BioQuant-96 Software included  Operating system  Windows 10/11  Minimum PC requirements  Memory: 512 MB, HD: 10 GB, CPU: Pentium® 4, Virtual Memory: ≥1,000 MB  Nominal operating voltage  Nominal operating voltage  Power consumption  600 W  Dimensions (W×D×H)  490 × 290 × 391 mm  Number of sockets  USB adapter, Bluetooth adapter	Temperature control mod	mode (automatic control
Hot-lid temperature range  Hot-lid temperature range  Hot-lid temperature range  Repeatability  Scan mode  Program  Max. 20 segments for each program, max. 99 cycles  Continuous operation  +  Scan period  Software  BioQuant-96 Software included  Operating system  Windows 10/11  Minimum PC requirements  HD: 10 GB, CPU: Pentium® 4, Virtual Memory: ≥1,000 MB  Nominal operating voltage  Power consumption  Dimensions (W×D×H)  Number of sockets  30~110°C (adjustable default: 105°C, Automatic Hot-lid)  Hot lid Segments for each program, max. 99 cycles  Footnaments  Horizonaments  Memory: 512 MB, HD: 10 GB, CPU: Pentium® 4, Virtual Memory: ≥1,000 MB  Nominal operating voltage  100-240 V; ~50/60 Hz  Power consumption  G00 W  USB adapter, Bluetooth adapter	Sample volume	1~100 µl
default: 105°C, Automatic Hot-lid)  Repeatability 5%  Scan mode Entire plate  Program Max. 20 segments for each program, max. 99 cycles  Continuous operation +  Scan period 5.5 s  Software BioQuant-96 Software included of the second system Windows 10/11  Minimum PC requirements Memory: 512 MB, HD: 10 GB, CPU: Pentium® 4, Virtual Memory: ≥1,000 MB  Nominal operating voltage 100-240 V; ~50/60 Hz  Power consumption 600 W  Dimensions (W×D×H) 490 × 290 × 391 mm  Number of sockets USB adapter, Bluetooth adapter	Gradient temperature ran	ge 1~36°C
Program  Max. 20 segments for each program, max. 99 cycles  Continuous operation  Scan period  Scan period  Software  BioQuant-96 Software included  Operating system  Windows 10/11  Minimum PC requirements  Memory: 512 MB, HD: 10 GB, CPU: Pentium® 4, Virtual Memory: ≥1,000 MB  Nominal operating voltage  Nominal operating voltage  Power consumption  600 W  Dimensions (W×D×H)  490 × 290 × 391 mm  Number of sockets  USB adapter, Bluetooth adapter	Hot-lid temperature rang	default: 105°C,
Program  Max. 20 segments for each program, max. 99 cycles  Continuous operation +  Scan period 5.5 s  Software BioQuant-96 Software included  Operating system Windows 10/11  Minimum PC requirements Memory: 512 MB, HD: 10 GB, CPU: Pentium® 4, Virtual Memory: ≥1,000 MB  Nominal operating voltage 100−240 V; ~50/60 Hz  Power consumption 600 W  Dimensions (W×D×H) 490 × 290 × 391 mm  Number of sockets USB adapter, Bluetooth adapter	Repeatability	5%
program, max. 99 cycles  Continuous operation +  Scan period 5.5 s  Software BioQuant-96 Software included  Operating system Windows 10/11  Minimum PC requirements Memory: 512 MB, HD: 10 GB, CPU: Pentium® 4, Virtual Memory: ≥1,000 MB  Nominal operating voltage 100-240 V; ~50/60 Hz  Power consumption 600 W  Dimensions (W×D×H) 490 × 290 × 391 mm  Number of sockets USB adapter, Bluetooth adapter	Scan mode	Entire plate
Scan period 5.5 s  Software BioQuant-96 Software included of the second	Program	
Software BioQuant-96 Software included Operating system Windows 10/11 Minimum PC requirements Memory: 512 MB, HD: 10 GB, CPU: Pentium® 4, Virtual Memory: ≥1,000 MB  Nominal operating voltage 100-240 V; ~50/60 Hz  Power consumption 600 W  Dimensions (W×D×H) 490 × 290 × 391 mm  Number of sockets USB adapter, Bluetooth adapter	Continuous operation	+
Operating system Windows 10/11  Minimum PC requirements HD: 10 GB, CPU: Pentium® 4, Virtual Memory: ≥1,000 MB  Nominal operating voltage 100–240 V; ~50/60 Hz  Power consumption 600 W  Dimensions (W×D×H) 490 × 290 × 391 mm  Number of sockets USB adapter, Bluetooth adapter	Scan period	5.5 s
Minimum PC requirements HD: 10 GB, CPU: Pentium® 4, Virtual Memory: ≥1,000 MB  Nominal operating voltage 100-240 V; ~50/60 Hz  Power consumption 600 W  Dimensions (W×D×H) 490 × 290 × 391 mm  Number of sockets USB adapter, Bluetooth adapter	Software	BioQuant-96 Software included
HD: 10 GB, CPU: Pentium® 4, Virtual Memory: ≥1,000 MB  Nominal operating voltage 100–240 V; ~50/60 Hz  Power consumption 600 W  Dimensions (W×D×H) 490 × 290 × 391 mm  Number of sockets USB adapter, Bluetooth adapter	Operating system	Windows 10/11
Power consumption 600 W  Dimensions (W×D×H) 490 × 290 × 391 mm  Number of sockets USB adapter, Bluetooth adapter	Minimum PC requirement	HD: 10 GB, CPU: Pentium <sup>®</sup> 4,
Dimensions (W×D×H) 490 × 290 × 391 mm  Number of sockets USB adapter, Bluetooth adapter	Nominal operating voltage	e 100-240 V; ~50/60 Hz
Number of sockets USB adapter, Bluetooth adapter	Power consumption	600 W
Bluetooth adapter	Dimensions (W×D×H)	490 × 290 × 391 mm
Interface Built-in Touchscreen	Number of sockets	
	Interface	Built-in Touchscreen







### **Enzyme-linked** immunosorbent assay (ELISA)

**Sample preparation** in laminar flow cabinet

**Biological Safety Cabinet class II** 







Incubation

PST-60HL

PST-60HL-4

PSU-2T















Washing (Manual)

FTA-2i

FTA-U





Reading and Analysis

HiPo MPP-96 Microplate Photometer with QuantAssay software

### IW-8, Intelispeed Washers



Intelispeed Washer **IW-8** is designed to wash standard flat-bottom (two point aspiration) and U-shape (only in single point aspiration) 96 well plates and microstrips. The unit is fully programmable, ensuring multi-step solution ripening, aspiration (aspiration, combination of aspiration/liquid dispensing and soaking, as well as soaking cycle during a particular time).

The unit has 100 user-defined programs. Standard version is supplied with an 8-channel washing head for dispensing/aspiration, three bottles for washing and rinsing solutions, a waste bottle and bottle with filter. Optional 4-channel washing solution weight logger, 4 CHW Logger is available.

#### The unit provides:

- Washing mode;
- · Rinsing mode;
- · Mixing mode;
- · Single point, two point aspiration;
- Possibility of additional solution mixing during time gap between two work cycles;
- Possibility to use microtest plates by different manufacturers, ensured by automated plate set up (adjusting to different depths of plate wells);
- · Plate and strip washing mode;
- · User-defined programs with adjustable parameters;
- · Saving work programs.





ORDERING INFORMATION:

Cat. number

IW-8 BS-060106-AAI

4 CHW Logger BS-060102-AK

Choice of 3 washing liquid bottles

Minimum dispense volume

Dimensions (W×D×H)

Weight with accessories

External power supply

Consumed power

25 µl

375 × 345 × 180 mm

9.6 kg

22 W

DC 12 V, 5 A

### **IW-8,** Intelispeed Washer



4-channel washing solution weight logger, 4 CHW Logger provides automatic control of rinsing solutions and waste volume. The washer shows the remaining volume for each bottle as a percentage and gives a warning message in case of low solution volume or full waste bottle when 4 CHW Logger is connected.



#### 4 CHW Logger Specifications:

Max. loading per scale cup	2 kg
Dimensions (W×D×H)	267 × 252 × 97 mm
Weight	3 kg

Maximum dispense volume	1,600 μΙ	
Dispense increment	25 μΙ	
Dispensing accuracy	±2.5%	
Allowed residual liquid volume not more than 2 µl in plate well		
Number of wells washed simultaneou	sly 8	
Number of washing cycles for each ch	annel 1–15	
Timer sound signal	yes	
Aspiration time	0.2-3 s	
Aspiration/dispensing speed	3 levels	
Max. number of channels in a program	1 2	
Soaking time	0-300 s (increment 10 s)	
Shaking time	0–150 s (increment 5 s)	
Number of washed rows	1-12	
Time of plate single wash (350 µl), not	more 45 s	
Number of programs	101	
Plate platform and washing head mov	vement automated	
Indication of operation modes	8-line LCD	

The unit is designed for use in closed laboratory rooms at temperatures from +4 to +40°C and relative humidity up to 80% at +31°C decreasing linearly to 50% relative humidity at 40°C

### 3D-IW8, Inteliwasher



Inteliwasher **3D-IW8** series microplate washer is designed to wash various types of standard 96-well microtitre plates, microstrips as well as microarrays on FastFRAME (rectangular well shape). It is suitable for washing wells with different bottom shapes: flat, U-shape and V-shape. The unit is fully programmable, ensuring multi-step solution ripening, aspiration (aspiration, combination of aspiration/liquid dispensing and soaking, as well as soaking cycle during a particular time). Dispense system of liquid dosage for each channel separately.

#### The unit provides:

- · Washing mode;
- · Rinsing mode;
- · Mixing mode;
- Single point, two point, circular (circle or rectangular path) aspiration;
- Possibility of additional solution mixing during time gap between two work cycles;
- Possibility to use microtest plates by different manufacturers, ensured by automated plate set up (adjusting to different depths of plate wells);
- · Round-bottom plate and strip washing mode;
- Possibility of user-defined programs with adjustable parameters.



#### **ORDERING INFORMATION:**


Cat. number

3D-IW8


BS-060102-AAI

**4 CHW Logger** 

BS-060102-AK




### 3D-IW8, Inteliwasher



# The unit has 50 programs divided into 5 following aspiration categories (see figure bellow):

- **1 Type 1** (1.0–1.9) **IPF96 U/V** is intended for round and V-shape immunoplates, 1 point aspiration.
- **2 Type 2** (2.0–2.9) **IPF96 FLAT-2** is intended for flat-bottom shape immunoplates, 2 point aspiration.
- **3 Type 3** (3.0–3.9) **IPF96 FLAT-C** is intended for rectangular shape immunoplates, full-circle aspiration direction.
- Type 4 (4.0-4.9) FastFRAME-2 is intended for multi-slide plate\* with rectangular wells, 2 point aspiration.
- Type 5 (5.0–5.9) FastFRAME-C is intended for multi-slide\* plate with rectangular wells, fullsquare aspiration direction.
- \* The **FastFRAME** multi-slide plate or analog plate of another manufacturer, that is compatible with standard  $25 \times 76$  mm (1 × 3 inch) glass slides.



Minimum dispense volume	25 μΙ	
Maximum dispense volume	1,600 μΙ	
Dispense increment	25 μΙ	
Dispensing accuracy	±2.5%	
Allowed residual liquid volume in	n plate well, not more 2 μl	
Number of wells washed simulta	aneously 8	
Number of washing cycles	1–15	
Timer sound signal	yes	
Aspiration time	1-3 s	
Final aspiration time	1-3 s	
Aspiration/dispensing speed	3 levels	
Max. number of channels in a program 2		
Choice of 3 washing liquid bottle	es	
Soaking time	0-300 s (increment 10 s)	
Shaking time	0–150 s (increment 5 s)	
Number of washed rows	1–12	
Time of one plate wash (300 µl),	not more 45 s	
Number of programs	50	
Plate platform and washing hea	d movement automated	
Indication of operation modes	LCD, 8-line	
Dimensions (W×D×H)	375 × 345 × 180 mm	
Weight with accessories	9.9 kg	
External power supply	Input AC 100–240 V 50/60 Hz, Output DC 12 V	
Input current/ power consumpti	ion 12 V, 1.8 A / 22 W	

The unit is designed for use in closed laboratory rooms at temperatures from +4°C to +40°C and relative humidity up to 80% at +31°C decreasing linearly to 50% relative humidity at 40°C.

4-channel washing solution weight logger, 4 CHW Logger, provides automatic control of rinsing solution and waste volumes. The washer shows the volume for each bottle as a percentage and gives a warning message in case of low solution volume or full waste bottle when 4 CHW Logger is connected.

#### 4 CHW LOGGER SPECIFICATIONS:

Max. loading per scale cup	2 kg
Dimensions	267 × 252 × 97 mm
Weight	3 kg

Absorbance

12 V. 5 A / 60 W

### HiPo MPP-96, Microplate Photometer

Microplate Photometer HiPo is a compact tabletop device for measuring optical density – results of ELISA and microbiological studies in 96-well microplates. Photometer is controlled and outputs data via computer. An extensive range of additional interference filters is available (with average increment of 10 nm).

The device is supplied with specialized software

#### QuantAssay. Features of QuantAssay software:

- ELISA assays of any complexity can be carried out via robust assay editor with help of Assay Wizard
- · Quantitative assay, includes up to 20 standards
- · Avidity/Affinity assays
- Multiplex assays with up to 7 assays on one plate
- · Qualitative assay includes up to 11 controls
- · BestFit function for selecting the best calibration curve
- · User-friendly interface: get your results in 3 clicks
- · Save, load and export results
- · LIMS export integration
- · Creates visual reports

Detection mode

SPECIFICATIONS

Light source	LED, self-calibrating
Photodetector	8 silicon photodiodes
Plate type	96-well microplates (including strip-well microplates)
Reading Speed	5–8 s per wavelength
Measurement modes	Endpoint, Kinetic
Measurement channels	8
Reference channel	1
Measurement range	0-4.3 OD
Resolution	0.0001 OD
Wavelength range	400-700 nm
Wavelength selection	up to 8* filters on wheel standard filters 405, 450, 492 and 620 nm
Shaking	4 amplitudes, 4 speeds
Software	QuantAssay / compatible with IDEXX xChekPlus™
PC system requirements	Intel/AMD Processor, 1 GB RAM, Windows 10/11, USB
Overall dimensions (W×D	0×H) 140 × 300 × 130 mm
Weight	4.6 kg
External power supply	Input AC 100–240 V 50/60 Hz, Output DC 12 V

<sup>\* –</sup> It is possible to install up to 4 additional filters on request. Additional filters are available in two specifications: optical absorption not less than 3.5 OD or 4.3 OD

Input current / power consumption





Accuracy (405, 450, 492, 620 nm)		
0.000-2.000 OD	≤ (0.5 % ±0.010 OD) typical	
2.000-3.000 OD	≤ (1 % ±0.010 OD) typical	
Precision / Reproducibility (405, 450, 492, 620 nm)		
0.000-2.000 OD	≤ (0.5 % ±0.005 OD)	
2.000-3.000 OD	≤ (1.0 % ±0.005 OD)	

ORDERING INFORMATION:	Cat. number
HiPo MPP-96	BS-050108-A02
Optional accessories:	

**OD Plate,** Verification tool BS-050108-AK
Additional filters\* On request

### Quant Assay, Software for HiPo MPP-96





Software video is available on the website

ELISA assays of any complexity can be carried out via robust assay editor with help of Assay Wizard:

Heasurement option	5	
Assay name:		
Nosay Namo (28.06	11:38:27)	
Assay type		Wavderigth
○ Quantitative	1 * Post control count	Channel 1
<ul> <li>Qualitative</li> </ul>	1 w Neg. control count	450 nm Channel 2
(i) Avidity	1 w Group count	90 nm Channel 3
○ Multiplex	v Standards count	
		Description
form		

Qualitative assay includes up to 11 controls;

Results can be outputted as Positive/Negative or Positive/Gray Zone/Negative;

Gray zone can be set as symmetric and non-symmetric; Positivity ratio can be outputted




Avidity/Affinity results be outputted as Positive/Negative or Positive/Gray Zone/Negative;

Avidity index margins can be easily set; Avidity Index can be outputted



User-friendly interface: get your results in 3 clicks:

Choose an assay, a template and press Play



Save, load and export results

Creates reports: Excel, PDF, CSV














Quantitative assay includes up to 20 standards; User can choose Standard/Reverse type of curves



BestFit function for selecting the best calibration curve from following models:

4/5 Parameters logistics, Piece-wise linear, Linear, Index/Logarithm/Exponent/Cubic regression models



Install up to 7 assays on one plate by using multiplex

	1	2	3	4	5	6	7	
A	Smp1							
	0	1	2	3	4	5	6	
В	Smp2							
	0	1	2	3	4	5	6	
c	Smp3							
	0	1	2	3	4	5	6	

#### Easy fill of the samples



PDF report contains: Experiment information, Results table, List of variables and it's calculations, Interpretation parameters

								Results					
OH!	Type	Sample Stane	am :	Correspo	00 etil un	Result 1	Result 2	Onen Concentration	Moon Concentration	Calculated Concentration	Meson (CIC)	(Nandard Caucation (CE)	Coefficient of Yorkshore(CO
At	190	580 500			0.000	04.		0.00	10480	1,04 (0)	0.000	0.000	0.10%
42	90	58d 500			0.000	06		0.00	1,24 (0)	1.34 (0)	0.000	0.000	0.00%
AD	71	Emg1			1.390	h-Sarge			49.21 83	44.00 (0.1	1.900	0.096	2.79%
AR	11	(Sing f		1	1.360	9-Range			45.21 (0)	40:30 fU	1.300	0.036	2.79%
AS	79	Smpli			1.915	p. Aurge			40.40 (0	04.30 (0	1,965	0.091	2.71%
AS	79	Smpli			1,016	in Kurge			62:62:60	60/95 RJ	1.865	0.001	2.71%
AT	THE	Sing 17		OT.	1.081	h-Sarge			04.14 (0.1	63-29-97	1.607	0.036	1.62%
Ab	717	(Impt)		17	1.00	9-Range			\$4.14.00	54.99 (0.1	1.967	0.626	1.82%
Ab	125	Simplifi.		25	3.942	Oat of Name			119:SFW/	4F 51 R/I	3.458	0.004	25.00%
A/10	T25.	Smp25		26.	4.300	Out of Range			119:57 W/	105.96 84	3.656	0.004	25.00%
A/I	T10	Smp10		23	889	9-Range			28.47 83	28.47 (0.1	1.00	0.000	0.00%
A/2	130	Simplify		10	480	9-Range			26.47 (0	20.47 (0	100	0.000	0.30%
	91	Stat 91			0.000	046		5.67	2:48 (0)	2:48 (0.)	9.000	0.000	0.10%
80	81	Stel 91			0.000	06.		5.60	2.48 (0.)	2:48.00	9.000	0.000	0.00%
80	13	Smg2		3	1.040	9-Range			38.00 0.1	35 (3.6)	1.110	0.030	2.79%
64	13	Smg2		2	1.140	9-Range			34.06 (0.1	20104 (U	1.710	0.030	2.79%
86	710	Sing18		10	1.080	b Surge			50:41 (6)	50-79 RJ	1.584	1142	2.79%
86	THE	Smp10		10	180	ti-Sarge			52:41.00	61.04 (0.	1.504	0.042	2.79%
87	718	Smg18		18	1.318	9-Range			45.40 %	44.70 (0)	1.560	0.033	1.67%
80	710	Smp18		18	1.361	9-Range			45-46 (U	46 15 (0)	1.940	0.022	1,81%
80	124	Smp26		24	3.160	to Surger			97 84 94	79:54 RJ	2.000	0.700	25.00%
810	T26	Smp26		26	3.980	ti-Surge			97 84 83	125.26 (4)	3.990	0.739	25 00%
in i	134	Simple		34	0.790	P-Range			21 81 82	2010	176	0.000	0.00%
612	T34	Simple		14	0.790	9-Range			21 81 RJ	27 69 FU	1.790	0.000	0.10%
ús.	102	Stel 552			0.160	06		10.00	Fig. 84	Figures.	8.768	0.000	0.00%

## **OD Plate, Verification Instrument for MPP-96 HiPo**




**OD Plate** is the measurement verification instrument for microplate photometer MPP-96 HiPo. The instrument is designed to verify the accuracy and precision of measurements of the photometer at six levels of nominal optical density: 0.3; 0.6; 1.0; 2.0; 3.0; 4.0 OD. The instrument is supplied with the following verification wavelength range: 405-700 nm.

### Instrument is provided in a shockproof container with an USB flash drive containing:

- · Copy of measurement results
- User manual



Nominal optical density levels	0.3; 0.6; 1.0; 2.0; 3.0; 4.0 OD (±0.1 OD)
Verification wavelength range	405, 414, 450, 480, 492, 515, 540, 550, 560, 568, 580, 594, 620, 630, 650, 690, 700 nm
Instrument dimensions	128 × 86 × 12 mm
Net weight	0.3 kg





SPECIFICATIONS

ORDERING INFORMATION:

Cat. number

**OD Plate**, Verification tool

BS-050108-AK



### **General Information**

#### Safety

All Biosan laboratory equipment meets the requirements of International Standard IEC 61010-1: Safety requirements for electrical equipment for measurement, control and laboratory use, and applicable specific parts e.g. IEC 61010-2-010: particular requirements for laboratory equipment for the heating of materials, IEC 61010-2-020: particular requirements for laboratory centrifuges, IEC 61010-2-051: particular requirements for laboratory equipment for mixing and stirring.

#### **CE Mark**

All Biosan laboratory equipment bears a CE mark to indicate that it meets the requirements of all applicable European Directives

Compliance with the Low Voltage Directive is demonstrated by meeting EN 61010 (as indicating in paragraph on safety) and the EMC Directive by meeting EN61326-1: EMC requirements for electrical equipment for measurement, control and laboratory use. Some products also fall within the scope of IVD Directive.

#### **Electrical Supplies**

All standard Biosan laboratory equipment is available for voltages within the range 220–240 V, 50 or 60 Hz. Most of the equipment is also available for voltages 100–120 V, 50 or 60 Hz.

#### Quality

Biosan's core values are reflected in a consistent focus on quality, reliability, and responsibility. The company's Quality Management System is certified according to LVS EN ISO 9001:2015 since 2004, with a scope covering the development, production, sales, and service of laboratory equipment. In 2022, Biosan expanded its certification to include LVS EN ISO 13485:2016. This certification covers the development, design, production, service, and distribution of medical devices, including optical density meters and automatic microplate washers, as well as the sales, storage, and distribution of active and non-active non-implantable medical devices. In 2025, Biosan received certification for its environmental management system according to LVS EN ISO 14001:2015. This certification supports the company's ongoing commitment to sustainability and responsible resource management throughout all operations.

#### **Environmental Conditions**

Biosan laboratory equipment is designed for operation in cold rooms, incubators (excluding  $CO_2$  incubators) and closed laboratory rooms at ambient temperature from  $+4^{\circ}C$  to  $+40^{\circ}C$  in a non-condensing atmosphere and maximum relative humidity 80% for temperatures up to 31°C decreasing linearly to 50% relative humidity at  $40^{\circ}C$ .

#### Guarantee and After Sales Service

Biosan equipment is reliable, designed and built to provide years of trouble-free service. Most Biosan equipment is guaranteed for two years against faulty materials and workmanship and premium product class up to 3 years upon registration in our website support section. Warranty terms and conditions are indicated in the product manual. All Grant standard laboratory equipment is guaranteed for three years against faulty materials and workmanship. Local distributors and service centres provide necessary technical assistance within and outside the warranty period. Biosan technical support team provides direct support offering the best solution for assistance upon receipt of request via e-mail support@biosan.ly or forms available at Technical Support section of Biosan web-site

### World Wide Availability and Support for Biosan Laboratory Equipment

Biosan laboratory equipment and specialist technical support is available world-wide. Please, visit multilingual (English, French, German, Italian, Latvian, Russian, Spanish) web-site <a href="http://www.biosan.ly">http://www.biosan.ly</a> for further product information (videos, brochures, manuals, articles), placing enquiries and locating your locally appointed distributor or contact customer service direct at <a href="mailto:support@biosan.ly">support@biosan.ly</a>.

Biosan has established branch in Dubai bringing products and support closer to customers in the Middle East. Company details: Biosan Middle East L.L.C.-FZ.

As Biosan is committed to a continuous program of improvement, specifications may be changed without notice.

Premium

**Smart Plus** 

**Basic Plus** 

PRODUCT CLASS FEATURES	Basic Plus	Premium	Smart Plus
Designed to complete basic sample preparation tasks	•	•	•
Designed to complete sophisticated sample preparation tasks		•	•
Advanced specifications and special features		•	•
PC interface for logging, control, programming, alarms, online monitoring functions			•
Modern Bioform design	•	•	•
Small footprint	•	•	•
Low power consumption	•	•	•
Safe 12V DC	•	•	•
High quality	•	•	•
2 year warranty + 3rd year purchased via distributors	•		•
2 year warranty + 3rd year for free upon product registration		•	

### **Applications and Articles**



Reverse-Spin® Technology - Innovative Principle of Microbial Cultivation.

biosan.lv/rts-tech





Investigation of the effect of aeration on growth dynamics, respiratory rate and pH changes of the aerobic bacterium *E. coli* BL21 cultivated in RTS-8 PLUS single-use bioreactor.

biosan.lv/rts-growth





Development and evaluation of DNA amplicon quantification.

biosan.lv/uv-box





UVR-M and UVR-Mi, UV Air Recirculators Test Report.

biosan.lv/uvr-test





Germicidal and Antiviral decontamination of air by UV irradiation and UV recirculator method.

biosan.lv/uv-effect





# **HOW TO CHOOSE**

A PROPER SHAKER, ROCKER, VORTEX



Sample volume 10<sup>3</sup> ... 10<sup>2</sup> ml

Erlenmeyer flask and Cultivation flask



Sample volume 101 ml



Sample volume 10º ... 10-3 ml

PCR plates, microtest plates and Eppendorf type tubes





ES-20/80, Orbital Shaker-Incubator



#### **Applications:**

- . Microbiology
- Extraction
- Cell cultivation



RTS-1 and RTS-1C. Personal bioreactor

### Applications:

Multi Bio RS-24,

Programmable rotator

- Microbiology
- Extraction
- Cell cultivation
- Hematology



V-1 plus, Vortex



PST-60HL-4,

PST-60HL,

Thermo-Shaker



PST-100HL, Thermo-Shaker

TS-DW,

Thermo-Shaker for deep well plates





PSU-10i Orbital Shaker

ES-20. Orbital

Shaker-Incubator

**Applications:** Agglutination · Gel staining/destaining



Mini Rocker-Shaker

### Applications:

Multi Bio 3D, Mini Shaker

MSV-3500.

- Nucleic acid Analysis
- Molecular Analysis
- Protein Analysis

Multi Speed Vortex

· Genomic Analysis



**Applications: ELISA Analysis** Genomic Analysis

Hybridization

Immunology

PSU-2T. Mini-Shaker







Centrifuge vortex for PCR plates

TS-100, TS-100C, TS-100C Smart Thermo-Shakers



V-32. Multi-Vortex





#### Applications:

- Agglutination
- Extraction
- Blot hybridisation
- Gel staining/destaining

### Manual and automatic DNA/RNA extraction

**Manual** preparation



**Automated** preparation

Working area

**DNA/RNA UV-cleaner box** 



Resuspension of probes and reagents

MPS-1 High-speed multi

MSC-3000 Centrifuge/vortex multispin



plate shaker

Sample lysis

Sample wash |



Dry block thermostat

Pipette series

Aspirator



Elution of NA



MagSorb 16 Magnetic rack for manual nucleic acid extraction

Microspin 12 Plus High-speed mini-centrifuge



TS-100C Smart

and PCR plates

Programable termo-shaker

with cooling for microtubes

Biomagpure 12 Plus Compact bench-top robotic workstation for automated nucleic acid purification (+ reagents)



**Aspiration** 



FTA-2i Aspirator



Universal aspirator



**Aspirator** 

**Incubation** 



Dry block thermostat



TS-100C **Programable** termo-shaker with cooling

### **PCR** Analysis

Sample resuspendation and droplet spin down



FVL-2400N (tubes)



MSC-2P (PCR-plates)



BioQuant-96





### **Enzyme-linked** immunosorbent assay (ELISA)

Sample preparation • in laminar flow cabinet

**Biological Safety Cabinets class II** 







Incubation •

PST-60HL

PST-60HL-4

PSU-2T







Washing (Manual)





MA-8



Reading and Analysis

HiPo MPP-96 Microplate Photometer with QuantAssay software

### **Microbial Cell Cultivation**



**Biological Safety Cabinet class II** 







CO<sub>2</sub> cultivation



S-Bt Smart Biotherm

Measurement



DEN-1

or



DEN-1B



**DEN-600** 

**Cultivation and Real-time OD Measurements** 

and Logging



RTS-1



RTS-1C



RTS-8 and RTS-8 plus





Biosan, SIA

Ratsupites iela 7 k-2, Riga, LV-1067, Latvia

Phone: +371 67 426 137

E-mail: marketing@biosan.lv

www.biosan.lv

