



Heating mantle w/stirring

Code JGB001, JGB002, JGB003, JGB004, JGB005, JGB006, JGB007, JGB008

656 series, analogical regulation

Code JGB001, JGB002, JGB003, JGB004



- Analogical speed and temperatura adjustment through independent knobs.
- Indicating lights for heating and stirring operation.
- Maximum heating temperature: 380 °C.
- Maximum stirring speed: 1400 rpm.
- Supporting rod and stirring magnet included.

Code	Capacity	Consumption	internal Ø	Depth
JGB001	100 mL	100 W	70	45 mm
JGB002	250 mL	150 W	90	55 mm
JGB003	500 mL	250 W	105	60 mm
JGB004	1000 mL	350 W	140	85 mm

658 series, digital regulation

Code JGB005, JGB006, JGB007, JGB008



- Analogical speed adjustment through potentiometer.
- Digital temperature adjustment through LED display.
- pT100 sensor for temperature control in solution.
- Indicating lights for heating and stirring operation.
- Maximum heating temperature: 380 °C.
- Maximum stirring speed: 1400 rpm.
- Supporting rod and stirring magnet included

Code	Capacity	Consumption	internal Ø	Depth
JGB005	250 mL	150 W	90	55 mm
JGB006	500 mL	250 W	105	60 mm
JGB007	1000 mL	350 W	140	85 mm
JGB008	2000 mL	450 W	180	105 mm