



► Digital oil bath



- 1 | Digital selection of parameters with display of actual temperature and time of operation.
- 2 | Operation in continuous mode or in programmed periods of time with acoustic signal at the end of the selected time.
- 3 | Acoustic signal indicating overheating of the bath when the bath temperature exceeds the selected temperature by 10°C.
- 4 | Stainless steel tank in a single piece free of welds, coupled to an exterior metal cabinet with epoxy coating.
- 5 | Including removable lid with holes and reducing rings of 11, 8 and 5 cm.
- 6 | Supporting tray covering the heating element.

Code	JDD005	JDD006
Model	602/3	602/5
Capacity	2,5 L	5 L
Temperature range	RT+5°C – 200°C	
Accuracy	0,1°C	
Timer	1' - 9999' min	
Approx useful dimensions	145 x 130 x 100 mm	295 x 145 x 100 mm
Holes	1	2

DIGITAL OIL BATH
[J] HEATING/COOLING

Note: manufacturer reserves the right to modify features and photographs of the goods