



AUXILAB BENCHTOP AUTOCLAVES

Code ZCK016, ZCK017, ZCK018, ZCK019

- Small, fully automatic class B autoclaves
- Easy operation and maintenance
- LCD screen
- Interface available in Spanish and English
- Six predefined sterilisation programmes for solid, porous and hollow materials
- Special predefined sterilisation programme for prion inactivation
- User-defined programme
- Programmes for vacuum, BD (Bowie-Dick) and Helix tests
- Cleaning programme
- Option to set a user password
- With printer and USB output
- Electrical and mechanical door locking system
- Electrical protection against short circuits



| Code | ZCK016 | ZCK017 | ZCK018 | ZCK019 |
|------------------------------|---------------------|----------------|----------------|----------------|
| Model | AXL8L | AXL12L | AXL18L | AXL22L |
| Capacity | 8 L | 12 L | 18 L | 22 L |
| Sterilisation temperature | 105-134°C | | | |
| Sterilisation time | 4-60 min | | | |
| Maximum working pressure | 0,21-0,23 MPa | | | |
| Vacuum cycles | 1-10 | | | |
| Drying time | 1-25 min | | | |
| Chamber dimensions | Ø170x320 mm | Ø200x360 mm | Ø250x350 mm | Ø250x450 mm |
| Chamber material | 304 stainless steel | | | |
| Equipment dimensions (LxWxH) | 535x365x320 mm | 535x365x320 mm | 610x500x430 mm | 650x500x430 mm |
| Net weight | 45 kg | 48 kg | 53 kg | 55 kg |
| Power supply | 220V, 50Hz | | | |
| Power | 1,95 kW | | | |