

CAPPRondo 4-Place Plate Shaker is the perfect solution for efficiently and thoroughly mixing up to four individual microplates of various types at once. Its precise 3 mm orbital motion delivers exceptional mixing performance, ensuring consistent results with every use.

### Features

- Shaker is compatible with all commonly used skirted microplates
- Pulse mode for continuous mixing in both clockwise and counterclockwise directions
- Microprocessor controlled for efficient performance even under high load
- Last run feature
- Rubber O-rings ensure an ultimate grip for all four plates
- Digital display gives the user an overview of the speed and time settings
- Easy access to all 4 plates for easy loading and unloading



**2 YEAR  
Warranty**

### Specifications

<b>Speed</b>	200 - 1200 RPM	<b>Digital Display</b>	Yes
<b>Orbital Diameter</b>	3 mm	<b>Platform Size</b>	274 x 188 mm
<b>Motion</b>	Orbital	<b>Protection Class</b>	IP 21
<b>Motor Type</b>	BLDC motor	<b>Noise Level</b>	<60 dB
<b>Time Setting</b>	1 to 99 mins and infinity mode	<b>Dimensions (W x D x H)</b>	267 x 393 x 105 mm
<b>Pulse Mode</b>	Yes, clockwise and counter clockwise	<b>Weight</b>	4.91kg
<b>Capacity</b>	Up to 4 microplates	<b>Power Consumption</b>	36 W

## Ordering Information

Cat No.	Description
CRP-412X	CAPP Rondo 4-Place Plate Shaker 1200 RPM

### Scope of Delivery

1 pc. 4-Place Plate Shaker
1 pc. Power Adapter
1 pc. O-Ring
1 pc. User Manual
1 pc. Warranty Card