

**NEW  
PRODUCT**



Now Portable  
With Optional  
Accessory

**2 YEAR  
Warranty**



**Rotor**

CAPP Rondo Clinical Centrifuge is a compact and reliable low-speed centrifuge, ideal for blood centrifuge applications involving volumes below 15 mL. It stands out in its category for having great features that ensure an easy-to-use, reliable, safe performance.

### Features

- Rotor capacity of up to 8 x 15 mL tubes on a fixed-angle rotor
- Smart airflow design reduces noise
- Precise RPM setting with increments from 100 RPM and setting from 500-6500 RPM
- Built with a long-life brushless DC motor and a microprocessor
  - Can accommodate widely used blood collection tubes due to included adaptors for 5 mL, 7mL, 1.5 mL and 2 mL tubes
  - Digital display with speed and timer indicator
  - Automatic imbalance detection
- Lid lock safety – lid opens automatically after the operation is completed
- Functions – countdown timer / RPM/RCF speed setting / error display

### Specifications

<b>Speed</b>	500 - 6500 RPM / 3873 x g (max. RCF)
<b>Speed Accuracy</b>	±100 RPM
<b>Motor Type</b>	Brushless DC Motor
<b>Time Setting</b>	1 to 30 min and infinite mode
<b>Imbalance Detection</b>	Yes
<b>Noise Level</b>	< 60 dB
<b>Acceleration Time</b>	43 ± 2 sec.
<b>Deceleration Time</b>	40 ± 2 sec.
<b>Display</b>	Yes
<b>Dimensions (W x D x H)</b>	260x 244 x 203 mm
<b>Weight</b>	4kg (with rotor)
<b>Power Consumption</b>	100 W

### Scope of Delivery

1 pc. Clinical Centrifuge
1 pc. 8-slots Fixed Angle Rotor (Pre-installed)
1 pc. Power Adapter
1 pc. T-Allen Key
8 pc. Reduction Adaptors
8 pc. Tube Holder and Rubber Cap
1 pc. User Manual
1 pc. Warranty Card

## Ordering Information

Cat. No.	Description
CRC-658	CAPP Rondo Clinical Centrifuge 6500 RPM / 3.873 x g for 8 x 15 mL
CRC-658-0006	Portable Cable