

CAPP Rondo Mini Vortex Mixer offers the ideal solution for thorough mixing and quick vortexing, delivering high performance in terms of speed, reliability and safety.

### Features

- Microprocessor-controlled ensures delivery of the maximum speed, even under varying load conditions
- Variable speed is adjusted by turning speed control button
- Highly touch-sensitive switch ensures real-time mixing
- Suspension design for longer switch life
- Low-profile design with a small footprint, making it ideal for laboratory benches
- In-built counterbalance system creates maximum vortexing action, while minimising noise and excessive vibration
- Brushless DC motor for safety, quick acceleration and low maintenance
- Ideal for cold rooms and incubators, it maintains a set speed between 5-40°C
- Chemically-resistant PC ABS molded exterior, with aluminium die-cast base



**2 YEAR**  
Warranty

### Specifications

<b>Speed Range</b>	2000- 4500 RPM	<b>Protection Class (DIN EN 60529)</b>	IP 43
<b>Speed Accuracy</b>	± 5 %	<b>Material</b>	Body- PC ABS / Platform- Die Cast
<b>Motor Type</b>	Brushless DC motor	<b>Operating Mode</b>	Touch
<b>Max. Load Capacity</b>	500 g	<b>Dimensions (W x D x H)</b>	100 x 100 x 70 mm
<b>Max Tube Diameter</b>	30 mm	<b>Weight</b>	0.55kg
<b>Motion</b>	Orbital	<b>Power Consumption</b>	15 W
<b>Orbital Diameter</b>	4.5 mm		

### Scope of Delivery

- 1 pc. Mini Vortex Mixer
- 1 pc. Inbuilt Power Adapter
- 1 pc. User Manual
- 1 pc. Warranty Card

## Ordering Information

Cat. No.	Description
CRV-45X	CAPP Rondo Mini Vortex Mixer 4500 RPM w/ adjust. speed