



Eye washer and shower eye washer

Code RDB004, RDB005, RDB006

Showers and eye wash of emergency are a necessary security measure with high utility in most of the companies and laboratories. They can reduce the damage caused by the splashing in the eyes, face and the body with acids, bases, organic compound and other toxic or corrosive components, as well in case of fire.

It is the easiest and cheapest way of decontamination in case of accident.

For a correct use and to obtain the higher benefit it is necessary to follow some basics rules of localization, installation and signposting.

Eye washer

Code RDB005, RDB006

- The eyewash is in a pedestal which base is fixed in the floor to obtain a big stability.
- It has two sprays with a big volume of water with low pressure that can bath the eyes and the face in case of accident.
- The valve has automatic opening when you push the kick starter.
- It keeps opening till you push the lever in the other direction to close it.

Code	Material
RDB005	Stainless steel
RDB006	ABS plastic tank and PVC-coated steel pipe



Shower eye washer

Code RDB004

- Shower to be use in case of emergency that is combined with the eyewash in a pedestal.
- It is made in stainless resistant to corrosion.
- The shower head as a big diameter, it gives a big volume of water and a circular water curtain.
The eyewash has two sprays with a big volume of water with low pressure that can bath the eyes and the face of the person that has been injured.
- The shower operation is by a rigid handle and the eye wash has a lateral manual opening that works pressing the lever.

Code	Material
RDB004	Stainless steel 304

