



Mini Centrifuge RS-1008

Code GDC007

- It cannot operate if the lid is open
- Ideal for microfiltration and quick spin-downs
- Easy-to-use with start/stop function
- Fixed speed, 7000 rpm
- Two rotors included; for microtubes and PCR strips
- Rotor change without tools
- Low noise level
- Acceleration and braking in seconds



Technical specifications

Code	GDC007
Speed	7000 rpm
Max. RCF	2680×g
Capacity	
Angle rotor for microtubes:	8 x 0,2/0,5/1,5/2 mL
Rotor for PCR strips:	4 strips of 8 x 0,2 mL
Run time	Continuous operation
Motor	DC brush motor
Power supply	External AC/DC adapter:
Input	AC 100-240V, 50/60Hz
Output	DC 24V
Noise level	≤45dB
Dimension [D×W×H]	160×170×122 mm
Weigh	0,5 kg

Compatible tubes

Tube	Use of adapter	Codes
Microtube 0,5 mL	With adapter	BGN004
Microtube 1,5 mL	Without adapter	BGN005, BGN011, BGN012, BGN013, BGN014, BGN015, BGN027
Microtube 2 mL	Without adapter	BGN025
Strip PCR 8x0,2 mL	-	BGN006
Microtube 0,2 mL	With adapter	BGN003, BGN026