

D1012UA

High Speed Palm Micro Centrifuge

Features

- Max.Speed 10000rpm (7490xg) ;
- OLED display, may switch speed and time;
- DC motor drive, smooth and quiet operation;
- Safe,endurable and robust rotor with clamp locking design, easy to operate;
- Wide voltage technology: the speed is not affected by voltage fluctuations, high speed accuracy;
- Compatible with 0.5/1.5/2mL × 12 tubes+4 × PCR8 strips, 0.2mL × 44 PCR tubes;
- Screen buttons can adjust the start/stop of the motor, supports both soft braking and instant stop modes;
- Simple appearance and structure, only weight 1.5Kg;



Application

This centrifuges is ideal for quick spins,microfiltration of samples,cell separations and many other routine laboratory procedures.



Specifications

D1012UA

Max.Speed	10000rpm
Max.RCF	6708 × g-1.5/2mL; 6596-7490 × g-PCR8 strips;
Rotor Capacity	0.2mL × 44 PCR tubes 0.5/1.5/2mL × 12+ 0.2mLPCR 8 strips × 4
Run Time	10s-20min
Display	OLED
Acceleration/Deceleration	1 ↑ / 2 ↓
Driving Motor	DC motor
Power	AC100-240V/50Hz/60Hz 40W
Noise Level	≤58dB
Dimension (D × W × H)	210 × 180 × 130mm
Weight	1.5kg



DLAB Scientific Inc.
18150 Rowland st. City of
Industry, CA 91748



DLAB Scientific Co., Ltd
Yu An Road 31, Beijing Airport
Economic core zone,
Shunyi District, Beijing 101318 China
info@dlabsci.com



info@dlabsci.com



facebook.com/dlabsci/



linkedin.com/company/dlab-scientific-inc/



tinyurl.com/DLABYouTube



twitter.com/DLabScientific



service@dlabsci.com