



- 1 | User-friendly and intuitive moisture analyzer with 5-inch high-definition touch screen.
- 2 | Silent double action opening. It has an RS232 port to connect a printer. Possibility to print up to 9 groups of data records stored in memory.
- 3 | Equipped with a level gauge and two adjustable leveling feet to compensate for slight unevenness on the operating table.

Halogen Moisture Analyzer



Code	KCD001
Capacity	100g/1mg
Readability	50g/5mg
Repeatability	1mg
Operation Temperature	5°C-40°C
Pan size	100mm
Heating Element	Halogen
Temp.Range	50°C-180°C
Temp.Accuracy	1°C
Moisture Range	0%-100%
Moisture Readability	0.01%
Dry weighing range	100%-0%
Dry weighing Readability	0.01%