



► Dry Bath RS-120S



- Temperature control up to 120°C
- With timer function
- The device is supplied without blocks. Eight optional blocks are available for tubes of different volumes and microplates
- Overheating protection
- Lid for protection and heat preservation

Technical specifications

Code	JDS007
Functions	Heating
Regulator temperature range	15°C - 120°C, increment: 1°C
Temperature control accuracy (at 37°C)	± 0,5°C
Temperature uniformity	± 0,5°C
Max. heating rate	5,5°C/min
Timer	1 min – 99 h 59 min, continuous operation (00:00)
Screen	LED
Overheating protection	140°C
Power supply	AC 100-120V/220-240V, 50/60Hz
Power	160W
Dimension [D×W×H]	175×290×85 mm
Weight	3 kg

Available heating blocks

Code	Description
JDV017	Block 54x0,2mL
JDV018	Block 40x0,5mL
JDV019	Block 40x1,5mL
JDV020	Block 40x2mL
JDV021	Block 28x5/15mL
JDV022	Block 8x50mL
JDV023	Block for 96/384 microplate
JDV024	Block for 0.2mL, 0.5mL, 1.5/2mL tubes, 18 holes each volume