



Microporous plate oscillation mixer, model 730

Code LGD003

- Microplate shaker with compact appearance, easy operation, safe, stable and quiet. 3-inch LCD screen, real-time display of system operating status and parameters
- Speed and time controlled by microprocessor. Capacity for 4 plates or microplates. Equipped with power outage recovery function, the instrument can automatically resume operation according to the original set program after power outage recovery
- Practical fixing bracket for quick and easy installation of microplates
- Maintenance free brushless motor.



Technical Specifications

Code	LGD003
Speed range	200-1500 rpm, (increments of 10)
Speed accuracy	±10 rpm
Operating modes	Continuous, pulse
Storage	Up to 10 programs
Plate diameter	3 mm
Timer	1s-999 min.
Maximum capacity	4 microplates
Time to maximum speed (upwards)	10 s
Time to maximum speed (downwards)	5 s