



Pipette controller

Code - MBB002, MBB004, MBB007, MBG018

- The pipette controllers collection let you choose for the most effective solution in every application, helping you to get the most comfortable and fast manipulated of liquids through the use of graduated or volumetric pipettes.
- The pipette controllers as well as the pipette pumps, ergonomically designed, it have an internal mechanism that allows suck and dispense different volumes of liquids without change the manual position just with your thumb, avoiding the hand fatigue.
- With flexible connection and manufactured in rubber, they can be used with pipettes of different capacities.
- They are resistant to the action of acids and bases, and easily removable for cleaning.

Manual controllers

Code - MBB002, MBB004, MBB007

- Ergonomic, lightweight and easy to handle with one hand
- With pipette fillers, suction and expulsion lever of liquid and button for expulsion of the last residual fluid
- For pipettes of 0.1 mL to 100 mL
- Includes replaceable hydrophobic 3 μ m filter to protect the internal mechanism.



Code	Color
MBB002	Blue
MBB004	Green
MBB007	Yellow

Accessories included

Code - MBG018

Code	Description
MBG018	Filter for DIGIPETTE manual pipette controller