



Ultra-flat magnetic stirrer RSLAB 1F - The new ultra-flat magnetic stirrer of our RSLAB range is particularly suitable for the following applications: physical/chemical, analysis and biotechnology laboratories.



#### — Features

- 1 | Thickness : only 15 mm.
- 2 | Automatic reverse rotation for better mixing results.
- 3 | High IP protection class : IP65.
- 4 | Set-up plate is made from chemically resistant materials.
- 5 | Slip-proof and stable.

#### — Technical specifications

Code	LBG004
Model	1F - Flat Design
Speed range	15-1500 rpm
Max. stirring volume	2 kg
Speed adjustment	By +/- buttons
Work light	Yes
Ø. work plate	90 mm
Heating	No
Dimensiones (LxIxH)	175X120X15 mm
Weight	0,30 Kg

## MAGNETS AND RODS [L] STIRRING

Note: manufacturer reserves the right to modify features and photographs of the goods.