



► Magnetic stirrers RSLAB-3C



- 1| Operating temperature from room temperature up to +550°C.
- 2| Heating stop system if the safety temperature is exceeded.
- 3| Precise temperature and speed control: Microprocessor technology guarantees temperature and speed accuracy.
- 4| Includes Pt1000 temperature probe (support kit not included).
- 5| LED display for temperature control.

Code	LBD003
Display	LED
Speed range	0-1500 rpm
Maximum volume	10 L
Temperature range	RT - 550°C
Security temperature	580°C
Plate material	Glass ceramic
Plate dimensions	184x184 mm
Approx. dimensions	360 x 215 x 112 mm
Approx. dimensions	4,5 Kg

► Accessories (not included)



LBM002 Temperature probe
PT1000



LBM012 Poly-block blue stand
with 4 holes for 30 mL tubes



LBM005 Temperature probe
PT1000 with glass coating



LBM013 Poly-block black stand
with 4 holes for 40 mL tubes



LBM003 Support kit for probe PT1000

LBM004 Support kit for probe PT1000



LBM014 Poly-block green stand
with 6 holes for 80 mL tubes



LBM010 Poly-block red stand with
11 holes for 4 mL tubes



LBM015 Poly block yellow stand with 4
holes with 16 mL tubes



LBM011 Poly-block purple stand
with 4 holes for 20 mL tubes