




## Stainless steel vacuum manifold kit - 1 position

Code. PJB018

Vacuum manifold kit consisting of:

1. Vacuum manifold with 1 position
2. Stainless steel filter holder
3. Stainless steel funnel 300mL
4. Adapter with stainless steel valve
5. Aluminium clamp
6. Filters PTFE 5 $\mu$ m

Code	Description
PJB018	Stainless steel vacuum manifold kit, 1 position




## Stainless steel vacuum manifold kit - 3 positions

Code. PJB019

Vacuum manifold kit consisting of:

1. Vacuum manifold with 3 positions
2. 3 stainless steel filter holders
3. 3 stainless steel funnels 300mL
4. 3 adapters with stainless steel valve
5. 3 aluminium clamps
6. Filters PTFE 5 $\mu$ m

Code	Description
PJB019	Stainless steel vacuum manifold kit 3 positions



## Stainless steel vacuum manifold kit - 6 positions

Code. PJB020

Vacuum manifold kit consisting of:

1. Vacuum manifold with 6 positions
2. 6 stainless steel filter holders
3. 6 stainless steel funnels 300mL
4. 6 adapters with stainless steel valve
5. 6 aluminium clamps
6. Filters PTFE 5 $\mu$ m

Code	Description
PJB020	Stainless steel vacuum manifold kit, 6 positions