



Tube roller mixers

Code **LJD003** - 683/20, **LJD004** - 683/30

Nahita-Blue roller mixers offer high performance and versatility, making them ideal for applications such as haematological analysis and the homogeneous mixing of different types of samples in clinical and research laboratories.

Features

- High-precision rollers made of silica gel ensure smooth and uniform operation.
- Optimised internal structure with exclusive anti-blocking design for quieter and more stable performance, minimising the risk of blockages or failures during use.



Code	LJD003	LJD004
Model	683/20	683/30
Type	Rotary agitation	Horizontal agitation
Rollers	6	4
Containers	Tubes	Reagent bottles
Maximum capacity	10x50 ml (0.5 kg)	2x1L (2kg)
Permitted container dimensions	Ø 6-34 mm	Ø 45-100 mm
Control mode	Knob	Digital
Motor	Brush motor	Brushless
Speed	20-80 rpm	
Time	-	1 min-23h59min
Power	25 W	
Voltage	100-220V 50/60 Hz	
Dimensions	440x255x105 mm	