



VORTEX STIRRER MODEL 740

Code LGB009

The Nahita-Blue 740 vortex mixer is compact and efficient, ideal for basic laboratory applications. Designed to provide uniform and reliable agitation in volumes of up to 50 ml, this device is suitable for tasks in various fields.



Features

- **Compact size:** small design, easy to handle
- **Control type:** Pressure operation (touch) for simple use
- **Agitation orbit:** 4.5 mm
- **Maximum agitation capacity:** 50 ml
- **Motor type:** Brushless, maintenance-free motor

Technical Specifications

Code	LGB009
Model	Model 740
Control	Constant speed
Base material	Zinc alloy
Orbital diameter	4.5 mm
Maximum agitation volume	50 ml
Speed	3000 rpm
Operating conditions	5-40 °C / 80% HR
Protection level	IP42
Power	12 W
Power supply	100-240 V 50/60 Hz
Dimensions	93x93x72 mm
Weight	0.40 Kg