



Vacuum pump B-10

Code PXB017

- Compact design with metal housing
- Diaphragm pump with overtemperature or overload protection sensors
- Suction cups on the base of the unit and carrying handle
- 8 mm diameter hose connector (olive)
- Includes 2 rubber tubes (ID 6 mm; OD 12 mm) and 2 microfilters



Technical specifications

Reference	PXB017
Model	B-10
Type	Negative pressure
Adjustable pumping speed	NO
Pumping speed	10 l/min
Ultimate vacuum/negative pressure	0.075 Mpa
Noise	≤40 db
Power	20 W
Diaphragm	Rubber
Overload protection	YES
Overtemperature protection	YES
Short circuit protection	YES
Voltage	100-240V 50/60 Hz
Dimensions	192x100x135 mm
Weight	1.58 kg