

VIAL RACKS

The stackable vial racks made of transparent acrylic or blue PP allow the easy handling and secure transport of sample vials.

Transparent acrylic allows the fill level of the sample vials to be easily checked. PP is not transparent, but more robust than acrylic.

VIAL RACKS, ACRYLIC



Dimensions W x H x D mm	Positions	For	PK	Art. no.
173 x 20 x 95	50	ND8*	1	7.621 486
173 x 20 x 95	50	ND8**, ND9, ND10, ND11	1	7.620 863
175.8 x 20 x 115.5	40	ND13, 4 ml shell vials	1	7.616 896
160 x 30 x 160	25	ND18, ND20	1	7.620 864

* ND8 crimp neck vials and micro-vials

** ND8 screw neck vials, small opening



VIAL RACKS, PP

Dimensions W x H x D mm	Positions	For	PK	Art. no.
200 x 17 x 105	50	ND8*, ND9, ND10, ND11	1	7.654 509
230 x 28 x 117	50	ND13, 4 ml shell vials	1	7.671 275

* ND8 screw neck vials

THE VIAL RACKS ARE SUITABLE FOR SAMPLE VIALS
WITH THE FOLLOWING ARTICLE NUMBERS:

ND8*	ND8**, ND9, ND10, ND11			ND13, 4 ml shell vials	ND18, ND20
7.612 017	6.088 871	7.620 898	7.660 024	7.603 252	6.204 710
7.614 045	6.205 647	7.621 171	7.660 048	7.613 421	7.612 926
7.616 829	6.258 862	7.621 337	7.670 623	7.616 808	7.613 328
7.616 830	6.401 175	7.621 467	7.672 235	7.616 870	7.613 394
7.616 831	6.901 405	7.622 228	7.672 236	7.616 880	7.615 808
7.616 832	6.901 955	7.626 843	7.970 595	7.632 226	7.615 908
7.616 878	7.608 132	7.626 899	7.677 377	7.648 254	7.616 139
7.617 518	7.608 141	7.629 622	7.677 368	7.648 518	7.616 552
7.618 925	7.608 160	7.630 175	7.663 231	7.648 602	7.616 883
7.620 436	7.612 960	7.631 401	7.663 233	7.654 504	7.616 895
7.621 750	7.613 087	7.631 402	7.663 232	7.654 505	7.616 899
7.622 387	7.613 330	7.631 774	7.645 571		7.620 147
7.631 599	7.613 388	7.631 798	7.677 365		7.620 148
7.622 388	7.615 163	7.632 401	7.663 226		7.620 798
7.630 552	7.615 291	7.636 093	7.677 379		7.620 815
7.615 704	7.615 715	7.639 156	7.663 230		7.621 127
7.677 195	7.616 003	7.639 476	7.663 229		7.621 813
7.654 498	7.616 019	7.639 477	7.664 225		7.630 303
	7.616 109	7.639 478	7.677 376		7.632 402
	7.616 839	7.643 512	7.645 592		7.648 101
	7.616 848	7.644 559	7.663 234		7.648 634
	7.616 849	7.647 475	7.645 633		7.648 635
	7.616 850	7.647 476	7.677 375		7.663 432
	7.616 851	7.647 477	7.648 512		
	7.616 859	7.647 478	7.648 513		
	7.616 860	7.647 479	7.648 514		
	7.616 861	7.647 480	7.648 515		
	7.616 879	7.648 146	7.648 516		
	7.616 932	7.648 519	7.648 517		
	7.616 934	7.648 520	7.648 521		
	7.616 935	7.648 597	7.648 522		
	7.618 897	7.648 599	7.648 523		
	7.618 914	7.651 116	7.648 524		
	7.620 828	7.654 554	7.648 525		
	7.620 829	7.655 281	7.648 526		

* ND8 crimp neck vials and micro vials

** ND8 screw neck vials, small opening