

# Heating mantle | beaker

for beaker, integrated controller, up to 450°C



## WHM for beaker, 450°C integrated controller

### Ideal for:

- for beaker 100 - 5,000 ml

### Features:

- 100°C - 450°C
- full choice of vessel, capacity, housing and control are offered
- nickel chrome heating element
- coil type: acid-resistant, long life time
- all kinds of beaker can be used (Duran®/Kimax/Pyrex/Iwaki)
- CE certified and unique serial number for tracing

### Controller:

- built-in *temperature controller*



1.000 ml beaker heating mantle

Application (beaker brand)	Beaker capacity (ml)	Beaker Ø (mm)	Watts (W)	Inside depth (mm)	Mantle outside Ø x H mm	Support rod Ø x H mm	Order number
							with 230V
common-use	100	50 - 52	80	40	Ø 145 x 105	Ø 13 x 440	DH.WHM12141
common-use	250	68 - 70	140	65	Ø 145 x 140	Ø 13 x 440	DH.WHM12142
for DU.	400	79 - 81	240	70	Ø 145 x 140	Ø 13 x 440	DH.WHM12143
for Iw./Ki./Py.	400	76 - 78	240	75	Ø 145 x 140	Ø 13 x 440	DH.WHM121431
common-use	600	88 - 90	240	75	Ø 145 x 140	Ø 13 x 440	DH.WHM12144
common-use	800	98 - 100	300	90	Ø 190 x 165	Ø 13 x 440	DH.WHM12145
common-use	1,000	105 - 108	350	90	Ø 190 x 165	Ø 13 x 440	DH.WHM12146
common-use	2,000	130 - 132	450	125	Ø 190 x 180	Ø 13 x 660	DH.WHM12147
for DU./Iw.	3,000	151 - 153	550	160	Ø 270 x 235	Ø 13 x 660	DH.WHM12148
for Ki./Py.	4,000	160 - 162	700	160	Ø 270 x 235	Ø 13 x 660	DH.WHM121481
for DU.	5,000	169 - 171	800	160	Ø 270 x 235	Ø 13 x 660	DH.WHM12149

Remarks: DU. = Duran®, Ki. = Kimax, Iw. = Iwaki, Py. = Pyrex