


**! Aparát KDA - CASSOFIX umožňuje
PLYNULÚ KOMPRESIU A DISTRAKCIU FRAGMENTOV
pri vzájomnej rovnobežnej polohe prvej a ľavej časti rámu !**

**! The apparatus KDA - CASSOFIX enables
CONTINUOUS COMPRESSION AND DISTRACTION OF FRAGMENTS
at the mutual concurrent position of the right and left part of the frame !**



**distrakčný rozsah - 150 mm
distraction range - 150 mm**

**Aparát KDA - CASSOFIX
Apparatus KDA - CASSOFIX**

Obj. č. **00045**
Order. No. **00045**



**Istiacie prvky - Kirschnerove drôty
Securing elements - Kirschner's wires**



FERRMAT
Kurská 20
040 22 Košice
Slovak republic

Tel/Fax +421 55 6714 178
Mobil +421 907 768 421
+421 905 837 841
E-mail cassofix@cassofix.sk

IČO 10694609
DIČ 1024244298
Http www.cassofix.sk

