

JB Academy, Basic unstirred water baths

DESCRIPTION

An ideal choice for schools and colleges requiring a basic simple-to-use quality water bath. Base tray included as standard. A great value range consisting of three models – 5L, 12L and 18L.

- Ambient +5°C to 95°C
- Unique Set and Forget® technology – fast heat-up, reliable temperature control
- Stability $\pm 0.5^{\circ}\text{C}$
- Simple, intuitive controls – quick and easy to set temperature
- Practical front panel lock – disables front panel controls preventing unintentional temperature changes

Applications

- Practical science demonstration, sample warming, media preparation
- Markets: Education, schools, colleges and industry

SPECIFICATIONS

Tank capacity	5 / 12 / 18 L
Temperature range	ambient +5 to 95°C
Temp. display and setting resolution	0.5°C
Temp stability (DIN 12876) @ 70°C	$\pm 0.5^{\circ}\text{C}$
Nominal operating voltage	120 or 230 V



CAT. NUMBER

JBA5	230VAC 50/60Hz Euro plug
JBA12	230VAC 50/60Hz Euro plug
JBA18	230VAC 50/60Hz Euro plug